

# **Trust but always Verify:**

---

**Best practices in spectrophotometric measurement**

**L S Theobald Lecture  
15<sup>th</sup> May 2003**

# A journey in analytical science

---

- Sourced in classical analysis (ABC)
- Inspired by qualitative analysis
- Progressed via colorimetry
- Matured in spectrophotometry
- Continued quest for correctness and traceability in spectroscopic measurements

## Sourced in classical analysis

---

- ICI General Chemicals Division Widnes 1962-1967
- 'Techniques Course' Widnes Tech. 12 months before you were 'acceptable' within the analytical laboratory
- Classical gravimetric and volumetric analysis (Over 150 procedures against 'unknown' samples from Vogel 3<sup>rd</sup> Edition)
- Trained analyst

# Traceability in classical analysis

---

## Gravimetric analysis


- Traceable masses
- Method of swings on 2 pan balances
- Substitution
- Buoyancy correction

## Volumetric analysis

- Calibration of pipettes
- Calibration of burettes with drainage correction
- Calibration of volumetric flasks

Chemical Analysis –The Working Tools  
Analytical Chemists' Committee of ICI  
2<sup>nd</sup> Edition 1961

# Traceability in volumetric analysis



# Qualitative analysis

---


- The flame test; low cost atomic emission (Flame and bead tests data have been restored in the 84<sup>th</sup> Edition, 2003-2004, CRC Handbook of Chemistry and Physics)
- Calibration of spectrometers
  - 'Raies Ultimes'
  - Dispersion curves

# Bunsen spectroscope

---



# The eye as a detector



Approximate wavelength range (nm)	Colour of reflected light	Colour of absorbed light
400-456	Violet	Yellow-green
465-482	Blue	Yellow
482-487	Greenish-blue	Orange
487-493	Blue-green	Red-orange
493-498	Bluish-green	Red
498-530	Green	Red-purple
530-559	Yellowish-green	Reddish-purple
559-571	Yellow-green	Purple
571-576	Greenish-yellow	Violet
576-580	Yellow	Blue
580-587	Yellowish-orange	Blue
587-597	Orange	Greenish-blue
597-617	Reddish-orange	Blue-green
617-780	Red	Blue-green

# Absorptiometry; the eye as a detector

---

- Nessler tubes
  - Iron by thioglycollic acid
  - Standard series
  - Visual interpolation
- Nessleriser Glass standards
  - 'Permanent' standard
  - Traceable (?)
  - Samples only
- Lovibond tintometer

# Nessleriser standards

---



# Absorptiometry; the photocell as a detector

---

- Hilger Spekker absorptiometer
- Measurement of transmittance
- Calibration curves and the problem of absorbance
- Filter selection
  - Measurement
  - Band width
  - Stray light


# Hilger Spekker Absorptiometer

---




# Hilger Spekker Absorptiometer

- Light from lamp G illuminates 2 photocells A & O via optical null
- F & H are heat absorbing filters
- M is the wavelength filter
- I is the variable window
- The sample is in the cell L



# Band-pass filters for filter photometers




- The chromogenic reagent yields an orange solution (D)
- Filter A will be satisfactory
- Filters B & C transmit considerable radiation where the solution has little absorption

# Loughborough University of Technology


---

- 1967 to 1972
- The Photoelectric Spectrometry Group Bulletin
- Basic dyes as chromogenic reagents
- Elwell Ward 1971
- Society for Analytical Chemistry

# Hilger & Watts Uvispeck



# Unicam SP 600



# Basics of spectrophotometry

---

- Theoretical aspects
- Nuts & bolts
- Confidence in the numbers!
- J R Edisbury (1966)  
‘Practical hints on absorption spectrometry 200 -800m $\mu$ ’

## Beer Lambert Law Assumptions

---

- ✓ The radiation is perfectly monochromatic
- ✓ There are no uncompensated losses due to scattering or reflection
- ✓ The radiation beam strikes the cuvette at normal incidence
- ✓ There are no molecular interactions between the absorber and other molecules in solution
- ✓ The temperature remains constant

## Two More Assumptions

---

- ✓ **The instrument works properly**
- ✓ **The person operating it is competent**

# Beer's Law and Instrument Linearity

---

The voltage generated by the detector is a function of the power of the incident radiation falling on it;  $V = I^\beta$

Beer's law

$$I = I_0 e^{-\varepsilon cb} \text{ or } A = \log\left(\frac{I_0}{I}\right) = \varepsilon cb$$

Combining these two equations we get

$$V = I^\beta = (I_0 e^{-\varepsilon cb})^\beta = I_o^\beta e^{-\varepsilon \beta cb}$$

and hence the measured absorbance  $A$  is

$$\log\left(\frac{V_0}{V}\right) = \varepsilon \beta cb$$


and therefore a proportional change in  $c$  or  $b$  only confirms that  $\beta$  is constant and does not confirm unity

## Nuts and bolts


---

- Know your instrument
- Importance of optical design


# Pye Unicam SP500



# Double beam recording design



# HP 8450 diode array spectrometer



# The rise of spectrophotometry

---

- Photoelectric Spectrometry Group
  - 1948
  - Beckman model DU & Unicam SP 500
  - UVSG
  - 50<sup>th</sup> anniversary in 1998

# UVSG monographs

---

1. Standards in Absorption Spectrometry  
Ed. C. Burgess & A. Knowles, 1981
2. Standards in Fluorescence Spectrometry  
Ed. J. N. Miller, 1981
3. Practical Absorption Spectrometry  
Ed. A. Knowles & C. Burgess 1984
4. UV Spectroscopy; Techniques, instrumentation,  
data handling  
Ed. B.J. Clark, T. Frost & M.A. Russell, 1993
5. Standards and Best Practice in Absorption  
Spectrometry  
Ed. C. Burgess & T. Frost, 1999

## UVSG/CORM conferences

---

1. Advances in Standards and Methodology  
in Spectrophotometry  
Ed. C. Burgess & K.D. Mielenz, 1987
2. Spectrophotometry, Luminescence and  
Colour;  
Science and Compliance  
Ed. C. Burgess & D.G. Jones, 1995
3. Optical Spectrometry: Applications and  
Instrumentation into the 21<sup>st</sup> Century  
1998


# Spectrophotometry in a regulated environment

---

Glaxo Laboratories Ltd 1974

- Monitoring the performance of UV-visible spectrophotometers;  
UV Spectrometry Group Bulletin 5 77, 1997
- Control charts
- Standards for accuracy and precision
  - Ho(III) & Sm(III)
  - Thompson's solution
  - Sap Green solutions

# Ho(III) as a wavelength standard



# Qualification and calibration of spectrophotometers

---

Basis for measurement


- **Wavelength scale**
  - Operational range
  - Accuracy & precision
- **Photometric scale**
  - Operational range
  - Accuracy & precision (Linearity)
  - Noise and drift
- **Resolution**
- **Stray light**

# Traceable standards



# Starna RM-RE

---




## Long term photometric stability

---

- How reliable is the instrument over its working life?
- How stable is the standard?
- How stable is the standardising laboratory?

# HP 8450 life time photometric accuracy


- NPL calibrated neutral density filter
- 546nm checked daily
- Nine years (110 months)
- 2161 data points
- Monthly averages
- 95% CLs



# Potassium dichromate


- 2 maxima/minima pairs
- One close to the visible
- One in 'bomb alley' of UV region
- Easy to prepare (60mg in 0.005M sulphuric acid)

European Pharmacopoeia 4.5 2003




# HP 8450 photometric accuracy

- acid dichromate
- Four years  
(47 months)
- Duplicates
- Satisfactory performance
- Reading slightly low at 350nm
- No correlation max/min




# HP 8450 photometric accuracy

- acid dichromate
- Four years  
(47 months)
- Duplicates
- Less satisfactory performance in UV
- Correlation max/min



# HP 8450 photometric accuracy

- acid dichromate
- Four years  
(47 months)
- Duplicates
- Maxima to maxima
- No correlation



# Technology transfer standards

---


- Stable
- Non hazardous
- UV & Visible range
- Easy to use

# Sap Green 1979

---



# Sap Green Linearity




# Comparison of scanning and DAS spectrophotometers 1990


---

- 7 DAS (HP8450 & HP8452)
- 6 Scanning instruments (PE Lambda 5, 6 and 15)
- All labs had the 'same' solution of Sap Green

# Collaborative trials with Sap Green



# Collaborative trials with Sap Green



# Spectrophotometric methods

---

The instrument is OK but what about the method?

Optimum conditions

- Best precision
- Best accuracy

# Precision and noise

---

The relative standard deviation of concentration

$$RSD_c = \frac{s_c}{c} = \frac{0.4343s_T}{T \log_{10} T}$$

however  $s_T$  is a function of  $T$

Three main sources;

1. Electronics, which is a constant (Johnson Noise)

$$s_T = k_1$$


2. Detector, photon statistics (Shot noise)

$$s_T = k_2 \sqrt{T^2 + T}$$


3. Cuvette positioning errors and source flicker

$$s_T = k_3 T$$

# Precision and noise



# Relative errors; HP 8450 diode array spectrometer




## Accurate values of absorbance

---


- Setting standard values
- Minimising bias
- Less out of specification results
- Forgotten lessons
  - Ringbom-Ayres plots
  - VFG method

# Vandenbelt, Forsyth and Garrett Ind. Eng. Chem. (Analytical Edition) 1945

- variation of molar extinction coefficient with absorbance
- Impact on standard values
- Beckman DU working absorbance range 0.5 to 1.9 for broad bands



# VFG plot



## In conclusion

---

- Spectrophotometry is alive if not well
- Standards are becoming more important
- Trust your instruments and your methods BUT always verify
- The journey continues...

# NIR wavelength for the future?

TS5 & SRM 2035 (offset)  
8<sup>th</sup> January 2003; Jasco V570 1nm SBW; 100nm min<sup>-1</sup>

