EXHIBITION CHEMISTRY

SPECTACULAR DEMONSTRATIONS FOR TEACHERS

DEMONSTRATIONS TO CAPTURE THE STUDENT'S IMAGINATION, BY DECLAN FLEMING OF PATE'S GRAMMAR SCHOOL, CHELTENHAM IN THIS ISSUE: CATALYTIC COPPER

Heterogeneous catalysis

Curricula for 14–16 and post-16 students usually involve some reference to heterogeneous catalysis in the gas phase. The 'go-to' demonstration for this is the oxidation of ammonia.¹ Ammonia is presumably chosen for the curriculum connection to the Haber process. However, the demonstration requires the use of platinum wire, which many schools have only in very small quantities due to its cost.

The oxidation of acetone (propanone) is a similar demonstration, which requires no additional oxygen supply. It involves readily available copper that can be used in larger quantities to dramatic effect.

Kit

- Copper foil (0.25 mm or thicker) or coin
- Acetone
- 250 cm³ conical flask
- Bunsen burner or blowtorch

Procedure

Place a 250 cm³ conical flask on a heat resistant mat. Add acetone to the flask to a depth of 5 mm (approximately 30 cm³ of acetone). Trim a sheet of copper foil and bend it to hook over a glass rod. Check that when the copper foil is lowered into the flask it is held approximately 2 cm above the base.

In front of the audience

Heat the copper directly until it is red hot, then lower it into the conical flask. Waves of colour will ripple across its surface as it catalyses the oxidation of the acetone. The effect is emphasised if the lights are turned off and the copper will continue to glow as long as there is a supply of acetone vapour.

Teaching goal

The objective in this experiment is to demonstrate heterogeneous catalysis. The copper catalyses





🔯 Catalytic copper – see how to perform this demonstration at www.rsc.li/EiC0512ec

the formation of ketene (ethenone) and methane.

 $(CH_3)_2C=O \rightarrow H_2C=C=O + CH_4$

The ketene and methane are oxidised (and presumably some acetone is directly oxidised as well), providing the necessary heat to keep the copper glowing.

The role of oxygen in the reaction can be shown by covering the top of the flask with a heat resistant mat: the copper will quickly cool and dim. The effluent gases have been shown by infrared spectroscopy² to mostly consist of carbon dioxide and water, with small amounts of ketene.

If you find the copper does not

REFERENCES

- 1. T Lister, Classic chemistry demonstrations, No. 100. Cambridge, UK: Royal Society of Chemistry, 1996 (http://bit.ly/NemGp4)
- 2. R Whitaker and T McGarian, J. Chem. Educ., 1976, **53**, 776 (DOI: 10.1021/ ed053p776



Safety

- Wear eye protection.
- The flask will get very hot during the reaction.
- The oxygen in the flask is quickly absorbed by the primary reaction. However, keep an additional heat resistant mat or damp cloth nearby in case the acetone auto-ignites.
- Ketene is a respiratory irritant: although most will be oxidised, the demonstration should be performed in a well ventilated area.

sustain the reaction, it is likely that its position in the flask is incorrect to allow efficient mixing of the acetone and oxygen – try raising or lowering the copper. Thin copper foil used for electrodes often cools too quickly to maintain the reaction. It is possible to leave the flask reacting for the whole lesson and weigh the copper at the end to show that none has been consumed in the reaction.

An interesting point to draw the audience towards is the yellow colour of the copper that is seen when it is heated in the reducing gas mixture. In the absence of any oxides on the copper surface, its colour is closer to that of gold.

Alternative methods

This demonstration can also be performed with a coil of copper wire or with a copper coin held in a cradle of nichrome wire. The latter has the advantage of a good heat capacity to ensure the reaction is sustained. However, it is less impressive for showing the waves of colour changes on the surface of the copper as it finds itself under alternating oxidising and reducing conditions.