

Supporting information for

Per- and polyfluoroalkyl substances and fluorine mass balance in cosmetic products from the Swedish market: Implications for environmental emissions and human exposure

Lara Schultes,^{*a} Robin Vestergren,^{*b} Kristina Volkova,^c Emelie Westberg,^b Therese Jacobson,^c
Jonathan P. Benskin^{*a}

a Department of Environmental Science and Analytical Chemistry (ACES), Stockholm University, Stockholm, Sweden

b Swedish Environmental Research Institute (IVL), Stockholm, Sweden

c Swedish Society for Nature Conservation (SSNC), Stockholm, Sweden

*Corresponding authors:

Lara.Schultes@gmail.com

Robin.vestergren@ivl.se

Jon.Benskin@aces.su.se

Contents

Reagents and standards	S3
Mass spectrometer settings	S3
Figure S1. Inter-laboratory comparison of target PFASs analysis	S4
Table S1. Target compounds and internal standards	S5
Table S2. Mobile phase gradient profile for PFASs, PFCAs, FOSA, FOSAA, FTSA and alternatives.	S6
Table S3. Mobile phase gradient program for PAPs.	S7
Table S4. LC-MS/MS method accuracy and method precision.....	S8
Table S5. Absolute recoveries of internal standards.....	S9
Table S6. CIC method accuracy and method precision	S10
Table S7. PFASs concentrations (ng/g) in cosmetic samples	S11
Table S8. PAPs concentrations (ng/g) in cosmetic samples	S12
Table S9. C_{F_TF} , C_{F_EOF} , ΣC_{F_PFAS} , $C_{F_extr.unknown}$ and $C_{F_non\ extr.}$ concentrations in $\mu\text{g/g}$	S13
References	S14

Reagents and standards

Recovery standards $^{13}\text{C}_8$ -PFOA (MS transition 421>376) and $^{13}\text{C}_8$ -PFOS (MS transition 507>80) were purchased from Wellington (Ontario, Canada). For all other standards refer to Table S1. Fluoride standard (1000 mg/l) was obtained from Thermo Scientific. Certified reference material (BCR®-461) was obtained from Sigma-Aldrich. Solvents and reagents used for extraction and analytical procedures: Methanol (LiChrosolv®, Merck), acetonitrile (Chromasolv™, Honeywell), ammonium acetate (Merck), 1-methylpiperidine (Merck), sodium hydroxide (Sigma Aldrich), hydrochloric acid (37 %, VWR), acetic acid (Sigma Aldrich), MilliQ water (TOC = 3 ppb, conductivity = 18.2 MΩ; Millipore, Merck) and Supelclean™ ENVI-Carb™ SPE Bulk Packing (Supelco). Argon and oxygen gases were of purity grade 5.0.

Mass spectrometer settings

The mass spectrometer was operated with the following settings adapted from Gebbink et al.¹: source temperature = 150 °C, desolvation temperature = 350 °C, capillary voltage = 3kV, cone gas (nitrogen) flow = 150 l/hr, desolvation gas (nitrogen) flow = 650 l/hr, nebuliser gas = 7 bar. Details on MS transitions can be found in table S1.




Figure S1. Inter-laboratory comparison of target PFASs analysis carried out at ACES and IVL. Overall correlation coefficient $r^2 = 0.9605$.

Table S1. Target compounds and internal standards, selected instrumental parameters for quantification by UPLC/MS/MS.

Abbreviation [□]	Precursor Ion	Quantitative Product ion	Qualitative product ion	Internal standard	IS transition	Native standard used for quantification	Data quality*	Supplier [#]
PFBA	213	169	149	¹³ C ₄ -PFBA	217>172	PFBA	1	W
PFPeA	263	219	169	¹³ C ₅ -PFPeA	266>223	PFPeA	1	W
PFHxA	313	269	119	¹³ C ₂ -PFHxA	315>270	PFHxA	1	W
PFHpA	363	319	169	¹³ C ₄ -PFHpA	367>322	PFHpA	1	W
PFOA	413	169	369	¹³ C ₄ -PFOA	417>372	PFOA	1	W
PFNA	463	419	219	¹³ C ₅ -PFNA	468>423	PFNA	1	W
PFDA	513	469	269	¹³ C ₂ -PFDA	515>470	PFDA	1	W
PFUnDA	563	519	269	¹³ C ₂ -PFUnDA	565>520	PFUnDA	1	W
PFDoDA	613	569	169	¹³ C ₂ -PFDoDA	615>570	PFDoDA	1	W
PFTTrDA	663	619	169	¹³ C ₂ -PFTTrDA	615>570	PFDoDA	2	W
ADONA	377	251	85	¹³ C ₄ -PFOA	417>372	ADONA	2	W
9Cl-PF3ONS	531	351	83	¹³ C ₄ -PFOS	503>80	9Cl-PF3ONS	2	W
11Cl-PF3OUdS	631	451	83	¹³ C ₄ -PFOS	503>80	11Cl-PF3OUdS	2	W
PFBS	299	80	99	¹⁸ O ₂ -PFHxS	403>84	PFBS	2	W
PFHxS	399	80	99	¹⁸ O ₂ -PFHxS	403>84	PFHxS	1	W
PFHpS	449	80	99	¹³ C ₄ -PFOS	503>80	PFHxS	3	
PFOS	499	80	99	¹³ C ₄ -PFOS	503>80	PFOS	1	W
PFNS	549	80	99	¹³ C ₄ -PFOS	503>80	PFOS	3	
PFDS	599	80	99	¹³ C ₄ -PFOS	503>80	PFDS	2	W
PFUnDS	649	80	99	¹³ C ₄ -PFOS	503>80	PFDS	3	
FOSA	498	78	478	¹³ C ₈ -FOSA	506>78	FOSA	1	W
4:2 FTSA	327	307	81	¹³ C ₂ -6:2 FTSA	429>409	4:2 FTSA	2	W
6:2 FTSA	427	407	81	¹³ C ₂ -6:2 FTSA	429>409	6:2 FTSA	1	W
8:2 FTSA	527	507	81	¹³ C ₂ -6:2 FTSA	429>409	8:2 FTSA	2	W
FOSAA	556	498	419	D ₃ -MeFOSAA	573>419	FOSAA	2	W
4:2 monoPAP	343	97	323	¹³ C ₂ -6:2 monoPAP	445>97	4:2 monoPAP	2	D
6:2 monoPAP	443	97	423	¹³ C ₂ -6:2 monoPAP	445>97	6:2 monoPAP	1	W
8:2 monoPAP	543	97	523	¹³ C ₂ -8:2 monoPAP	545>97	8:2 monoPAP	1	W
10:2 monoPAP	643	97	623	¹³ C ₂ -8:2 monoPAP	545>97	10:2 monoPAP	2	C
4:2/4:2 diPAP	589	343	97	¹³ C ₄ -6:2/6:2 diPAP	793>445	4:2/4:2 diPAP	2	D
4:2/6:2 diPAP	689	443	343	¹³ C ₄ -6:2/6:2 diPAP	793>445	6:2/6:2 diPAP	3	
6:2/6:2 diPAP	789	443	97	¹³ C ₄ -6:2/6:2 diPAP	793>445	6:2/6:2 diPAP	1	W
6:2/8:2 diPAP	889	443	543	¹³ C ₄ -6:2/6:2 diPAP	793>445	6:2/8:2 diPAP	2	W
8:2/8:2 diPAP	989	543	97	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	1	W
6:2/10:2 diPAP	989	443	643	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	3	
8:2/10:2 diPAP	1089	543	643	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	3	
6:2/12:2 diPAP	1089	443	743	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	3	
8:2/12:2 diPAP	1189	543	743	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	3	
6:2/14:2 diPAP	1189	443	843	¹³ C ₄ -8:2/8:2 diPAP	993>545	8:2/8:2 diPAP	3	

□ For full chemical names refer to Buck et al²

* 1 = native standard and exact matched mass labelled standard, 2 = native standard but no exact matched mass labelled standard, 3 = no native standard, no exact matched mass labelled standard

W = Wellington Laboratories (Guelph, ON, Canada), C = Chiron Chemicals (Hawthorn, VIC, Australia), D = Donated by Dr. Xenia Trier (National Food Institute, Denmark)

Table S2. Mobile phase gradient profile for PFSAs, PFCAs, FOSA, FOSAA, FTSA and alternatives.

Time [min]	LC Gradient Program		LC Flow Rate
	Mobile phase A [%]¹	Mobile Phase B [%]²	[ml/min]
0.0	90	10	0.40
0.5	90	10	0.40
5.0	20	80	0.40
5.1	0	100	0.40
6.6	0	100	0.40
8.0	0	100	0.55
10.0	90	10	0.40

¹ Mobile phase A: 90 % water and 10 % acetonitrile containing 2 mM ammonium acetate.

² Mobile phase B: 100 % acetonitrile containing 2 mM ammonium acetate.

Table S3. Mobile phase gradient program for PAPs.

Time [min]	LC Gradient Program		LC Flow Rate
	Mobile phase A [%]¹	Mobile Phase B [%]²	[ml/min]
0.0	80	20	0.30
4.0	0	100	0.30
7.5	0	100	0.30
9.0	80	20	0.30

¹ Mobile phase A: 95 % water and 5 % methanol containing 2 mM ammonium acetate and 5 mM 1-methylpiperidine (1-MP).

² Mobile phase B: 5 % water, 75 % methanol and 20 % acetonitrile containing 2 mM ammonium acetate and 5 mM 1-methylpiperidine (1-MP).

Table S4. LC-MS/MS method accuracy as % recovery (spike recovery), method precision as % RSD (n = 3). Low fortification amount was 1 ng per analyte and sample, high fortification 48 ng.

Analyte	Low fortification		High fortification	
	% recovery	RSD	% recovery	RSD
PFBA	100	9%	90	7%
PFPeA	93	17%	86	6%
PFHxA	93	17%	92	11%
PFHpA	161	14%	162	11%
PFOA	97	11%	92	11%
PFNA	109	4%	98	3%
PFDA	99	17%	101	24%
PFUnDA	78	25%	85	4%
PFDoDA	87	8%	89	12%
PFTTrDA	71	46%	47	7%
ADONA	114	31%	90	59%
9Cl-PF3ONS	280	27%	216	32%
11Cl-PF3OUdS	188	23%	186	34%
PFBS	101	16%	93	5%
PFHxS	92	17%	90	15%
PFOS	111	63%	112	32%
PFDS	83	26%	93	38%
FOSA	120	38%	94	17%
4:2 FTSA	142	22%	82	16%
6:2 FTSA	120	37%	77	13%
8:2 FTSA	124	49%	126	10%
FOSAA	29	57%	40	10%
4:2 monoPAP	-		80	14%
6:2 monoPAP	87	15%	72	9%
8:2 monoPAP	222	26%	119	37%
10:2 monoPAP	-		70	99%
6:2/6:2 diPAP	85	15%	102	17%
6:2/8:2 diPAP	59	21%	57	21%
8:2/8:2 diPAP	92	48%	186	65%

Table S5. Absolute recoveries of internal standards, determined by recovery standards M8PFOA and M8PFOS. Recoveries are given as averages from triplicate high and low spike experiments. n.d. = not determined.

Recovery standard	M8PFOA		M8PFOS	
	Average recovery (%)	RSD	Average recovery (%)	RSD
MPFBA	83%	11%	86%	21%
M3PFPeA	125%	16%	129%	25%
MPFHxA	109%	13%	113%	24%
MPFHpA	63%	10%	65%	20%
MPFHxS	109%	15%	111%	14%
MPFOA	104%	14%	106%	18%
MPFNA	105%	7%	108%	16%
MPFOS	80%	27%	80%	18%
MPFDA	53%	21%	54%	18%
MPFUnDA	85%	7%	88%	17%
MPFDoDA	52%	23%	53%	23%
M8FOSA	47%	27%	48%	24%
13C 6:2 monoPAP	150%	15%	173%	12%
13C 8:2 monoPAP	80%	27%	92%	29%
13C 6:2/6:2 diPAP	n.d.	n.d.	78%	39%
13C 8:2/8:2 diPAP	n.d.	n.d.	21%	118%

Table S6. CIC method accuracy as % recovery (spike recovery), method precision as % RSD (n = 3).

Sample	Fortification [µg]		Expected recovery [%]	Apparent recovery [%]	PFOS recovery*	RSD
	PFOS	NaF				
FOUN07	7.6		100	68	68	12%
FOUN07		5.0	0	-	-	-
FOUN07	7.6	5.0	50	36	72	14%
Blank	7.6		100	60	60	13%
Blank		5.0	0	-	-	-
Blank	7.6	5.0	50	38	76	16%

*PFOS recovery assuming no NaF was extracted.

Table S7. PFASs concentrations (ng/g) in cosmetic samples.

	PFBA	PFPeA	PFHxA	PFHpA	PFOA	PFNA	PFDA	PFUnDA	PFDoDA	PFTrDA	6:2 FTS	8:2 FTS
<i>LOD</i>	0.80	5.35	3.35	0.50	4.38	3.45	3.73	0.45	4.38	3.10	8.85	3.38
CRE01	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE02	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE03	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE04	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE05	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE06	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
CRE07	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN01	382	542	2380	1590	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	21.9	<3.38
FOUN02	437	178	1850	1200	2330	651	1410	269	634	77.6	132	156
FOUN03	609	699	4640	2520	8.38	<3.45	6.34	1.46	<4.38	<3.10	60.6	<3.38
FOUN04	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN05	93.8	69.6	1290	157	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN06	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN07	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN08	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
FOUN09	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
PEN01	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW01	51.0	41.0	50.2	72.0	50.8	47.2	56.2	70.1	102	103	<8.85	<3.38
POW02	<0.80	<5.35	<3.35	1.55	<4.38	<3.45	<3.73	1.96	<4.38	<3.10	<8.85	<3.38
POW03	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW04	5.18	<5.35	3.90	8.99	6.52	7.70	6.32	9.68	13.2	<3.10	<8.85	<3.38
POW05	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW06	38.0	33.3	418	57.3	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW07	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW08	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW09	66.3	34.8	358	44.0	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
POW10	4.78	<5.35	<3.35	5.15	<4.38	<3.45	<3.73	5.19	6.92	<3.10	<8.85	<3.38
POW11	8.57	6.00	8.13	15.4	7.46	10.7	11.4	14.9	23.6	11.2	<8.85	<3.38
POW12	84.9	52.8	447	63.5	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
SHAV01	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38
SHAV02	<0.80	<5.35	<3.35	<0.50	<4.38	<3.45	<3.73	<0.45	<4.38	<3.10	<8.85	<3.38

Table S8. PAPs concentrations (ng/g) in cosmetic samples; n.d. = not detected.

	4:2 monoPAP	6:2 monoPAP	8:2 monoPAP	10:2 monoPAP	4:2/4:2 diPAP	4:2/6:2 diPAP	6:2/6:2 diPAP	6:2/8:2 diPAP	8:2/8:2 diPAP	6:2/10:2 diPAP	8:2/10:2 diPAP	6:2/12:2 diPAP	8:2/12:2 diPAP
<i>LOD</i> *	70.3	3.75	10.0	188	70.0	10.3	10.3	188	190	190	190	190	190
CRE01	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE02	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE03	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE04	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE05	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE06	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
CRE07	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
FOUN01	847	55490	<10.0	<188	691	1380	255570	<188	<190	<190	<190	<190	<190
FOUN02	1640	50450	23450	24330	<70.0	202	63700	38330	11260	5180	891	473	21.6 *
FOUN03	562	61980	38.3	<188	799	1130	405450	682	<190	68.8 *	<190	<190	<190
FOUN04	<70.3	<3.75	<10.0	<188	<70.0	<10.3	11.2	<188	<190	<190	<190	<190	<190
FOUN05	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
FOUN06	<70.3	<3.75	22.2	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
FOUN07	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
FOUN08	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
FOUN09	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
PEN01	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW01	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW02	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW03	<70.3	<3.75	282	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW04	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW05	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW06	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW07	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW08	<70.3	<3.75	145	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW09	<70.3	26.6	<10.0	<188	<70.0	<10.3	27.0	<188	<190	<190	<190	<190	<190
POW10	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW11	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
POW12	<70.3	<3.75	26.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
SHAV01	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190
SHAV02	<70.3	<3.75	<10.0	<188	<70.0	<10.3	<10.3	<188	<190	<190	<190	<190	<190

* Samples FOUN01, FOUN02 and FOUN03 were diluted 1:5 – 1:500 times for PAPs analysis (LOD refers to undiluted samples), therefore concentrations marked with asterisk are reported with lower concentrations than the LOD.

Table S9. C_{F_TF} , C_{F_EOF} , ΣC_{F_PFAS} , $C_{F_extr.unknown}$ and $C_{F_non\ extr.}$ concentrations (\pm standard deviation) in $\mu\text{g/g}$. n.d. = not detected.

sample	C_{F_TF} [$\mu\text{g/g}$]	C_{F_EOF} [$\mu\text{g/g}$]	ΣC_{F_PFAS} [$\mu\text{g/g}$]	$C_{F_EOF} /$ C_{F_TF}	$\Sigma C_{F_PFAS} /$ C_{F_EOF}	$C_{F_extr.unknown}$ [$\mu\text{g/g}$]	$C_{F_non\ extr.}$ [$\mu\text{g/g}$]
CRE01	11100 \pm 2000	<6.65	n.d.	0.0%	-	-	11098
CRE02	4790 \pm 862	<20.44	n.d.	0.0%	-	-	4789
CRE03	4040 \pm 728	<6.65	n.d.	0.0%	-	-	4043
CRE04	2720 \pm 490	<20.44	n.d.	0.0%	-	-	2723
CRE05	2490 \pm 447	11.1 \pm 0.44	n.d.	0.4%	0.0%	11.1	2473
CRE06	<91.1	31.9 \pm 1.27	n.d.	-	0.0%	31.9	-
CRE07*	<91.1	<1.02	n.d.	-	-	-	-
FOUN01	3120 \pm 562	1720 \pm 68.8	196	55.1%	11.4%	1524	1404
FOUN02	2900 \pm 522	1380 \pm 55.1	140	47.5%	10.2%	1237	1522
FOUN03	2570 \pm 463	1050 \pm 42.0	296	40.8%	28.2%	754	1520
FOUN04	2000 \pm 360	552 \pm 22.1	0.01	27.6%	0.0%	552	1449
FOUN05	1660 \pm 299	374 \pm 14.9	1.07	22.5%	0.3%	373	1289
FOUN06	450 \pm 81	5.12 \pm 0.20	0.01	1.1%	0.3%	5.10	445
FOUN07*	326 \pm 58.6	1.29 \pm 0.05	n.d.	0.4%	0.0%	1.29	324
FOUN08*	<91.1	5.98 \pm 0.24	n.d.	-	0.0%	5.98	-
FOUN09*	<91.1	30.2 \pm 1.21	n.d.	-	0.0%	30.2	-
PEN01	438 \pm 78.9	5.76 \pm 0.23	n.d.	1.3%	0.0%	5.76	433
POW01	19200 \pm 3456	<6.65	0.47	0.0%	-	-	19201
POW02	8870 \pm 1596	12.3 \pm 0.49	0.00	0.1%	0.0%	12.3	8854
POW03	4240 \pm 763	289 \pm 11.5	0.17	6.8%	0.1%	289	3948
POW04	3430 \pm 617	4.91 \pm 0.20	0.04	0.1%	0.9%	4.86	3424
POW05*	3330 \pm 600	0.83 \pm 0.03	n.d.	0.0%	0.0%	0.83	3333
POW06	3020 \pm 544	<6.65	0.36	0.0%	-	-	3024
POW07*	2600 \pm 467	<20.44	n.d.	0.0%	-	-	2596
POW08	2570 \pm 463	296 \pm 11.8	0.09	11.5%	0.0%	296	2278
POW09	2460 \pm 443	187 \pm 7.49	0.37	7.6%	0.2%	187	2276
POW10	1770 \pm 318	6.97 \pm 0.28	0.02	0.4%	0.2%	6.96	1758
POW11	1510 \pm 272	6.96 \pm 0.28	0.08	0.5%	1.2%	6.88	1502
POW12	547 \pm 98.4	35.3 \pm 1.41	0.44	6.5%	1.3%	34.8	511
SHAV01	837 \pm 151	<1.02	n.d.	0.0%	-	-	837
SHAV02*	<91.1	3.53 \pm 0.14	n.d.	-	0.0%	3.53	-

References

- 1 W. A. Gebbink, A. Glynn and U. Berger, *Environ. Pollut.*, 2015, **199**, 166–173.
- 2 R. C. Buck, J. Franklin, U. Berger, J. M. Conder, I. T. Cousins, P. De Voogt, A. A. Jensen, K. Kannan, S. A. Mabury and S. P. J. van Leeuwen, *Integr. Environ. Assess. Manag.*, 2011, **7**, 513–541.