

## *Supporting Information*

### **Exploring the utility of PVP@Au-Polyamide-Triton X-114 for SERS tracking of extracellular senescence associated-beta-galactosidase activity**

Shaofei Li,<sup>abc</sup> Yizhuang Cheng<sup>ab</sup>, Miao Qin<sup>ab</sup>, Siyu Chen<sup>ab</sup>, Pan Li<sup>\*ad</sup>, Liangbao Yang<sup>\*abd</sup>


<sup>a</sup> Institute of Health and Medical Technology, Hefei Institutes of Physical Science, Chinese Academy of Sciences, Hefei, Anhui 230031, China

<sup>b</sup> University of Science and Technology of China, Hefei, Anhui 230026, China


<sup>c</sup> School of Life Science, Anhui University, Hefei, Anhui 230601, China

<sup>d</sup> Cancer Hospital, Chinese Academy of Sciences, Hefei, Anhui 230031, China


\* To whom correspondence should be addressed. Email: [lipan2011@iim.ac.cn](mailto:lipan2011@iim.ac.cn); [lbyang@iim.ac.cn](mailto:lbyang@iim.ac.cn)




**Fig. S1** X-gal staining for *E. coli* DH5 $\alpha$  containing pUC18 vector. Colonies were cultured on the Luria-Bertani (LB) solid plate medium containing 1 mmol L<sup>-1</sup> IPTG and 0.1 mmol L<sup>-1</sup> X-gal at 37°C overnight, and they appeared blue by X-gal staining.




**Fig. S2** Phenol-chloroform-isoamyl alcohol extraction. Colonies of *E. coli* DH5 $\alpha$  containing pUC18 vector were cultured in LB liquid medium containing 1 mmol L<sup>-1</sup> IPTG and 0.1 mmol L<sup>-1</sup> X-gal at 37°C overnight. After centrifuging the culture at high speed, the supernatant was collected. (a) The blue supernatant of the culture for *E. coli* DH5 $\alpha$  containing pUC18 vector. (b) Phenol-chloroform-isoamyl alcohol mixture was added in the blue supernatant. (c) The proteins with the blue X-gal hydrolysates were kept at interface between aqueous and organic layer.




**Fig. S3** Adsorption performance of polyamide powder for X-gal hydrolysates. Polyamide powder and triton X-114 were investigated. Blue supernatant by adding amide powder became transparent in more than one hour, and blue supernatant by adding triton X-114 and polyamide powder could effectively adsorb X-gal hydrolysates in 5 minutes without the heating process of cloud point extraction.




**Fig. S4** X-gal staining in liquid medium. Colonies of *E. coli* DH5 $\alpha$  containing pUC18 vector were cultured in LB liquid medium containing 1 mmol L<sup>-1</sup> IPTG and 0.1 mmol L<sup>-1</sup> X-gal at 37°C for 6 hours. Both bacteria and medium could be dyed to the distinct blue color.




**Fig. S5** The relationship between the intensity of the SERS characteristic peak at  $599 \text{ cm}^{-1}$  and the concentration of X-gal hydrolysates for samples b1-b7 in Fig 2. The numbers on the abscissa represent the exponent of 10-fold dilutions.




**Fig. S6** Evaluation on SERS performance of Au@PVP-Polyamide-Triton X-114 for urine SA-β-gal activity. a1-a8 denote eight urine specimens from four healthy children and four patients with chronic kidney disease diagnosed by clinic. The color of X-gal hydrolysates was absent for all specimens in the beginning and gradually appeared at 4 hours. Moreover, the feature Raman band of X-gal hydrolysates at 599 cm<sup>-1</sup> was absent for all specimens in the beginning.




**Fig. S7** SERS spectrum analyzed by Au@PVP for the stained urine specimens at half an hour. The feature Raman band at 599 cm<sup>-1</sup> was also absent for all specimens.




**Fig. S8** Evaluation on SERS performance of Au@PVP-Polyamide-Triton X-114 for malignant effusion SA- $\beta$ -gal activity. a1-a8 denote eight specimens from eight patients with different types of malignancies diagnosed by clinic. The color and feature Raman band of X-gal hydrolysis at 599 cm<sup>-1</sup> were absent for all specimens in the beginning.



**Fig. S9** SERS spectrum analyzed by Au@PVP for the stained malignant effusions SA- $\beta$ -gal activity. The feature Raman band of X-gal hydrolysis at 599 cm<sup>-1</sup> was absent in all at 2 hours.



**Fig. S10** Evaluation on SERS performance of Au@PVP-Polyamide-Triton X-114 for throat secretions SA-β-gal activity. c1-c3 denote three specimens from three patients with respiratory tract infection diagnosed by clinic. The blue band indicates the marker of X-gal hydrolysates at 599 cm<sup>-1</sup>. (A) The feature Raman band of X-gal hydrolysates at 599 cm<sup>-1</sup> was absent for all specimens at the beginning. (B) The color of X-gal hydrolysates was absent for all specimens at 1 hour, and that became obvious after 8 hours. The feature Raman band of X-gal hydrolysates at 599 cm<sup>-1</sup> could be obviously distinguished at 1 hour.



**Fig. S11** Evaluation on SERS performance of Au@PVP-Polyamide-Triton X-114 for nasal secretions SA- $\beta$ -gal activity. d1-d3 denote three specimens from three patients with respiratory tract infection. The blue band indicates the marker of X-gal hydrolysates at  $599 \text{ cm}^{-1}$ . (A) The feature Raman band of X-gal hydrolysis at  $599 \text{ cm}^{-1}$  was absent for all specimens at the beginning. (B) The color of X-gal hydrolysates was absent for all specimens at 1 hour, and that became obvious after 12 hours. The feature Raman band of X-gal hydrolysis at  $599 \text{ cm}^{-1}$  could be obviously distinguished at 1 hour.

**Table S1** Summary table of approximately detected times of X-gal hydrolysates by different methods.

<b>Target molcules</b>	<b>X-gal hydrolysates</b>		
<b>Types of body fluids</b>	Urines	Malignant effusions	Throat and nasal secretions
<b>Developed method</b>	Half an hour	2 hours	1 hour
<b>Colorimetric</b>	4-12 hours	12 hours	8 hours