Electronic Supplementary Material (ESI) for Analytical Methods. This journal is © The Royal Society of Chemistry 2015

Supplementary Material

One-step electrochemical preparation of reduced graphene oxide/poly (sulfosalicylic acid) nanocomposite film for detection of acetaminophen and its application in human urine and serum studies

Bo Zheng, Cong Li, Yaru Li, Yue Gu, Xiaoyi Yan, Liu Tang, Ruixue Chen,

Zhiquan Zhang*

College of Chemistry, Jilin University, Changchun 130012, P. R. China.

* Corresponding author. Tel: +86-431-85168352-7; Fax: +86-431-85168399; *E-mail address*: zzq@jlu.edu.cn (Z. Zhang)




Fig. S1 DPVs of AC with different concentrations (0, 1, 5, 10, 50, 100, 200, 300 μ M) at the PSA/rGO/GCE in 0.1 M PBS (pH 6.0). Inset: the different concentrations of PCT (0, 1, 5 μ M). A step increment = 4 mV; pulse amplitude = 50 mV; pulse width = 60 ms; pulse periods = 200 ms.