

Supplementary data

Manuscript title:

Simultaneous detection of endocrine disrupting chemicals including conjugates in municipal wastewater and sludge with enhanced sample pretreatment and UPLC-MS/MS

Authors:

Bing Zhu, Xiangjuan Yuan, Weiwei Ben, Yu Zhang, Min Yang, Zhimin Qiang*

List of supplementary material:

Material	Title	Page
Table S1	Optional parameters of tandem MS and calibration curves of target EDCs.	S2
Table S2	Intra- and inter-day precisions evaluated at three concentration levels of target EDCs.	S3
Table S3	Major characteristics of the influent and effluent of a local WWTP.	S4
Fig. S1	Sampling points along the A/A/O treatment process in a local WWTP.	S5
Fig. S2	(a) Illustration of E2 fragmentation to produce the 145.18 m/z ion; and (b) mass chromatograms of the product ions of E2 in the MS2 scan mode.	S6
Fig. S3	Total-ion MRM chromatograms of 12 target EDCs (100 µg L ⁻¹ each) in standard solution (prepared in MQ water), unpurified sludge sample, and purified sludge sample.	S7

Table S1

Operational parameters of tandem MS and calibration curves of target EDCs.

Analyte	Formula	Mother ion	Fragmentor (V)	Quantifier (cone energy/eV)	Qualifier (cone energy/eV)	Linear slope ^a	R ²	LOD ^b (pg)
E1	C ₁₈ H ₂₂ O ₂	269.2	160	145.2 (35)	183.2 (45)	478.9	0.9991	2.1
E2	C ₁₈ H ₂₄ O ₂	271.3	170	145.2 (40)	183.2 (40)	514.3	0.9990	8.4
E3	C ₁₈ H ₂₄ O ₃	287.2	170	145.2 (50)	170.9 (34)	466.9	0.9987	9.1
EE2	C ₂₀ H ₂₄ O ₂	295.3	170	145.2 (38)	199.1 (35)	611.0	0.9992	10.5
BPA	C ₁₅ H ₁₆ O ₂	227.3	130	133.0 (30)	211.0 (10)	284.0	0.9982	12.7
NP	C ₁₅ H ₂₄ O	219.3	125	133.1 (25)	117.1 (55)	568.4	0.9966	20.0
E1-3G	C ₂₄ H ₂₉ O ₈ Na	445.2	150	269.1 (35)	113.0 (20)	113.6	0.9997	10.6
E2-3G	C ₂₄ H ₃₁ O ₈ Na	447.2	155	113.0 (20)	271.1 (40)	121.5	0.9999	10.0
E2-17G	C ₂₄ H ₃₁ O ₈ Na	447.2	165	85.0 (35)	271.1 (40)	231.5	0.9991	5.2
E1-3S	C ₁₈ H ₂₁ O ₅ NaS	349.1	110	269.1 (35)	79.9 (55)	751.4	0.9998	1.6
E2-3S	C ₁₈ H ₂₃ O ₅ NaS	351.0	170	271.1 (35)	79.9 (35)	451.5	0.9996	2.7
E3-3S	C ₁₈ H ₂₃ O ₆ NaS	367.0	165	287.1 (40)	79.9 (35)	709.7	0.9966	1.7

^a For each estrogen and BPA, a 10-point calibration curve was established (0.5 to 500 ng L⁻¹); while the concentrations of NP were raised by 10 folds.

^b The instrumental limit of detection was determined by directly injecting serially diluted standards (with MQ water) into the mass spectrometer and then multiplying the concentration at 3:1 of signal to noise ratio by the injection volume (10 µL).

Table S2

Intra- and inter-day precisions evaluated at three concentration levels of target EDCs.

Analyte	Intra-day precision			Inter-day precision		
	Recovery (RSD ^a) (%/%)			Recovery (RSD ^b) (%/%)		
	5 ng L ⁻¹	50 ng L ⁻¹	500 ng L ⁻¹	5 ng L ⁻¹	50 ng L ⁻¹	500 ng L ⁻¹
E1	89.2 (4.8) ^a	86.3 (5.9)	98.8 (2.1)	96.6 (13.0)	86.5 (6.1)	93.0 (10.5)
E2	109.9 (0.8)	110.7 (0.3)	100.7 (0.5)	105.1 (5.9)	91.2 (18.1)	99.8 (2.0)
E3	93.7 (1.5)	95.7 (6.6)	101.7 (1.8)	96.0 (7.6)	97.9 (4.6)	99.2 (7.3)
EE2	90.3 (8.7)	88.1 (8.0)	98.3 (3.7)	87.1 (13.7)	91.3 (10.6)	89.7 (4.6)
BPA	94.1 (2.0)	98.7 (2.9)	106.4 (5.5)	92.6 (3.7)	105.0 (5.6)	107.7 (6.3)
NP (/10)	99.0 (2.5)	99.7 (6.7)	101.4 (3.3)	95.7 (11.3)	88.1 (10.3)	102.2 (9.2)
E1-3G	99.0 (2.9)	89.9 (1.9)	98.6 (4.5)	95.0 (5.6)	94.8 (9.6)	96.4 (2.3)
E2-3G	99.8 (0.9)	94.4 (13.2)	96.7 (2.1)	104.1 (7.1)	84.3 (5.9)	101.7 (4.9)
E2-17G	93.6 (9.5)	101.6 (14.3)	92.0 (8.2)	100.6 (7.4)	105.1 (9.0)	92.7 (2.1)
E1-3S	95.0 (2.9)	98.1 (6.5)	98.5 (1.2)	98.1 (5.7)	95.0 (12.3)	94.3 (8.3)
E2-3S	96.9 (1.3)	94.1 (4.4)	93.6 (0.8)	89.4 (8.7)	100.3 (8.6)	95.1 (1.7)
E3-3S	102.8 (2.5)	87.0 (6.6)	102.3 (2.4)	97.3 (4.5)	89.7 (13.9)	101.7 (4.3)

^a Intra-day repeatability ($n = 3$).^b Inter-day repeatability ($n = 3$).

Table S3

Major characteristics of the influent and effluent of a local WWTP.

Wastewater	COD (mg L ⁻¹)	NH ₄ ⁺ -N (mg L ⁻¹)	TN (mg L ⁻¹)	TP (mg L ⁻¹)
Influent	330.4 (2.8) ^a	18.8 (2.1)	33.5 (5.1)	6.2 (0.7)
Effluent	35.8 (4.4)	<LOD ^b	17.5 (0.3)	0.8 (0.1)

^a Mean concentration (SD), $n = 3$.

^b LOD = 0.4 mg L⁻¹.




Fig. S1. Sampling points along the A/A/O treatment process in a local WWTP.




Fig. S2. (a) Illustration of E2 fragmentation to produce the 145.18 m/z ion; and (b) mass chromatograms of the product ions of E2 in the MS2 scan mode.




Fig. S3. Total-ion MRM chromatograms of 12 target EDCs ($100 \mu\text{g L}^{-1}$ each) in standard solution (prepared in MQ water), unpurified sludge sample, and purified sludge sample.