## **Electronic Supporting Information**

## A fluorescence nanosensor for lipase activity: enzyme-regulated

## quantum dots growth in situ

Wei Zhang<sup>‡a</sup>, Jia Liu<sup>‡a</sup>, Lei Zhang<sup>b</sup>, Jian-hong Gan<sup>a</sup>, Yaqi Ding<sup>a</sup>, Wei Huang<sup>c</sup>, Fengwei Huo<sup>c</sup>, and Danbi Tian<sup>a</sup>\*

<sup>a</sup> College of Chemistry and Molecular Engineering, Nanjing Tech University, 30 South Puzhu Road, Nanjing 211816, China.

<sup>b</sup> Key Laboratory for Organic Electronics and Information Displays & Institute of Advanced Materials (IAM), National Synergistic Innovation Center for Advanced Materials (SICAM), Nanjing University of Posts & Telecommunications, 9 Wenyuan Road, Nanjing 210023, China.

<sup>c</sup> Key Laboratory of Flexible Electronics (KLOFE) & Institute of Advanced Materials (IAM), National Synergistic Innovation Center for Advanced Materials (SICAM), Nanjing Tech University (NanjingTech), 30 South Puzhu Road, Nanjing 211816, China.

\* To whom correspondence should be addressed. Phone/Fax: +86 25 86210899; E-mail: danbi@njtech.edu.cn.

*‡* These authors have contributed equally.




Fig. S1 Distribution of CdS QDs diameter determined by TEM image analysis.






Fig. S2 The zeta potential of the assay system, before and after addition of 2.0 mg mL<sup>-1</sup> lipase.



**Fig.S3** Fluorescence intensity of the assay system in different buffer solutions including a. Tris-HCl; b. KH<sub>2</sub>PO<sub>4</sub>-NaOH; c. Na<sub>2</sub>HPO<sub>4</sub>-citric acid; d.PBS; e. H<sub>2</sub>O (pH=8.0) in the presence of 2 mg mL<sup>-1</sup> lipase.



**Fig.S4** (A) Effects of pH; (B) Effects of temperature; (C) Effects of the substrate concentration(METG); (D) Effects of sodium sulfide. The experiments were performed with 2.0 mg mL<sup>-1</sup> lipase.

**Table.S1** Use of METG as substrate in pH-stat method for detecting the hydrolytic activity of commerical lipases. All assays were performed with 0.02 g lipase at 35°Cand pH 8.0 for 10min.

Lipase	CRL	PPL	SBE-01Li	RNL
V <sub>NaOH</sub> (mL)	7.82	5.60	5.20	4.05
Activity(U/g)	189.24	135.52	125.84	98.01