Electronic Supplementary Material (ESI) for RSC Advances. This journal is © The Royal Society of Chemistry 2015

## Label-free detection of nicotinamide adenine dinucleotide based on ligation-triggered exonuclease III-assisted signal amplification

Haiyan Zhaoa, Lei Wangb, Xingti Liua, Zhiyue Gaoa, Wei Jianga,\*

<sup>a</sup> Key Laboratory for Colloid and Interface Chemistry of Education Ministry, School of Chemistry and Chemical Engineering, Shandong University, 250100 Jinan, P.R. China

<sup>b</sup> School of Pharmaceutical Sciences, Shandong University, 250012 Jinan, P. R. China

Corresponding author: Tel: 86-531-88363888; fax: 86-531-88564464.

E-mail: wjiang@sdu.edu.cn

Table S1. Sequences of oligonucleotides used in this study

Name	Sequence $(5' \rightarrow 3')$
Oligo1	p-TACAAGACACAATAAC
Oligo2	GACGGGAAG
cDNA	GTGTCTTGTACTTCCCGTC
G4 DNA	CTCTTCGAGGGTTTTGGGTTTTTGGGAGCTA
template (40nt)	AAAACCCAAAACCCAAAACCCGTGTCTTGTACTTCCC
	GTC
template (43nt)	CCCAAAACCCAAAACCCGTGTCTTGTACTT
	CCCGTC
template (34nt)	CAAAACCCAAAACCCGTGTCTTGTACTTCCCGTC
Lower-case letter 'p' indicates phosphate group.	

## Optimization of the reaction conditions

To achieve the best performance of our signal amplification strategy, various experimental conditions were investigated, including the concentration of oligo1/oligo 2 DNA, cDNA, and G4-template DNA, the dosage of Exo III and E. *coli* DNA ligase, different pH, and the concentration of NMM.




Fig. S1. The influence of Oligo DNA concentration on the NAD<sup>+</sup> amplified detection system.




Fig. S2. The influence of cDNA concentration on the NAD<sup>+</sup> amplified detection system.




Fig. S3. The influence of G4/template concentration on the NAD<sup>+</sup> amplified detection system.







Fig. S4. The influence of E. *coli* DNA ligase dosage on the NAD<sup>+</sup> amplified detection system.



**Fig. S5.** The influence of Exo III ligase dosage on the NAD<sup>+</sup> amplified detection system.



**Fig. S6.** The influence of different pH on the NAD<sup>+</sup> amplified detection system.



**Fig. S7.** The influence of NMM concentration on the NAD<sup>+</sup> amplified detection system.