

Supporting Information

Quantum dot cluster (QDC)-loaded phospholipid micelles as FRET probe for phospholipase A₂ detection

Junling Li¹, Yonghua Zhang¹, Junjie Ai¹, Qiang Gao^{1*}, Honglan Qi¹, Chengxiao
Zhang¹, Zhiliang Cheng^{2*}

¹Key Laboratory of Applied Surface and Colloid Chemistry, Ministry of Education,
School of Chemistry and Chemical Engineering, Shaanxi Normal University, Xi'an,
710062, China, ²Department of Bioengineering, University of Pennsylvania,
Philadelphia, PA 19104, USA

Corresponding Author

***Qiang Gao**

Tel: 86-29-81530726

Fax: 86-29-81530727

E-mail: gaoqiang@snnu.edu.cn

***Zhiliang Cheng**

Tel: 215-573-0290

Fax: 215-573-2071

Email: zcheng@seas.upenn.edu

Experimental

Synthesis of single quantum dot-loaded phospholipid micelles without NBD


The single quantum dot-loaded phospholipid micelles without NBD was prepared using an oil-in-water emulsion-based self-assembly method. A mixture containing HSPC (0.88 mg) and QD (1 mg) in 150 μ L chloroform was injected into a glass vial containing 3 mL of water, and the sample was sonicated until a homogenous mixture was obtained. The chloroform was then allowed to evaporate overnight. Following that, QDC-loaded phospholipid micelles samples were centrifuged at 1000 rpm for 30 minutes to remove large aggregates. To obtain the QDC-loaded micelle without NBD, the resulting supernatant was centrifuged at 3000 rpm for half hour, and the pellet was resuspended in water (1 mL). To obtain the single QD-loaded micelle without NBD, the resulting supernatant was then centrifuged at 10000 rpm for half hour, and the pellet was resuspended in water (1 mL). The collected samples of QDC- and QD-loaded phospholipid micelles were stored in the dark at 4 °C and used for UV-vis (Figure S2) and fluorescence (Figure S3) measurement.

Synthesis of NBD-containing liposomes


The NBD-containing liposome was prepared using an oil-in-water emulsion-based self-assembly method. A mixture containing HSPC (0.88 mg) and NBD PC (0.1 mg) in 150 μ L chloroform was injected into a glass vial containing 3 mL of water, and the sample was sonicated until a homogenous mixture was obtained. The chloroform was then allowed to evaporate overnight. Following that, NBD-containing liposomes samples were centrifuged at 3000 rpm for 30 minutes to remove large aggregates, the

resulting supernatant was stored in the dark at 4 °C and used for NBD fluorescence measurement (Figure S4).


Chemical information



HSPC



NBD PC



NBD



Hydrophobic CdSe/ZnS QD




Figure S1 TEM image of single quantum dot-loaded phospholipid micelles.




Figure S2 Absorption spectra of single quantum dot-loaded phospholipid micelles
without NBD.




Figure S3. Emission spectra of QDC-loaded phospholipid micelles without NBD. Ex, 460 nm. Aliquots (900 μ L) of 10.0 mM HEPES buffer (pH 7.4) containing 2.0 mM CaCl_2 and samples (50 μ L).




Figure S4. Emission spectra of NBD-containing liposome. Ex, 460 nm. Aliquots (900 μ L) of 10.0 mM HEPES buffer (pH 7.4) containing 2.0 mM CaCl_2 and NBD-containing liposomes (50 μ L).




Figure S5 Linear range of the probe for the PLA₂ detection. Concentration of PLA₂, 5,

10, 25, 50, 100, 200, 300, 400 U/L. Ex, 460 nm, Em, 545 nm.




Figure S6 Dependence of fluorescence intensity on inhibitor LY311727 concentrations. Concentrations of inhibitor LY311727: 0, 10 μM , 25 μM , 50 μM , 80 μM , 100 μM , 200 μM . The concentration of PLA_2 , 50 unit/L.