

checkCIF/PLATON report

You have not supplied any structure factors. As a result the full set of tests cannot be run.

THIS REPORT IS FOR GUIDANCE ONLY. IF USED AS PART OF A REVIEW PROCEDURE FOR PUBLICATION, IT SHOULD NOT REPLACE THE EXPERTISE OF AN EXPERIENCED CRYSTALLOGRAPHIC REFEREE.

No syntax errors found. CIF dictionary Interpreting this report

Datablock: I

Bond precision: Ga- P = 0.0047 Å Wavelength=0.71073

Cell: a=7.294(9) b=18.03(2) c=6.557(8)
 alpha=90 beta=90 gamma=90

Temperature: 296 K

	Calculated	Reported
Volume	862.3(18)	862.0(18)
Space group	I b a 2	Iba2
Hall group	I 2 -2c	?
Moiety formula	Ag Ga5 P8, 4(Ba)	?
Sum formula	Ag Ba4 Ga5 P8	Ag Ba4 Ga5 P8
Mr	1253.55	1253.59
Dx,g cm ⁻³	4.828	4.830
Z	2	2
Mu (mm ⁻¹)	18.472	18.478
F000	1092.0	1092.0
F000'	1091.89	
h,k,lmax	9,23,8	9,23,8
Nref	1002[549]	966
Tmin,Tmax	0.343,0.574	0.215,0.607
Tmin'	0.105	

Correction method= MULTI-SCAN

Data completeness= 1.76/0.96 Theta(max)= 27.510

R(reflections)= 0.0450(882) wR2(reflections)= 0.1211(966)

S = 1.122 Npar= 43

The following ALERTS were generated. Each ALERT has the format
test-name_ALERT_alert-type_alert-level.
Click on the hyperlinks for more details of the test.

● Alert level C

PLAT125_ALERT_4_C	No '_symmetry_space_group_name_Hall' Given	Please Do !
PLAT148_ALERT_3_C	su on the a - Axis is (Too) Large	0.009 Ang.
PLAT148_ALERT_3_C	su on the b - Axis is (Too) Large	0.0200 Ang.
PLAT148_ALERT_3_C	su on the c - Axis is (Too) Large	0.008 Ang.

● Alert level G

PLAT004_ALERT_5_G	Polymeric Structure Found with Dimension	2 Info
PLAT005_ALERT_5_G	No _iucr_refine_instructions_details in the CIF	Please Do !
PLAT083_ALERT_2_G	SHELXL Second Parameter in WGHT Unusually Large.	42.01 Why ?
PLAT301_ALERT_3_G	Main Residue Disorder	Percentage = 11 Note
PLAT328_ALERT_4_G	Possible Missing H on sp3? Phosphorus	P1 Check
PLAT899_ALERT_4_G	SHELXL97 is Deprecated and Succeeded by SHELXL	2014 Note

0 **ALERT level A** = Most likely a serious problem - resolve or explain
0 **ALERT level B** = A potentially serious problem, consider carefully
4 **ALERT level C** = Check. Ensure it is not caused by an omission or oversight
6 **ALERT level G** = General information/check it is not something unexpected

0 ALERT type 1 CIF construction/syntax error, inconsistent or missing data
1 ALERT type 2 Indicator that the structure model may be wrong or deficient
4 ALERT type 3 Indicator that the structure quality may be low
3 ALERT type 4 Improvement, methodology, query or suggestion
2 ALERT type 5 Informative message, check

checkCIF publication errors

🔴 Alert level A

PUBL004_ALERT_1_A The contact author's name and address are missing,
_publ_contact_author_name and _publ_contact_author_address.
PUBL005_ALERT_1_A _publ_contact_author_email, _publ_contact_author_fax and
_publ_contact_author_phone are all missing.
At least one of these should be present.
PUBL006_ALERT_1_A _publ_requested_journal is missing
e.g. 'Acta Crystallographica Section C'
PUBL008_ALERT_1_A _publ_section_title is missing. Title of paper.
PUBL009_ALERT_1_A _publ_author_name is missing. List of author(s) name(s).
PUBL010_ALERT_1_A _publ_author_address is missing. Author(s) address(es).
PUBL012_ALERT_1_A _publ_section_abstract is missing.
Abstract of paper in English.

● Alert level G

PUBL017_ALERT_1_G The _publ_section_references section is missing or
empty.

7 **ALERT level A** = Data missing that is essential or data in wrong format
1 **ALERT level G** = General alerts. Data that may be required is missing

Publication of your CIF

You should attempt to resolve as many as possible of the alerts in all categories. Often the minor alerts point to easily fixed oversights, errors and omissions in your CIF or refinement strategy, so attention to these fine details can be worthwhile. In order to resolve some of the more serious problems it may be necessary to carry out additional measurements or structure refinements. However, the nature of your study may justify the reported deviations from journal submission requirements and the more serious of these should be commented upon in the discussion or experimental section of a paper or in the "special_details" fields of the CIF. *checkCIF* was carefully designed to identify outliers and unusual parameters, but every test has its limitations and alerts that are not important in a particular case may appear. Conversely, the absence of alerts does not guarantee there are no aspects of the results needing attention. It is up to the individual to critically assess their own results and, if necessary, seek expert advice.

If level A alerts remain, which you believe to be justified deviations, and you intend to submit this CIF for publication in a journal, you should additionally insert an explanation in your CIF using the Validation Reply Form (VRF) below. This will allow your explanation to be considered as part of the review process.

Validation response form

Please find below a validation response form (VRF) that can be filled in and pasted into your CIF.


```
# start Validation Reply Form
_vrf_PUBL004_GLOBAL
;
PROBLEM: The contact author's name and address are missing,
RESPONSE: ...
;
_vrf_PUBL005_GLOBAL
;
PROBLEM: _publ_contact_author_email, _publ_contact_author_fax and
RESPONSE: ...
;
_vrf_PUBL006_GLOBAL
;
PROBLEM: _publ_requested_journal is missing
RESPONSE: ...
;
_vrf_PUBL008_GLOBAL
;
PROBLEM: _publ_section_title is missing. Title of paper.
RESPONSE: ...
;
_vrf_PUBL009_GLOBAL
;
PROBLEM: _publ_author_name is missing. List of author(s) name(s).
RESPONSE: ...
;
_vrf_PUBL010_GLOBAL
;
PROBLEM: _publ_author_address is missing. Author(s) address(es).
RESPONSE: ...
;
_vrf_PUBL012_GLOBAL
;
```

PROBLEM: _publ_section_abstract is missing.
RESPONSE: ...
;
end Validation Reply Form

If you wish to submit your CIF for publication in Acta Crystallographica Section C or E, you should upload your CIF via the web. If your CIF is to form part of a submission to another IUCr journal, you will be asked, either during electronic submission or by the Co-editor handling your paper, to upload your CIF via our web site.

PLATON version of 20/08/2014; check.def file version of 18/08/2014

Datablock I - ellipsoid plot



checkCIF/PLATON report

You have not supplied any structure factors. As a result the full set of tests cannot be run.

THIS REPORT IS FOR GUIDANCE ONLY. IF USED AS PART OF A REVIEW PROCEDURE FOR PUBLICATION, IT SHOULD NOT REPLACE THE EXPERTISE OF AN EXPERIENCED CRYSTALLOGRAPHIC REFEREE.

No syntax errors found. CIF dictionary Interpreting this report

Datablock: I

Bond precision: As-Ga = 0.0009 A Wavelength=0.71073

Cell: a=7.4769(8) b=18.5766(19) c=6.7590(7)
 alpha=90 beta=90 gamma=90

Temperature: 296 K

	Calculated	Reported
Volume	938.79(17)	938.79(17)
Space group	I b a 2	Iba2
Hall group	I 2 -2c	?
Moiety formula	Ag As8 Ga5, 4(Ba)	?
Sum formula	Ag As8 Ba4 Ga5	Ag As8 Ba4 Ga5
Mr	1605.15	1605.19
Dx,g cm-3	5.678	5.679
Z	2	2
Mu (mm-1)	30.277	30.277
F000	1380.0	1380.0
F000'	1379.83	
h,k,lmax	9,24,8	9,24,8
Nref	1080[589]	1078
Tmin,Tmax	0.350,0.546	0.065,0.583
Tmin'	0.002	

Correction method= MULTI-SCAN

Data completeness= 1.83/1.00 Theta(max)= 27.500

R(reflections)= 0.0207(1056) wR2(reflections)= 0.0499(1078)

S = 1.172 Npar= 44

The following ALERTS were generated. Each ALERT has the format
test-name_ALERT_alert-type_alert-level.
Click on the hyperlinks for more details of the test.

● Alert level C

PLAT125_ALERT_4_C No '_symmetry_space_group_name_Hall' Given Please Do !

● Alert level G

PLAT004_ALERT_5_G Polymeric Structure Found with Dimension 2 Info
PLAT005_ALERT_5_G No _iucr_refine_instructions_details in the CIF Please Do !
PLAT232_ALERT_2_G Hirshfeld Test Diff (M-X) Ag -- As1 .. 12.5 su
PLAT232_ALERT_2_G Hirshfeld Test Diff (M-X) As1 -- Ga1 .. 5.0 su
PLAT232_ALERT_2_G Hirshfeld Test Diff (M-X) As1 -- Ga2 .. 12.5 su
PLAT301_ALERT_3_G Main Residue Disorder Percentage = 11 Note
PLAT899_ALERT_4_G SHELXL97 is Deprecated and Succeeded by SHELXL 2014 Note

- 0 **ALERT level A** = Most likely a serious problem - resolve or explain
 - 0 **ALERT level B** = A potentially serious problem, consider carefully
 - 1 **ALERT level C** = Check. Ensure it is not caused by an omission or oversight
 - 7 **ALERT level G** = General information/check it is not something unexpected
-
- 0 ALERT type 1 CIF construction/syntax error, inconsistent or missing data
 - 3 ALERT type 2 Indicator that the structure model may be wrong or deficient
 - 1 ALERT type 3 Indicator that the structure quality may be low
 - 2 ALERT type 4 Improvement, methodology, query or suggestion
 - 2 ALERT type 5 Informative message, check
-

checkCIF publication errors

● Alert level A

PUBL004_ALERT_1_A The contact author's name and address are missing,
_publ_contact_author_name and _publ_contact_author_address.
PUBL005_ALERT_1_A _publ_contact_author_email, _publ_contact_author_fax and
_publ_contact_author_phone are all missing.
At least one of these should be present.
PUBL006_ALERT_1_A _publ_requested_journal is missing
e.g. 'Acta Crystallographica Section C'
PUBL008_ALERT_1_A _publ_section_title is missing. Title of paper.
PUBL009_ALERT_1_A _publ_author_name is missing. List of author(s) name(s).
PUBL010_ALERT_1_A _publ_author_address is missing. Author(s) address(es).
PUBL012_ALERT_1_A _publ_section_abstract is missing.
Abstract of paper in English.

● Alert level G

PUBL017_ALERT_1_G The _publ_section_references section is missing or
empty.

- 7 **ALERT level A** = Data missing that is essential or data in wrong format
 - 1 **ALERT level G** = General alerts. Data that may be required is missing
-

Publication of your CIF

You should attempt to resolve as many as possible of the alerts in all categories. Often the minor alerts point to easily fixed oversights, errors and omissions in your CIF or refinement strategy, so attention to these fine details can be worthwhile. In order to resolve some of the more serious problems it may be necessary to carry out additional measurements or structure refinements. However, the nature of your study may justify the reported deviations from journal submission requirements and the more serious of these should be commented upon in the discussion or experimental section of a paper or in the "special_details" fields of the CIF. *checkCIF* was carefully designed to identify outliers and unusual parameters, but every test has its limitations and alerts that are not important in a particular case may appear. Conversely, the absence of alerts does not guarantee there are no aspects of the results needing attention. It is up to the individual to critically assess their own results and, if necessary, seek expert advice.

If level A alerts remain, which you believe to be justified deviations, and you intend to submit this CIF for publication in a journal, you should additionally insert an explanation in your CIF using the Validation Reply Form (VRF) below. This will allow your explanation to be considered as part of the review process.

Validation response form

Please find below a validation response form (VRF) that can be filled in and pasted into your CIF.

```
# start Validation Reply Form
_vrf_PUBL004_GLOBAL
;
PROBLEM: The contact author's name and address are missing,
RESPONSE: ...
;
_vrf_PUBL005_GLOBAL
;
PROBLEM: _publ_contact_author_email, _publ_contact_author_fax and
RESPONSE: ...
;
_vrf_PUBL006_GLOBAL
;
PROBLEM: _publ_requested_journal is missing
RESPONSE: ...
;
_vrf_PUBL008_GLOBAL
;
PROBLEM: _publ_section_title is missing. Title of paper.
RESPONSE: ...
;
_vrf_PUBL009_GLOBAL
;
PROBLEM: _publ_author_name is missing. List of author(s) name(s).
RESPONSE: ...
;
_vrf_PUBL010_GLOBAL
;
PROBLEM: _publ_author_address is missing. Author(s) address(es).
RESPONSE: ...
;
_vrf_PUBL012_GLOBAL
;
```

PROBLEM: _publ_section_abstract is missing.
RESPONSE: ...
;
end Validation Reply Form

If you wish to submit your CIF for publication in Acta Crystallographica Section C or E, you should upload your CIF via the web. If your CIF is to form part of a submission to another IUCr journal, you will be asked, either during electronic submission or by the Co-editor handling your paper, to upload your CIF via our web site.

PLATON version of 20/08/2014; check.def file version of 18/08/2014

Datablock I - ellipsoid plot

