Electronic Supplementary Material (ESI) for Analytical Methods. This journal is © The Royal Society of Chemistry 2016

Disulfide Bond Characterization of Endogenous IgG3 Monoclonal Antibodies Using LC-MS: An Investigation of IgG3 Disulfide-mediated Isoforms.


Jude C. Lakbub¹, Daniel F. Clark¹, Ishan S. Shah², Zhikai Zhu¹, Eden P. Go¹, Thomas J. Tolbert², and Heather Desaire¹

¹Department of Chemistry, University of Kansas, Lawrence, KS, 66047


²Department of Pharmaceutical Chemistry, University of Kansas, Lawrence, KS, 66047


Supplementary Material

Supplementary Figures 1-6: Data supporting the assignment of disulfide bonds in the constant region of IgG3 monoclonal antibodies. Data from the Fitzgerald antibody is shown. The Sigma antibody produced the same results.





Supplementary Figure 1. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the CH3 domain disulfide bond. Asterisks indicate z and $c \pm 1$ ions







Supplementary Figure 2. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the lambda constant light chain disulfide bond






Supplementary Figure 3. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the disulfide bond between the heavy chain and the lambda light chain







Supplementary Figure 4. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the kappa constant light chain disulfide bond




Supplementary Figure 5. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the disulfide bond between the heavy chain and the kappa light chain



Supplementary Figure 6. XIC's (a and b) and ETD spectrum (c) supporting the assignment of the Hinge-4 disulfide bond which links two identical tryptic peptides



Supplementary Figure 7. CID spectrum supporting the assignment of the Hinge-2 disulfide bond



Supplementary Figure 8. CID spectrum supporting the assignment of the Hinge-3 disulfide bond

Supplementary Table 1. Summary of Mass Assignments of Disulfide Bonds in Fitzgerald and Sigma IgG3 mAbs Using High Resolution MS Data.

				Fitzgerald IgG3		Sigma IgG3	
Tryptic dipeptides	Position	Theor. m/z	cs	Exptl. m/z	Mass Error (ppm)	Exptl. m/z	Mass Error (ppm)
DYFPEPVTVSWNSGAL		1581.9839	5+	1581.9953	5	1581.9974	9
TSGVHTFPAVLQSSGL		1318.4878	6+	1318.4954	6	1318.4954	5
YSLSSVVTVPSSSLGT QTYTCNVNHKPSNTK	CH1	1130.2763	7+	1130.2854	8	1130.2886	11
QITICNVNHKPSNIK		989.1177	8+	989.1251	7	989.1251	8
STSGGTAALGCLVK		879.3276	9+	879.3309	4	879.3309	8
CPEPK CPEPK	Hinge-2	572.2623	2+	572.2564	5	572.2602	4
CPEPK			_		_		
SCDTPPPCPR		1070.4394	2+	1070.4419	2	1070.4434	4
	Hinge-3	713.9620	3+	713.9645	4	713.9641	< 1
SCDTPPPCPR		535.7234	4+	535.7264	6	535.7250	3