

Versatile UHPLC-MSMS method for simultaneous quantification of various alcohol intake related compounds in human urine and blood

Rastislav Monošík and Lars Ove Dragsted

Department of Nutrition, Exercise and Sports, Faculty of Life Science, University of Copenhagen, Nørre Alle' 51,
DK-2200 Copenhagen N, Denmark

Supporting Information:

Table S-1. Table overview of LC gradient	2
Table S-2. Overview of mass spectrometry parameters of target compounds with related internal standards.....	3
Figure S-1. Illustrative chromatograms of target compounds and related internal standards	4–6
Table S-3. Coefficients of variation (%) of repeatability (n=6) and intermediate precision (n=9).	7
Figure S-2a. Illustrative chromatograms of isoxanthohumol and isocohumulone in 10× diluted pilsner beer	8
Figure S-2b. Illustrative chromatograms of 3-nitro-tyrosine in a standard solution (1 $\mu\text{g mL}^{-1}$) and pooled sample (n=80).....	9
Figure S-2c. Illustrative chromatogram of indole-3-acetic acid in a standard solution.....	9
Figure S-2d. Illustrative chromatograms of cortisol and cortisol sulphate in standard solution.....	10

Table S–1. Table overview of LC gradient.

Time (min.)	Flow (mL min ⁻¹)	A:B	Curve
0	0.55	98:2	6
1	0.55	98:2	6
1.5	0.43	75:25	6
1.9	0.43	50:50	6
4.1	0.6	20:80	6
4.5	0.6	0:100	3
4.9	0.6	0:100	6
5.1	0.55	98:2	3
6	0.55	98:2	6

Table S–2. Overview of mass spectrometry parameters of target compounds with related internal standards (span = ± 0.2 ; dwell = 0.100 s).

Compound	RT	RT window	Primary transition (quantifier)	Secondary transition (qualifier)	Cone voltage	Collision E
L– (+)–tartaric acid	0.56	0–0.85	149.1→87.1	149.1→73	25	^a 10 / ^b 15
L– (+)–tartaric acid–d ₂	0.55	0–0.85	151.1→74	151.1→88.1	25	15
ethyl sulphate	0.80	0.65–1.15	125→97	125→80	25	20
ethyl sulphate–d ₅	0.79	0.65–1.15	130→98	130→80	25	20
ethyl-β-D-glucuronide	1.39	1–1.65	221.2→85.1	221.2→75.1	25	15
ethyl-β-D-glucuronide–d ₅	1.34	1–1.65	226.2→85.1	226.2→75.1	25	15
indoxyl sulphate	2.54	2.35–2.75	212→80	212→132	30	20
indoxyl sulphate–d ₄	2.53	2.35–2.75	216→80	216→136	30	20
p–cresol sulphate	2.68	2.65–3	186.7→107	186.7→80	30	20
p–cresol sulphate–d ₇	2.80	2.25–2.8	194.1→114.1	194.1→80	30	20
resveratrol	^c 3.04/ 3.26	2.85–3.55	227.1→143.2	227.1→185	30	^a 25 / ^b 20
resveratrol– ¹³ C ₆	^c 3.04/ 3.25	2.85–3.55	233.2→149.2	233.2→191.2	35	^a 25 / ^b 20
estrone 3–sulphate	3.75	3.6–4	349.2→269.4	349.2→80	35	30
estrone–d ₄ –3–sulphate	3.76	3.6–4	353.2→273.4	353.2→80	35	30
DHEAS	4.12	3.9–4.45	367.2→97	–	30	30
DHEAS–d ₅	4.13	3.9–4.45	372.2→98	–	30	30

RT – retention time; ^avalue for primary transition; ^bvalue for secondary transition; ^ccommercial resveratrol is a mixture of *cis*– and *trans*–resveratrol




Figure S-1. Illustrative chromatograms of target compounds and related internal standards.




Figure S-1. Illustrative chromatograms of target compounds and related internal standards.




Figure S-1. Illustrative chromatograms of target compounds and related internal standards.

Table S–3. Coefficients of variation (%) of repeatability (n=6) and intermediate precision (n=9).

Analyte	Repeatability			Intermediate precision		
	Standard solution	Urine	Plasma	Standard solution	Urine	Plasma
L-(+)-tartaric acid	1.7	1.4	2.0	3.5	2.0	3.3
ethyl sulphate	8.7	1.7	11	12	8.2	26
ethyl-β-D-glucuronide	3.3	1.6	2.1	3.9	3.2	2.5
indoxyl sulphate	2.5	2.1	2.7	3.4	1.8	1.9
p-cresol sulphate	4.2	2.1	3.8	4.5	2.3	19
resveratrol	2.1	1.7	1.6	7.5	3.3	3.7
estrone 3-sulphate	1.5	1.8	1.1	2.5	2.0	4.5
DHEAS	0.51	0.97	1.1	0.90	1.3	1.7




Figure S-2a. Illustrative chromatograms of isoxanthohumol and isocohumulone in 10× diluted pilsner beer.




Figure S-2b. Illustrative chromatograms of 3-nitro-tyrosine in a standard solution ($1 \mu\text{g mL}^{-1}$) and an unknown isomer in pooled sample (n=80).




Figure S-2c. Illustrative chromatogram of indole-3-acetic acid in a standard solution.




Figure S-2d. Illustrative chromatograms of cortisol and cortisol sulphate in standard solution.