

## Electronic Supplementary Information

# Fully automated and colorimetric foodborne pathogen detection on an integrated centrifugal microfluidic device

Seung Jun Oh, Byung Hyun Park, Goro Choi, Ji Hyun Seo, Jae Hwan Jung, Jong Seob Choi, Do Hyun Kim and Tae Seok Seo<sup>a\*</sup>

*Department of Chemical and Biomolecular Engineering, Korea Advanced Institute of Science and Technology (KAIST), 291 Daehak-ro, Yuseong-gu, Daejeon 305-701, Republic of Korea.*

**Table S1.** Primer sequence information for the multiplexed LAMP reaction

Target gene / Accession #	Primer	Sequence (5' to 3')	Length [bp]
<i>E. coli</i> O157:H7 fliC gene (CP002967.1 : 2147323-2149017)	FIP	CTACCAACCTGTCTGAAGCGCATCGACATGTTGGACACTTCG	42
	BIP	GGTGAUTGCGGAATCCAGACGTGAGGCAATTGCATCCATCG	41
	F3	GCCTGCTGGATGATCTGC	18
	B3	CGGCTCTGCAACCAAAGA	18
<i>S. Typhimurium</i> invA gene (M90846.1)	FIP	TCCCGCACACGTTCTGAACCCTCTATTGTCACCGTGGTCC	40
	BIP	TGCCGATTGAAGGCCGGTACAGTACGCTTCGCCGTT	38
	F3	GCGGTGGGTTTGTGTCT	19
	B3	CGTAAAGCTGGCTTCCCTT	20
<i>V. parahaemolyticus</i> toxR gene (L11929.1)	FIP	TCGACTCCACATTCACTCGATTACTGATAACTTGCCAGACG	41
	BIP	CCTGCCGAATGGCGATTACAATTATTGTCATTAGCCG	40
	F3	CAACCATGGTGACTGTGA	18
	B3	CCAGCGACCTTCTCTGA	18
<i>L. monocytogenes</i> iap gene (JQ015300.1)	FIP	ATCAAATGTAGTTGGCCGTTACCAAGCTGAAGCTAAAAACACC	44
	BIP	ATGTATTGCTAAAGCGGAATCTGTAGTGCTAGCGTATTGTGC	44
	F3	AATTCAAGTGCAAGTGCTAT	20
	B3	GTTTGCTTGAGATTAGAGA	21



**Fig. S1.** Singleplex *Listeria monocytogenes* detection on the lab-on-a-disc. A real sample containing  $10^4$  cells of *L. monocytogenes* was used. Colour change was only observed in the reaction chamber #4, where the LAMP primer set targeting *L. monocytogenes* had been pre-stored, which confirmed specific detection of that bacteria species.



**Fig. S2.** Negative control experiments using a real sample spiked with non-pathogenic *Escherichia coli* K-12. No colour change was shown in the reaction chambers. This result proves the ability of this assay to distinguish bacteria sub-types within *E. coli* species.