


Supplementary Information for:  
**Synthesis and Photophysical Characteristics of 2,3,12,13-Tetraalkylbacteriochlorins**

Shaofei Zhang, Han-Je Kim, Quan Tang, Eunkyung Yang,  
David F. Bocian,\* Dewey Holten,\* and Jonathan S. Lindsey\*

**Table of Contents**

Section	Page
1. Acid catalysis survey of the self-condensation of a dihydrodipyririn-carboxaldehyde	S1
2. X-ray structural data for <b>BC-MM</b>	S3
3. Spectral data	S3

**1. Acid catalysis survey of the self-condensation of a dihydrodipyririn-carboxaldehyde**



**Scheme S1.** Self-condensation of dihydrodipyririn-carboxaldehyde **S1-CHO**.

**1.1. Preparation of dihydrodipyririn-carboxaldehyde S1-CHO.**

**1-Formyl-2,3-dihydro-3,3-dimethyl-7-(4-iodophenyl)dipyririn (S1-CHO).** A solution of **S1-Me**<sup>S1</sup> (50 mg, 0.13 mmol) in 1,4-dioxane (2.6 mL) was treated with SeO<sub>2</sub> (21 mg, 0.19 mmol) at room temperature. The reaction mixture was stirred at room temperature and monitored via UV-Vis spectroscopy. After completion (in this case 2 h), the reaction mixture was treated with saturated aqueous NaHCO<sub>3</sub> solution and extracted with ethyl acetate. The organic extract was washed (brine), dried (Na<sub>2</sub>SO<sub>4</sub>), concentrated, and chromatographed [silica, hexanes/ethyl acetate (3:1)] to afford a red solid (30 mg, 57%): <sup>1</sup>H NMR δ 1.23 (s, 6H), 2.73 (s, 2H), 6.30 (s, 1H), 6.33 (m, 1H), 7.01 (m, 1H), 7.18 (d, *J* = 5.1 Hz, 2H), 7.75 (d, *J* = 5.1 Hz, 2H), 10.00 (s, 1H), 10.81 (br, 1H); <sup>13</sup>C NMR δ 29.4, 41.3, 46.2, 91.9, 110.0, 112.5, 122.2, 127.3, 130.8, 136.0, 137.9, 161.3, 169.4, 190.2 (14 out of 15 distinct carbons were found); ESI-MS obsd 405.0465, calcd 405.0458 [(M + H)<sup>+</sup>, M = C<sub>18</sub>H<sub>17</sub>N<sub>2</sub>OI].

**1.2. Survey of acidic conditions.** The reactions of **S1-CHO** in the presence of different acids were examined. Each reaction was carried out (1.0 mL scale) with 10 mM dihydrodipyrin-carboxaldehyde **S1-CHO** and the indicated amount of acid in CH<sub>2</sub>Cl<sub>2</sub> at room temperature in a microreaction vial. Each reaction was monitored by UV-Vis absorption spectroscopy. Once the starting material was consumed and the absorption spectrum no longer changed, the reaction mixture was quenched by the addition of saturated aqueous NaHCO<sub>3</sub> solution. The mixture was extracted with CH<sub>2</sub>Cl<sub>2</sub>. The organic extract was washed, dried and concentrated to give the crude product. The resulting **BC-S1** and **TDC-S1** were checked via absorption spectroscopy and MALDI-MS. The relative yield of each component was determined by absorption spectroscopy.

For most reactions (Entries 2, 3-7, and 9), the absorption spectrum of each crude product showed a typical tetrahydrocorrin absorption pattern. Also, the MALDI-MS analysis of reaction mixture showed a major peak at  $m/z = 762.9$ , consistent with a tetrahydrocorrin wherein X = H (likely a product rather than precursor ion on the basis of consideration of the reaction of **S1-CHO**). If not a product ion, such a composition would be expected to tautomerize to give the tetrahydrocorrole. Regardless, the resulting putative **TDC-S1** was unstable on silica column chromatography and was not isolated; hence, the assignment remains unproven.

**Table S1.** Survey conditions for self-condensation of **S1-CHO**.<sup>a</sup>

entry	acid	[acid], (mM)	time	BC-S1	TDC-S1
1	- <sup>b</sup>	-	48 h	- <sup>c</sup>	- <sup>c</sup>
2	TsOH·H <sub>2</sub> O in CH <sub>3</sub> OH	10	16 h	<1% <sup>d</sup>	major <sup>e</sup>
3	TFA	10	1 h	- <sup>c</sup>	- <sup>c</sup>
4	ZnCl <sub>2</sub>	10	16 h	- <sup>c</sup>	major <sup>f</sup>
5	In(OTf) <sub>3</sub>	10	1 h	<1% <sup>d</sup>	major <sup>e</sup>
6	In(OTf) <sub>3</sub>	50	10 min	<1% <sup>d</sup>	major <sup>e</sup>
7	Sc(OTf) <sub>3</sub>	10	1 h	<1% <sup>d</sup>	major <sup>e</sup>
8	Zn(OTf) <sub>2</sub>	10	90 min	- <sup>c</sup>	trace
9	ZnBr <sub>2</sub>	50	30 min	<1%	major <sup>f</sup>

<sup>a</sup>All reactions were carried out (1.0 mL scale) with 10 mM **S1-CHO** and the indicated amount of acid in CH<sub>2</sub>Cl<sub>2</sub> at room temperature. Each reaction was monitored by UV-Vis absorption spectroscopy. Products were detected with MALDI-MS. <sup>b</sup>No acid was added. <sup>c</sup>Not detected. <sup>d</sup>**BC-S1** was isolated by chromatography [silica, hexanes/CH<sub>2</sub>Cl<sub>2</sub> (1:1)], characterized by UV-Vis absorption spectroscopy and MALDI-MS. <sup>e</sup>Observed by MALDI-MS and UV-Vis absorption spectroscopy but was not isolated. M = H, H. <sup>f</sup>Observed by MALDI-MS and UV-Vis absorption spectroscopy but was not be isolated. M = Zn.

**Reference.**

(S1) Balasubramanian, T.; Strachan, J.-P.; Boyle, P. D.; Lindsey, J. S. *J. Org. Chem.* **2000**, *65*, 7919–7929.



**2. X-ray structural data for BC-MM****Table S2. Summary of Crystal Data for BC-MM**

CCDC registry	
Formula	C <sub>28</sub> H <sub>34</sub> N <sub>4</sub>
Formula Weight (g/mol)	426.59
Crystal Dimensions (mm)	0.136 × 0.241 × 0.483
Crystal System	Monoclinic
Space Group	P2 <sub>1/c</sub>
Temperature, K	100(2)
<i>a</i> , Å	7.6659(13)
<i>b</i> , Å	15.982(3)
<i>c</i> , Å	9.5235(15)
α, deg	90
β, deg	92.263(8)
γ, deg	90
<i>V</i> , Å <sup>3</sup>	1165.9(3)
Number of reflections to determine final unit cell	158
Min and Max 2θ for cell determination (deg)	7.262, 59.29
<i>Z</i>	2
<i>F</i> (000)	460
ρ (g/cm)	1.215
λ, Å, (MoKα)	0.71073
μ, (mm <sup>-1</sup> )	0.072
Max 2θ for data collection (deg)	2.49, 31.20
Measured fraction of data	0.996
Number of reflections measured	25895
Unique reflections measured	3766
R <sub>merge</sub>	3.83%
Number of parameters in least-squares	153
R <sub>1</sub>	0.0436
wR <sub>2</sub>	0.1184
R <sub>1</sub> (all data)	0.0547
wR <sub>2</sub> (all data)	0.1251

### 3. Spectral data.



STANDARD 1H OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.935 sec  
Width 8000.0 Hz  
32 repetitions  
OBSERVE H1, 400.1246889 MHz  
DATA PROCESSING  
FT size 32768  
Total time 1 min, 40 sec





13C OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.139 sec  
Width 25000.0 Hz  
10000 repetitions  
OBSERVE C13, 100.6114219 MHz  
DECOUPLE H1, 400.126927 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 3 hr, 34 min, 55 sec





STANDARD 1H OBSERVE

Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-300BB "ncsumerc300"  
Relax. delay 1.000 sec  
Pulse 36.0 degrees  
Acq. time 1.995 sec  
Width 4506.5 Hz  
S2 repetitions  
OBSERVE H1, 299.7918551 MHz  
DATA PROCESSING  
F1 size 32768  
Total time 1 min, 39 sec



13C OBSERVE



Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Filter: CNRR\_DNDP-ME  
Mercury-300BB "ncsumerc300"  
Pulse 39.0 degrees  
Acq. time 1.815 sec  
Width 18761.7 Hz  
10000 repetitions  
OBSERVE C13, 75.3827175 MHz  
DECOUPLE H1, 299.7932667 MHz  
Power 40 dB  
Continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
F1 size 131072  
Total time 5 hr, 52 min, 1 sec



STANDARD 1H OBSERVE



Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient Temperature  
Mercury-400BB "ncsumerc400"

Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.335 sec  
Width 6005.0 Hz  
32 repetitions  
OBSERVE H1, 400.1245660 MHz  
DATA PROCESSING  
F1 size 32759  
Total time 1 min, 40 sec





STANDARD 1H OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 78.7 degrees  
Acq. time 1.993 sec  
Width 6000.0 Hz  
32 F2 Acquisitions  
OBSERVE H1, 400.1245662 MHz  
DATA PROCESSING  
FT size 32768  
Total time 1 min, 40 sec





13C OBSERVE



Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 46.8 degrees  
Acq. time 1.159 sec  
Width 25000.0 Hz  
10000 F2 Acquisitions  
OBSERVE C13, 100.6113782 MHz  
DECUPLE H1, 400.1266627 MHz  
Power 99 dB  
continuously on  
MAGTC-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 3 hr, 34 min, 55 sec




STANDARD 1H OBSERVE  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-300BB "ncsumerc338"  
Relax. delay 1.000 sec  
Pulse 50.4 degrees  
Acq. time 1.395 sec  
Width 6506.5 Hz  
16 repetitions  
OBSERVE H1, 300.1683398 MHz  
DATA PROCESSING  
FT size 32748  
Total time 0 min, 49 sec




13C OBSERVE  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
10000 repetitions  
OBSERVE C13, 100.6113929 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65550  
Total time 3 hr, 34 min, 55 sec




STANDARD 1H OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-300B5 "nsumerc338"  
Relax. delay 1.000 sec  
Pulse 50.4 degrees  
Acq. time 1.995 sec  
Width 4506.5 Hz  
32 repetitions  
OBSERVE H1, 300.1683390 MHz  
DATA PROCESSING  
FT size 32768  
Total time 0 min, 0 sec




10-CHO



13C OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400B5 "nsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
10000 repetitions  
OBSERVE C13, 100.6113889 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 3 hr, 34 min, 55 sec




10-CHO




STANDARD 1H OBSERVE

Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.958 sec  
Width 7867.0 Hz  
286 repetitions  
OBSERVE H1, 400.1245737 MHz  
DATA PROCESSING  
F1 size 32760  
Total time 15 min, 4 sec




HBC-MM




13C OBSERVE



Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
45804 repetitions  
OBSERVE C13, 100.6113684 MHz  
DECOUPLE H1, 400.1268027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
F1 size 65536  
Total time 17 hr, 54 min, 26 sec





HBC-MM



STANDARD 1H OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-300BB "ncsumerc300"  
Relax. delay 1.000 sec  
Pulse 26.0 degrees  
Acq. time 1.995 sec  
Width 4516.5 Hz  
32 repetitions  
OBSERVE H1, 299.7918152 MHz  
DATA PROCESSING  
F1 size 32768  
Total time 1 min, 39 sec





13C OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-300BB "ncsumerc300"  
Pulse 39.0 degrees  
Acq. time 1.815 sec  
Width 18761.7 Hz  
10000 repetitions  
OBSERVE C13, 75.3826926 MHz  
DECOUPLE H1, 299.792667 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
F1 size 131072  
Total time 5 hr, 52 min, 1 sec





STANDARD 1H OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.995 sec  
Width 6006.0 Hz  
16 repetitions  
OBSERVE H1, 400.1245671 MHz  
DATA PROCESSING  
FT size 32768  
Total time 0 min, 50 sec





13C OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.198 sec  
Width 25006.0 Hz  
1920 repetitions  
OBSERVE C13, 100.6113960 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 3 hr, 34 min, 55 sec





STANDARD 1H OBSERVE

Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient Temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.995 sec  
Width 6000.0 Hz  
32 repetitions  
OBSERVE H1, 400.1245662 MHz  
DATA PROCESSING  
F1 size 32768  
Total time 1 min, 40 sec





13C OBSERVE



Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient Temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
2304 repetitions  
OBSERVE C13, 100.6113782 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
F1 size 65536  
Total time 1 hr, 47 min, 29 sec





STANDARD 1H OBSERVE  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient Temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.985 sec  
Width 6000.0 Hz  
32 repetitions  
OBSERVE H1, 400.1245440 MHz  
DATA PROCESSING  
F1 size 32768  
Total time 1 min, 40 sec





13C OBSERVE  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient Temperature  
Mercury-400BB "ncsumerc400"  
Pulse 29.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
1600 repetitions  
OBSERVE C13, 100.6113983 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
F1 size 65536  
Total time 3 hr, 34 min, 55 sec





2,3,4-triester-5-nitroethylpyrrole  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.395 sec  
width 6006.0 Hz  
64 repetitions  
OBSERVE H1, 400.1245495 MHz  
DATA PROCESSING  
FT size 32768  
Total time 3 min, 20 sec





13C OBSERVE  
Pulse Sequence: s2pu1  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 23.0 degrees  
Acq. time 1.199 sec  
width 25000.0 Hz  
7104 repetitions  
OBSERVE C13, 100.6113782 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 7 hr, 9 min, 47 sec





STANDARD 1H OBSERVE  
 Pulse Sequence: s2pul  
 Solvent: CDCl3  
 Ambient Temperature  
 Mercury-400BB "ncsumerc400"  
 Relax. delay 1.000 sec  
 Pulse 45.0 degrees  
 Acq. time 1.995 sec  
 Width 6000.0 Hz  
 64 repetitions  
 OBSERVE H1, 400.1245645 MHz  
 DATA PROCESSING  
 FT size 32768  
 Total time 3 min, 20 sec



13C OBSERVE  
 Pulse Sequence: s2pul  
 Solvent: CDCl3  
 Ambient Temperature  
 Mercury-400BB "ncsumerc400"  
 Pulse 23.0 degrees  
 Acq. time 1.199 sec  
 Width 25000.0 Hz  
 11456 repetitions  
 OBSERVE C13, 100.6114028 MHz  
 DECOUPLE H1, 400.1266027 MHz  
 Power 40 dB  
 continuously on  
 WALTZ-16 modulated  
 DATA PROCESSING  
 Line broadening 1.0 Hz  
 FT size 65536  
 Total time 10 hr, 44 min, 38 sec



STANDARD 1H OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Relax. delay 1.000 sec  
Pulse 45.0 degrees  
Acq. time 1.495 sec  
Width 6000.0 Hz  
32 repetitions  
OBSERVE H1, 400.1245634 MHz  
DATA PROCESSING  
FT size 32768  
Total time 1 min, 40 sec



13C OBSERVE  
Pulse Sequence: s2pul  
Solvent: CDCl3  
Ambient temperature  
Mercury-400BB "ncsumerc400"  
Pulse 28.0 degrees  
Acq. time 1.199 sec  
Width 25000.0 Hz  
8832 repetitions  
OBSERVE C13, 100.6113891 MHz  
DECOUPLE H1, 400.1266027 MHz  
Power 40 dB  
continuously on  
WALTZ-16 modulated  
DATA PROCESSING  
Line broadening 1.0 Hz  
FT size 65536  
Total time 10 hr, 44 min, 38 sec

