Supporting Information

Label-free detection of fibrinogen based on fibrinogenenhanced peroxidase activity of fibrinogen-hemincomposite

Tingting Hou, Yuanfu Zhang*, Tao Wu, Meifeng Wang, Yinghong Zhang, Rui Li, Lei Wang, Qingwang Xue, Shuhao Wang *

Shandong Provincial Key Laboratory of Chemical Energy Storage and Novel Cell Technology, School of Chemistry and Chemical Engineering, Liaocheng University, Liaocheng 252000, China



Fig. S1. UV–vis spectra of hemin (black), Fib (red) and the hemin-Fib (blue). Experimental conditions: hemin 25 μ M, Fib 0.1 μ M.



Fig. S2. CD spectra of hemin (red), Fib (blue) and the hemin-Fib (black). Experimental conditions: hemin 25 μ M, Fib 0.1 μ M.



Fig. S3. TEM images of Fib (A), hemin (B), and hemin-Fib (C).



Fig. S4. SDS-polyacrylamide gel electrophoresis of Fib.



Fig. S5. Langmuir isotherms obtained from the Fib interactions with various concentrations of hemin at room temperature. (a) Non-linear regression between ΔA and [Hemin], (b) linear regression between [Hemin]/ ΔA and [Hemin]. Experimental conditions: Fib 50nM, ABTS 5 mM, H₂O₂ 2.5 mM, 20 mM Tris-HCl (pH 8.0).

Hemin + Fib \implies Hemin · Fib (1)

$$K_{b} = \frac{[\text{Hemin}\cdot\text{Fib}]}{[\text{Hemin}][\text{Fib}]}$$
(2)

$$\triangle A = \triangle A^{\max} \frac{[\text{Hemin}]K_b}{1 + [\text{Hemin}]K_b}$$
(3)

$$\frac{[\text{Hemin}]}{\triangle A} = \frac{1}{\triangle A^{\max}K_b} + \frac{[\text{Hemin}]}{\triangle A^{\max}}$$
(4)



Fig. S6. The effect of catalytic activity on hemin. Experimental conditions: 100 μ L of 5 mM ABTS, 100 μ L of 2.5 mM H₂O₂, 100 μ L of 20 mM Tris-HCl (pH 8.0), 7 min.



Fig. S7. The effect of catalytic activity on the interaction between hemin and Fib with difference pH. Experimental conditions: 100 μ L of 5 mM ABTS, 100 μ L of 2.5 mM H₂O₂, 100 μ L of 20 mM Tris-HCl (pH 8.0), 50 μ L of 25 μ M hemin or hemin-Fib (Fib 50 nM, hemin 25 μ M), 7 min.



Fig. S8. The effect of catalytic activity on ABTS. Experimental conditions: 100 μ L of 2.5 mM H₂O₂, 100 μ L of 20 mM Tris-HCl (pH 8.0), 50 μ L of 25 μ M hemin or hemin-Fib (Fib 50 nM, hemin 25 μ M), 7 min.



Fig. S9. The effect of catalytic activity on pH. Experimental conditions: 100 μ L of 5 mM ABTS, 100 μ L of 2.5 mM H₂O₂, 100 μ L of 20 mM Tris-HCl, 50 μ L of 25 μ M hemin or hemin-Fib (Fib 50 nM, hemin 25 μ M), 7 min.



Fig. S10. The effect of catalytic activity on the reaction time of ABTS- H_2O_2 system. Experimental conditions: 100 µL of 5 mM ABTS, 100 µL of 2.5 mM H_2O_2 , 100 µL of 20 mM Tris-HCl (pH 8.0), 50 µL of 25 µM hemin or hemin-Fib (Fib 50 nM, hemin 25 µM), 7 min.