

1 **Electronic Supplementary Information**

2

3

4 **Optimisation of an extraction/leaching procedure for the characterization and quantification**
5 **of titanium dioxide (TiO₂) nanoparticles in the aquatic environment using SdFFF-ICP-MS**
6 **and SEM EDX analysis.**

7

8 *Florian Dutschke^{a,b}, Johanna Irrgeher^a, Daniel Pröfrock^{a*}*

9

10

11 ^a Helmholtz-Zentrum Geesthacht, Zentrum für Material- und Küstenforschung, Department of
12 Marine Bioanalytical Chemistry, Max-Planck Str. 1, 21502 Geesthacht,
13 Germany.

14

15 ^b University of Hamburg, Institute of Applied and Inorganic Chemistry, 20146 Hamburg,
16 Germany.

17

18

19 *Corresponding Author E-mail: daniel.proefrock@hzg.de

20

21

22

23 **Used equations**

24

25 **Calculation of the Settling time (t_s) for the size cut-off using a centrifuge**

26

$$t_s = \ln\left(\frac{r_{out}}{r_{in}}\right) \frac{18 \eta}{4\pi^2 \Delta\rho d^2 \left(\frac{rpm}{60}\right)^2} \quad (\text{ESI 1})$$

27

28

29 t_s : settling time (min)

30 r_{out} : outer diameter of the centrifuge tube (cm)

31 r_{in} : inner diameter of the centrifuge tube (cm)

32 η : viscosity kg/(s·m)

33 $\Delta\rho$: delta density (water and TiO₂, kg/m³)

34 rpm: rounds per min

35 d: diameter of the particles (m)

36

37

38 **Calculation of the particle size distribution**

39

40 The following equations have been used to perform different calculations related with the analysis of nano
41 particles via SdFFF. The SdFFF theoretical basis for the conversion of retention time data into size
42 distribution is described in detail elsewhere.^{1,2} Briefly, the retention ratio R obtained experimentally by
43 SdFFF is a function of the layer thickness λ :

44

$$R = \frac{t_0}{t_r} = 6\lambda \left[\coth\left(\frac{1}{2\lambda}\right) - 2\lambda \right] \quad (\text{ESI 2})$$

45

46

47 Being t_0 the void time and t_r the retention time.

48

49 In SdFFF, the layer thickness λ is expressed as:

50

$$\lambda = \frac{kT}{m\omega^2rw} \quad (\text{ESI 3})$$

51

52

53 k the Boltzmann constant

54 T the absolute temperature in K

55 w the channel thickness (m)

56 ω the angular velocity

57 r the radius of the centrifuge from rotation axis to the channel (m).

58 m mass of the particle

59 t_r retention time (min)

60

61 The particle size can be calculated from the retention time obtaining a size distribution profile by applying
62 the equation:

63

$$d = \left[\left(\frac{6kT}{\pi G w \Delta \rho t_0} \right) t_r \right]^{1/3} \quad (\text{ESI 4})$$

64

65

66 d the diameter of the particles according to FFF theory

67 k the Boltzmann constant

68 T the absolute temperature in K

69 w the channel thickness (m)

70 t_0 void time (min)

71 r the radius of the centrifuge from rotation axis to the channel (m)

72 t_r retention time (min)

73 G the centrifugal acceleration

74 $\Delta \rho$ the difference between the density of the particle and the one of the carrier liquid.

75

76

77 **Tables**

78

79 **ESI Table 1** Detailed information related with the different sampling sites

Sample name	Sampling date	Air- / Water temperature (°C)	pH-Value	Conductivity (µS/cm)
December 2016 (west)	20.12.2016	3 / 5.2	7.83	673
December 2016 (south)	20.12.2016	2 / 5.3	7.68	699
September 2016 (west)	23.09.2016	17 / 19.6	8.15	690
September 2016 (south)	23.09.2016	12 / 19.6	8.06	689
June 2016 (west)	01.06.2016	24 / 19.2	8.08	684
June 2016 (south)	01.06.2016	21 / 19.1	7.98	685
June 2015 (west)	01.06.2015	22 / 18.5	8.14	674

June 2015 (south)	01.06.2015	20 / 18.5	8.09	677
-------------------	------------	-----------	------	-----

80

- 81 **ESI Table 2** Instrumental parameters used for the operation of SdFFF (CF2000, PostNova Analytics) and
 82 ICP-MS/MS (Agilent 8800, Agilent Technologies) coupling

SdFFF

Tip to tip channel length	89.5 cm
Channel breadth	2.1 cm
Channel Thickness (<i>w</i>)	0.0231 cm
Injection volume	21.54 μ L
Void volume (V^0)	2.0 mL
Relaxation time	5 min
Channel flow rate	0.5 mL/min
Initial rotation speed	2500 rpm
Power field programming	$t_1 = 8$ min (focus time at the initial field), and $t_a = -64$ (field decay parameter)
Final field	100 rpm
Carrier solution	0.15 mmol/L NaOH

ICP-MS/MS

RF-Power	1600 W
Torch/Injector	Quartz/Quartz
Spray chamber	Quartz
Sampler and skimmer cones	Cu-Ni, Ni
Carrier gas flow	1.12 L/min
Makeup gas flow	0.11 L/min
Cell gas flow rate (He)	4.5 mL/min
Q1 bias	-2 V
Qtopole bias	-5 V
Dwell times	0.3 s
Tuning solution	10 μ g/L (Li, Y, Tl, Ce in 2 % (w/w) HNO ₃)

83

- 84 **ESI Table 3** Details of the used microwave assisted acid digestion power, temperature and time program

Energy / Temperature	Time
1600 W / 180 °C	60 min
1600 W / 180 °C	180 min
1600 W / 180 °C	60 min

85

86

87 **ESI Table 4** Instrumental parameter used during the measurement of the sediment digests via ICP-
88 MS/MS (Agilent 8800, Agilent Technologies)


ICP-MS/MS

RF-Power	1550 W
Torch/Injector	Quartz/Sapphire
Spray chamber	PFA
Sampler and skimmer cones	Cu-Ni, Ni
Carrier gas flow	1.12 L/min
Makeup gas flow	0.11 L/min
Cell gas flow rate (H ₂ ,He,O ₂)	6 mL/min, 4.5 mL/min, 2.5 mL/min (30 %)
Q1 bias	-2 V
Qctopole bias	-5 V
Dwell times	0.3 s
Tuning solution	10 µg/L (Li, Y, Tl, Ce in 2 % (w/w) HNO ₃)

89
90
91
92
93
94
95
96
97
98
99
100
101
102
103
104
105
106
107
108
109
110
111

112 **Figures**

113




114

115 **ESI Figure 1** Decay program used for the separation of the particles with the SdFFF system

116

117




118

119 **ESI Figure 2** Comparison of the intensity data from the ICP-MS measurement and the size-information
120 obtained based on mathematical calculations using the MALS (dots) signal for the reference material
121 NIST SRM 1898 before (black) and after (grey) the application of the developed sample pretreatment
122 procedure.

123

124




125

126 **ESI Figure 3** Fractograms and size-distribution profiles after spiking and extraction of a natural sediment
 127 test sample to compare the effects of different sample preparation procedures (freeze drying vs. direct
 128 preparation of the wet fresh sample) including data on the measured gyration diameter dg .

129


130



131

132 **ESI Figure 4** Different fractograms obtained for test mixtures of NIST SRM 2702 spiked with 100 mg/kg
 133 NIST SRM 1898 using different extraction solvents (grey cross (water), grey circle (0.05% Tego W752),
 134 grey box (0.05% FI-70))

135



136

137 **ESI Figure 5** Compilation of different fractograms showing the particle size-distribution of the
 138 investigated lake sediment samples. (A) December 2016 south, (B) December 2016 west, (C) September
 139 2016 south, (D) September 2016 west, (E) June 2016 south, (F) June 2016 west, (G) June 2015 south,
 140 (H) June 2015 west.