

Supplementary Information

Scaling the effect of hydrophobic chain length variations on gene transfer properties of di-alkyl, di-hydroxy ethylammonium chloride based cationic amphiphiles

Ankita A. Hiwale^{1†}, Chandrashekhara Voshavar^{2†}, Priya Dharmalingam³, Ashish Dhayani¹, Rajesh Mukthavaram⁴, Rasajna Nadella³, Omprakash Sunnapu⁵, Sivaraman Gandhi¹, VGM Naidu⁵, Arabinda Chaudhuri⁶, Srujan Marepally^{1,3*}, Praveen Kumar Vemula^{1*}

1. Institute for Stem Cell Biology and Regenerative Medicine (inStem), GKVK-post, Bellary Road, Bengaluru 560065, India.

2. BioSatva Technologies, Golnaka, Hyderabad 500013, India.

3. Centre for Stem Cell Research, Christian Medical College Campus, Bagayam, Vellore 632002, India.

4. Translational Neuro-oncology Laboratories, Moores Cancer Center, University of California San Diego, La Jolla, CA, USA.

5. National Institute for Pharmaceutical Education and Research, Balanagar, Hyderabad 500018, India.

6. Biomaterials Group, CSIR-Indian Institute of Chemical Technology, Hyderabad 500 007, India.

† Equally contributing authors

* Corresponding authors: Srujankm@cmcvellore.ac.in, Praveenv@instem.res.in

Contents	Page number
DLS data (Size & Zeta Potentials) for Liposome (Figure S1)	S3
DLS data (Size & Zeta Potentials) for Lipoplexes (Figure S2 and S3)	S4-S5
Cytotoxicity of lipids in HeLa cells (Figure S4)	S6
Transfection efficiency of lipids in CHO and B16F10 cells (Figure S5)	S7
SAXS diagrams of lipids (Figure S6)	S8
Characterization of lipids	S9-S23

Figure S1. Hydrodynamic diameters (top) and surface potentials (bottom) of liposomes of lipids, Lip1818-Lip1810





Figure S2. Sizes of the lipid-DNA complexes (lipoplexes) measured using dynamic light scattering.





Figure S3. Surface charge of the lipid-DNA complexes (lipoplexes) measured using dynamic light scattering.

Zeta potentials






Figure S4. Cytotoxicity evaluation of Lipids (Lip1818-Lip1810) in HeLa cells for 24h incubation. HeLa cells were incubated with Lipids (Lip1818-Lip1810) for 24h at varying Lipid:DNA charge ratios along with Lipofectamine 3K. Control cells received DMEM complete media with 0.01% DMSO. MTT assay was performed as per the protocol described in Materials & Methods. Values shown are means \pm SEM of three independent experiments performed in six to eight replicates.




Figure S5. In vitro transfection efficiencies of Lipids (Lip1818-Lip1810) in CHO and B16F10 cell lines at lipid:DNA charge ratios of 8:1-1:1. Absorption at 405 nm was converted to β -galactosidase units using a calibration curve constructed with pure commercial β -galactosidase enzyme. Units of β -galactosidase activity were plotted against the varying lipid-to-DNA charge ratios (8:1-0.5:1). The transfection values shown are the average of triplicate experiments performed with 4-6 data points.






Figure S6. SAXS diagrams of Lipids (Lip1818-Lip1810). Radially integrated scattering intensities of the lipids, with and without the Model membrane (MM) at 1:1 molar ratio. The Model membrane is a phospholipid mixture of DOPC:DOPE:DOPS:Chol at 45:20:20:15, w/w content that closely mimics mammalian cell membrane. The symbols correspond to the measured data and the solid curves represent the smoothed curve with weighted average 9 neighbours to the data. Curve smoothing was done using GraphPad Prism software.


1H NMR of Lip1818




¹³C NMR of Lip1818




ESI-Mass of Lip1818




1H NMR of Lip1816




¹³C NMR of Lip1816




ESI-Mass of Lip1816




1H NMR of Lip1814




13C NMR of Lip1814




ESI-Mass of Lip1814




1H NMR of Lip1812




13C NMR of Lip1812




ESI-Mass of Lip1812



1H NMR of Lip1810




13C NMR of Lip1810




HPLC Profiles of lipids, Lip1818-Lip1810


Lip1818




Lip1816




Lip1814



Lip1812



Lip1810



HPLC Conditions:

System:	Agilent 1100 Series HPLC Value System
Column:	LiChrospher® 100A, RP-18e (5µm), 125 x 4 mm
Mobile Phase:	100% Methanol.
Flow Rate:	2.0 mL/min
Column Pressure:	60-65 Bars
Detection:	Evaporative Light Scattering Detector (ELSD)