Electronic Supplementary Material (ESI) for Analytical Methods. This journal is © The Royal Society of Chemistry 2018

Supplementary material

Deep eutectic solvent-based extraction coupled with green two-dimensional HPLC-DAD-ESI-MS/MS for the determination of anthocyanins from *Lycium ruthenicum* Murr. fruit

Jun Sang, Bing Li, Ya-ya Huang, Qun Ma, Kang Liu, Cui-qin Li\*

Key Laboratory of the Ministry of Education for Medicinal Resources and Natural Pharmaceutical Chemistry, National Engineering Laboratory for Resource Development of Endangered Crude Drugs in Northwest of China, College of Life Sciences, Shaanxi Normal University, Xi'an, Shaanxi 710062, China

E-mail address: licuiqin16@snnu.edu.cn (C.-Q. Li), junsang@foxmail.com (J. Sang)

<sup>\*</sup> Corresponding author. Tel.: +86 29 85310266; fax: +86 29 85310546.




Fig. S1 Effects of liquid/solid ratio (A) and extraction time (B) on extraction yield of total anthocyanins from L. ruthenicum fruit. Means with different letters were significantly different at the level of p < 0.05.




Fig. S2 HPLC-DAD chromatograms of anthocyanins from *L. ruthenicum* fruit using ethanol and fruit acid aqueous solutions (0.1 mol/L) at a ratio of 15:85 (v/v) as mobile phases, including mandelic (A), ascorbic (B), lactic (C), glycolic (D), malic (E), citric (F), and tartaric (G) acids; peak 1, petunidin-3-*O-(trans-p-*coumaroyl)-rutinoside-5-*O-*glucoside.

Supplementary Table 1
Physical characteristics of studied macroporous resins.

Resins	Polarity	Surface area (m <sup>2</sup> /g)	Pore size (Å)
NKA-9	Polar	250-290	155-165
AB-8	Weak-polar	480-520	130-140
XDA-6	Weak-polar	$\geq 650$	125
XDA-7	Weak-polar	$\geq 800$	45
HPD-80	Non-polar	350-400	80-85
D101	Non-polar	400-600	100-110
HPD-100A	Non-polar	600-700	90-100
X-5	Non-polar	500-600	290-300
LS-305	Non-polar	$\geq 800$	230
LX-68	Non-polar	$\geq$ 980	45-50
LX-32	Non-polar	$\geq 1000$	_ a
LX-31B	Non-polar	$\geq 1200$	_ a

<sup>&</sup>lt;sup>a</sup> Date are unpublished by resin manufacturer.