

Microfluidic On-demand Engineering of Exosomes towards Cancer Immunotherapy

Zheng Zhao^a, Jodi McGill^b, Pamela Gamero-Kubota^c, and Mei He^{a,d*}

^a. Department of Chemical and Petroleum Engineering, Bioengineering Program, University of Kansas, Lawrence, Kansas, USA

^b. Department of Veterinary Microbiology and Preventive Medicine, Iowa State University, Iowa, USA

^c. Department of Molecular Biosciences, University of Kansas, Lawrence, Kansas, USA

^d. Department of Chemistry, University of Kansas, Lawrence, Kansas, USA

Supplemental Information




Figure s1. Illustration of 3D printing approach for one-step producing 3D mold and replicating PDMS microfluidic device integrated with cell culture and downstream exosome isolation, surface engineering, and on-demand photo release.




Figure s2. Investigation of the side-effect of UV exposure on exosome molecular contents in terms of proteins, DNAs and RNAs.





Figure S3. The fluorescence intensity analysis for showing the cellular uptake rate of gp-100 engineered exosomes and native exosomes.



Negative Control PWM Stimulator Engineered Exosomes

Figure s4. Human leukocytes culture under different stimulation conditions: 1) negative control is the leukocytes without any stimulation; 2) PWM protein stimulation as the positive control; 3) The gp-100 engineered exosome stimulation.




Figure s5. ex vivo testing of surface-engineered exosomes for activating transgenic mice spleen-derived CD8+ T cells. a) depicts representative flow plots from wells containing T cells + Activated JAWS cells with increasing concentrations of the gp100-engineered exosomes. The spiking gp-100 in μM serves as positive control. b) depicts the cumulative data from all four culture conditions showing the CD8+ T cell dividing rate under stimulation. The results are representative of 3 independent experiments with three duplicate wells for each culture condition (RSD < $\sim 5\%$).


Tumor peptide synthesis and characterization

The protocols follow standard Fmoc chemistry. The peptides were cleaved using a solution of 92.5:2.5:2.5:2.5 TFA:TIPS:H₂O:DODt and the crude peptides were purified using preparative HPLC (gradients of water/ acetonitrile (90:10 to 0:100 containing 0.1% TFA over 40 min) and lyophilized to obtain white powder. Analytical HPLC traces were acquired using an Agilent 1100 quaternary pump and a Hamilton PRP-1 (polystyrene-divinylbenzene) reverse phase analytical column (7 μm particle size, 4 mm x 25 cm) with UV detection at 210 nm. The eluents were heated to 45 $^{\circ}\text{C}$ to reduce separation of rotational isomers, and elution was achieved with gradients of water/ acetonitrile (90:10 to 0:100 containing 0.1% TFA) over 20 min. Low-resolution mass spectra (LRMS) were obtained using a Waters Micromass ZQ 4000 instrument with ESI+ ionization.


Peptide gp-100 Sequence: RLMKQDFSV

Chemical Formula: C₄₉H₈₂N₁₄O₁₄S₁

Molecular Weight: 1123.33



Analytical HPLC profile of synthesized peptide gp-100. Retention time = 7.51min (monitored at 210 nm). Purity > 90% by HPLC.




Low-resolution mass spectrum of peptide synthesized gp-100


Peptide MART-1 Sequence: ELAGIGILTV

Chemical Formula: C₄₅H₈₀N₁₀O₁₄

Molecular Weight: 985.18



Analytical HPLC profile peptide MART-1. Retention time = 9.80 min (monitored at 210 nm). Purity > 90% by HPLC



Low-resolution mass spectrum of peptide MART-1