Electronic Supplementary Material (ESI) for New Journal of Chemistry. This journal is © The Royal Society of Chemistry and the Centre National de la Recherche Scientifique 2018

Supplementary Information

Selective recognition of atropine in biological fluids and leaves of Datura


stramonium employing carbon nanotubes – chitosan film based biosensor

Suyash Mane, Rushda Narmawala and Sanghamitra Chatterjee*

Department of Chemistry, Institute of Chemical Technology, Matunga, Mumbai 400019, India

^{*}Corresponding author. Tel.: +91-22-33611144

Email address: sk.chatterjee@ictmumbai.edu.in (S. Chatterjee)



Scheme S-1 Synthetic route of the developed biosensor.




Fig. S-1 Voltammetric response obtained for 100 μM atropine in PBS of pH 10.0 with (A) (a) 0.25, (b) 0.5; (B) (c) 1.0 and (d) 2.0 mg mL⁻¹ suspension of SWCNTs in 4% chitosan coated on GCE.




Fig. S-2 Linear calibration plot of atropine in the range **(A)** 100 to 1000 nM and **(B)** 5 to 150 μM observed at SWCNT/Chit/GCE.




Fig. S-3 Effect of interferents of 500 μ M each of glucose, uric acid, ascorbic acid, paracetamol, xanthine and hypoxanthine on 5 μ M of atropine at SWCNT/Chit/GCE.




Fig. S-4 (A) Stability and (B) Sensor reproducibility of the SWCNT/Chit/GCE for 50 μM atropine in PBS of pH 10.0. (C) Intra day and (D) Inter day reproducibility of the SWCNT/Chit/GCE for 50 μM atropine in PBS of pH 10.0.




Fig. S-5 (A) Voltammetric response obtained at SWCNT/Chit/GCE for leaves of *Datura* stramonium sample 1 (a) before spiking and after (b) 15, (c) 30 and (d) 45 μM spiking of atropine.

(B) Voltammetric response obtained at SWCNT/Chit/GCE for pharmaceutical sample

1 (a) before spiking and after (b) 10, (c) 20 and (d) 30 µM spiking of atropine.

	Spiked (µM)	Detected (µM)	Recovery (%)
Human Urine			
Sample 1	0.00		
	20.00	19.63	98.15
	40.00	39.06	97.65
	60.00	59.12	98.53
Sample 2	0.00		
	20.00	19.31	96.55
	40.00	40.75	101.88
	60.00	57.87	96.45
Human Serum			
Sample 1	0.00		
	30.00	30.60	102.00
	60.00	58.81	98.02
	90.00	87.02	96.69
Sample 2	0.00		
	30.00	29.97	99.90
	60.00	60.56	100.93
	90.00	87.65	97.39

Table S-1 Concentration of atropine in human urine and human serum samples atSWCNT/Chit/GCE at pH 10.0.

Technique	Sensor Used	Analytical Range (M)	Limit of Detection (M)	Sensitivity (µA/µM)	Reference
DPV	Co ₃ O ₄ -reduced graphene oxide modified carbon paste electrode	0.1 × 10 ⁻⁶ – 3.2 × 10 ⁻⁶	0.03 × 10 ⁻⁶	2.59	28
LSV	Poly (3, 4-ethylene- dioxythiophene) film modified Pt electrode	0.1 × 10 ⁻⁶ – 100 × 10 ⁻⁶	27 × 10 ⁻⁹		31
CV	Screen printed graphite electrode	5 × 10 ⁻⁶ – 50 × 10 ⁻⁶	3.9 × 10 ⁻⁶	0.01	13
DPV	Modified pencil graphite electrode	$\begin{array}{c} 0.6 \times 10^{\text{-6}} - \\ 600 \times 10^{\text{-6}} \end{array}$	$30 imes 10^{-9}$	0.20	29
SWV	SWCNT/Chit/GCE	100 × 10 ⁻⁹ - 150 × 10 ⁻⁶	16.5 × 10 ⁻⁹	2.70	This work

Table	S-2	Comparison	of	the	SWCNT/Chit/GCE	with	the	reported	sensors	for	the
electrochemical determination of atropine.											