


Smart Electrochemical Sensing Platform for the Simultaneous Determination of Psychotic Disorders Drugs Isopropamide Iodide and Trifluoperazine Hydrochloride

Mona A. Mohamed ^{a,b*}, Ahmed S. Saad ^{c*}, Sara H. Koshek ^a and Mohamed R. El-Ghobashy ^c

^a Pharmaceutical Chemistry Department, National Organization for Drug Control and Research, Pyramids Ave, P.O. Box 29, Giza, Egypt

^b Institute of Electronics, Microelectronics and Nanotechnology (IEMN, UMR CNRS 8520), Université Lille 1, Cité Scientifique Avenue Poincaré-BP60069, 59652 Villeneuve d'Ascq, France

^c Cairo University, Faculty of Pharmacy, Analytical Chemistry Department, Kasr Elaini St., P.O. box 11562, Cairo, Egypt.



Scheme S1. Suggested electrochemical oxidation of ISP and TFP




Figure S1. Nyquist plots for MC-CPE, MWCNT/CPE, and MC-MWCNT-CPE using the redox probe 5.0 mmol L⁻¹ K₃Fe(CN)₆ (1:1) solution.




Figure S2. Chronoamperograms for the oxidation of different concentrations ISP at MCP-MWCNT/CPE in B-R buffer (pH 4.0), for a potential step of +0.51 V vs. Ag/AgCl. The numbers 1 to 4 in Cottrell's plot (inset *i*) correspond to 62.50, 100.00, 117.60, and 134.60 $\mu\text{mol L}^{-1}$ of ISP, respectively. Inset *ii* shows the variation of chronoamperometric currents at $t = 30$ s vs. ISP concentration.




Figure S3. Chronoamperograms for the oxidation of different concentrations TFP at MCP-MWCNT/CPE in B–R buffer (pH 4.0), for a potential step of +0.794 V vs. Ag/AgCl. The numbers 1 to 4 in Cottrell's plot (inset *i*) correspond to 81.63, 100.00, 117.60, and 134.60 $\mu\text{mol L}^{-1}$ of TFP, respectively. Inset *ii* shows the variation of chronoamperometric currents at $t = 30$ s vs. TFP concentration.

Table S1. Determination of ISP and TFP in human plasma sample using the new proposed sensing protocol.

| Parameters | Conc. Added ($\mu\text{mol L}^{-1}$) | Conc. found ($\mu\text{mol L}^{-1}$) | %Recovery* |
|-------------------|--|--|-------------------|
| | 12.50 | 12.36 | 98.90 |
| <u>ISP</u> | 15.50 | 15.75 | 101.61 |
| | 20.90 | 20.72 | 99.13 |
| | 23.40 | 23.47 | 100.32 |
| %Mean \pm %RSD | | | 99.99 \pm 1.24 |
| | 62.24 | 63.30 | 101.75 |
| <u>TFP</u> | 77.50 | 75.80 | 97.75 |
| | 105.00 | 105.80 | 100.73 |
| | 117.00 | 116.90 | 99.90 |
| %Mean \pm %RSD | | | 100.03 \pm 1.70 |

*Average of three determinations.