

Supporting Information

Dual-emissive gold nanoclusters for separation-free and label-free ratiometric fluorescent sensing of 4-nitrophenol based on inner filter effect

Tong Shu,^a Jianxing Wang,^a Xiangfang Lin,^a Ziping Zhou,^a Feng Liang,^b Lei Su*,^a and Xueji Zhang*,^a

^a Beijing Key Laboratory for Bioengineering and Sensing Technology, Research Center for Bioengineering and Sensing Technology, School of Chemistry and Biological Engineering, University of Science and Technology Beijing, Beijing 100083, P. R. China.

^b The State Key Laboratory of Refractories and Metallurgy, Wuhan University of Science and Technology, Wuhan 430081, P.R. China.

*Email - sulei@iccas.ac.cn




Figure S1. a) TEM image of the GNCs@BSA taken on a JEOL JEM-2010 microscope. b) The corresponding histogram of size distribution of the AuNCs@BSA.




Figure S2. PL emission (solid line) and excitation (dash line) spectra of the 4-fold diluted GNCs@BSA stock solution.




Figure S3. PL emission spectra of the 4-fold (red) and 10-fold (black) diluted GNCs@BSA stock solution excited at 330 nm.




Figure S4. Photostability of the 410-nm emission of the GNCs@BSA in tri-HCl buffer solution at pH 8.0.




Figure S5. Bar plot of varied pH values on the 410-nm emission intensity of the GNCs@BSA.




Figure S6. UV-vis absorption spectra of 3.3 μM 4-NP at varied pH values (left arrow: from up to bottom, pH = 4, 5, 6, 7, 8, 9, 10, 11, 12, respectively).




Figure S7. The molecular structures of of phenol, pAP, pIP, DA, BA, HQ, pCP and 2-NP.

Table S1. Spiked recoveries and RSDs (%; n=3) for the determination of 4-NP in tap water and local river water samples using the ratiometric sensor.

Samples	Spiked (μM)	Found (μM , average value)	Recovery ^a (%)	RSD (%)
Tap water	0	ND ^b	/	/
	0.680	0.682	100.2	1.3
River water	0	ND	/	/
	0.680	0.676	99.4	1.6

^a Average value from three individual experiments.

^b ND: non detected.