Supplementary material

Quadruple-labeling luminescence strategy for multiplexed immunoassay of 51

drugs in milk with automated pretreatment system

Yaqing Zhang^{1a}, Xiaoxi Chang,^{1a} Xin Wang,¹ Xiaoqi Tao^{1, 2*}

¹College of Food Science, Southwest University, Chongqing, 400715, China

² Key Laboratory of Luminescence and Real-Time Analytical Chemistry, Ministry of

Education, Chongqing, 400715, China

The preparation of immunomagenetic beads

The preparation of immunomagnetic beads was according to the manufacture's product specifications. Briefly, streptavidin magnetic beads were added to centrifuge tubes, 0.02 M PB was added to and made the concentration of streptavidin magnetic beads to be 1 mg·mL⁻¹ (1 mL). The biotinylated mAb (30.0 μ g) was putted into centrifuge tubes and incubated for 30 min at room temperature, the immunomagnetic beads were formed and separated from the solution by magnetic grate, and washed for 4 times with 0.01 M PBS (containing 0.1% BSA). The immunomagnetic beads were resuspended with 0.1 mL of 0.02 M PB and placed at 4 °C before use.

^{*}Address correspondence to Xiaoqi Tao, College of Food Science, Southwest University, NO. 2 Tiansheng Road, Beibei District, Chongqing, 400715, China. E-mail: <u>77179000@qq.com</u>

^a These authors contributed equally to this work and should be considered co-first authors.



Supplementary Figure 1 Gross standard inhibition curves and revised standard inhibition curves in two detection steps



Supplementary Figure 2 Inhibition curves of NOR, AMP, SMX and CAP in buffer and milk extract

Compoud		LOD	MRL/MPL (ng·L ⁻¹)		
		$\mu g \cdot L^{-1}$	EU ^{20, 21}	USA ²²	China ³
FQs	CIP, ENR, NOR, PEF	≪0.04			
	DAN, LEV, OFL, AMI, ENO, ORB, FLER, SPA, MAR, LOM	< 0.08	10.0-100.0	10.0-100.0	30.0-100.0
	FLU, PAZ	< 0.20			
	PRU, SAR, DIF	<1.60			
	TRO	<4.00			
β-lactams	PEN, AMP, AMO, NAF, CEF, CEQ	≤0.06	4.0-125.0	5.0-10.0	
	DIC, OXA, CEM, CEZ, CEP, CET	< 0.24			4.0-100.0
	CLO, CER	<1.20			
	CEL	<6.00			100.0
SAs	SMX, SDM, SDT, SIM, SMD, SMM, SMP, SPY, SNT	≤0.06			25.0-100 0
	SCP, SMR, SMZ, acetyl-SMZ, SQX	< 0.30	25.0-100.0 10.0-100.0		
	SDZ	<1.20			
CAP	САР	0.006	-a (MRPL	_	_
			$0.3 \ \mu g \cdot L^{-1})$		

Supplementary Table 1 Limits of detection for FQs, β -lactams, SAs and CAP, and MRL (MRPL) set by EU, USA and China in milk

^a No level of chloramphenicol is acceptable