

Supporting Information

Immunosensing of prostate cancer in human plasma samples using immobilization of antibody on the surface of mesoporous silica modified silver nanoparticles and its immunocomplex with prostate specific antigen

Fatemeh Farshchi ^a, Mohammad Hasanzadeh ^{a*}, Elham Solhi ^b

^a Pharmaceutical Analysis Research Center, Tabriz University of Medical Sciences, Tabriz, Iran.

^b Biotechnology Research Center, Tabriz University of Medical Sciences, Tabriz, Iran, Tabriz 51664, Iran.

Corresponding Author

E-mail address: (*) hasanzadehm@tbzmed.ac.ir

Tabriz University of Medical Sciences, Tabriz.

Tel: +98(41) 33363311; Fax: +98(41)33363231




Fig S1. Chronoamperogram of Au electrode in the solution of (4.4 mM TEOS+0.1 M NaNO_3 +1.4 mM CTAB). Time: 100 s and $E = -1.3\text{V}$.




Fig. S2. CVs of Au electrode modified by SiO₂-NH₂ on the solution of AgNO₃ (0.1 M) + HNO₃ (0.1 M). Scan rate is 100 mV/s, number of cycle = 30 cycle.




Fig. S3. A) FE-SEM imaging of SiO₂-Ag NPs/Au in different magnification. Inset is FESEM image of bare Au electrode.




Fig. S3B) FE-SEM imaging of SiO₂- AgNPs /Ab/Au in different magnification.








Fig. S3C) FE-SEM imaging of SiO₂-AgNPs /Ab/Ag/Au in different magnification.








Fig. S5. A) DPVs of the engineered immunosensor for the detection of PSA biomarker in the concentration range of 0.03-0.001 $\mu g L^{-1}$ in human plasma samples. **B)** Calibration curve.




Fig. S6. A) CVs of the immunosensor in different cycle number (2, 5, 10, and 50 cycles). **B)** Histogram of the peak currents *versus* cycle number. Scan rate: 100 mV/s.





Fig. S7. DPVs (**A**) and histograms (**B**) for the stability study of $\text{SiO}_2\text{-AgNPs /Ab/BSA/Ag/Au}$ in the potential range of 0 to 0.8 and sweep 50 mV/s in 1mM of $\text{Fe (CN)}_6^{-3/-4}$ periodic time (24,48,72, and 96 h).




Fig. S8. A) DPVs SiO₂-AgNPs/Ab/BSA/Ag/Au in the presence of serum proteins (BSA, CA15.3, CA125). (B and C) Histogram of peak potential and currents *versus* type of interfering species.