

Supplementary Information

A sample-to-answer, portable platform for rapid detection of pathogens with a smartphone interface

Yu-Dong Ma¹, Kuang-Hsien Li², Yi-Hong Chen¹, Yung-Mao Lee², Shang-Ta Chou¹,

Yue-Yuan Lai¹, Po-Chiun Huang², Hsi-Pin Ma^{2*} and Gwo-Bin Lee^{1, 3, 4*}

¹Department of Power Mechanical Engineering, National Tsing Hua University,
Hsinchu, Taiwan

²Department of Electrical Engineering, National Tsing Hua University, Hsinchu,
Taiwan

³Institute of NanoEngineering and Microsystems, National Tsing Hua University,
Hsinchu, Taiwan

⁴Institute of Biomedical Engineering, National Tsing Hua University, Hsinchu,
Taiwan

Table S1. The sequences of the LAMP primers and aptamer for detection of H1N1 viruses and MRSA

Primer name	Sequence (‘5→ 3’)	Length (bp)
Inf A- F3	TCAGACAATGGAACGTGTT	19
Inf A- B3	CTTCCCTTTATCATTAATGTAGGA	25
Inf A- FIP	GTCTTGGGGAATATCTCAAACCTTTATTATGAGG AGCTAAGAGAGC	46
Inf A- BIP	GTGTAACGGCAGCATGTCCTGAATTTCTTTTTT AACTAGCCA	43
Inf A- LF	TGATGACACTGAGCTCAATT	20
Inf A- LB-1	CATGCTGGAGCAAAAAGC	18
Inf A- LB-2	GCTGGAGCAAAAAGCTTCT	19
Inf A- LB-3	CTCATGCTGGAGCAAAAAG	19
<i>mecA</i> -F3	TGATGCTAAAGTTCAAAAAGAGT	22
<i>mecA</i> -B3	GTAATCTGGAACCTTCTTGAGC	21
<i>mecA</i> -FIP	TGAAGGTGTGCTTACAAGTGCTAATCAACATGA AAAATGATTATGGCTC	49
<i>mecA</i> -BIP	TGACGTCTATCCATTTATGTATGGCAGGTTCTTT TTTATCTTCGGTTA	48
<i>mecA</i> -LF-1	AGGGTGGATAGCAGTAC	17
<i>mecA</i> -LF-2	GAGGGTGGATAGCAGTAC	18
<i>mecA</i> -LF-3	TCACCTGTTTGAGGGTGGA	19
<i>mecA</i> -LB	ATGAGTAACGAAGAATATAATAAAT	25
Inf A-aptamer	GGCAGGAAGACAAACAGCCAGCGTGACAGCGA CGCGTAGGGACCGGCATCCGCGGGTGGTCTGTG GTGCTGT	72




Figure S1. The effect of betaine on the real-time LAMP reaction for detection of MRSA.




Figure S2. The fabrication process of hydrophilic pattern slide. (a) Spin-coating of PDMS. (b) Surface activation by oxygen plasma. (c) Patterning. (d) Spin-coating of a hydrophilic UV glue film. (e) Nitrogen gas purging. (e) UV irradiation.




Figure S3. Schematic illustration of mobile APP interface.




Figure S4. The heating profile of the temperature control module, which demonstrated high precision across a range of temperatures (<1.5°C fluctuation).




Figure S5. The specificity of the colorimetric LAMP assay for detection of (a) influenza A (H1N1) virus and (b) MRSA. (c) Agarose gel electrophoresis analysis of the colorimetric LAMP products from (a). (d) Agarose gel electrophoresis analysis of the colorimetric LAMP products from (b). L=50-bp DNA ladder, N=negative control (distilled water).