

Supporting Information

A FRET-based ratiometric fluorescent probe for highly selective detection of cysteine based on a coumarin-rhodol derivative

Yu Bai^a, Ming-Xia Wu^a, Qiu-Juan Ma^{a,*}, Chun-Yan Wang^a, Jing-Guo Sun^a, Mei-Ju Tian^a,
Jian-Sheng Li^{b,*}

^a *School of Pharmacy, Henan University of Chinese Medicine, Zhengzhou 450046, PR China*

^b *Collaborative Innovation Center for Respiratory Disease Diagnosis and Treatment &*

Chinese Medicine Development of Henan Province, Henan University of Chinese Medicine,

Zhengzhou 450046, PR China

*Corresponding author, E-mail: maqiujuan104@126.com (Q. J. Ma); mxwu711@163.com (M.

X. Wu); Tel: +86-371-65676656; Fax: +86-371-65680028.

Table of contents

Figure S1	S3
Figure S2	S4
Figure S3	S5
Figure S4	S6
NMR and MS data for compounds	S7




Figure S1. The effect of CH₃CN volume fraction on fluorescence intensity ratio ($I_{543\text{nm}}/I_{470\text{nm}}$) of probe **1** (10 μM) in the absence (filled circles) and presence of 100 μM cysteine (clear circles).




Figure S2. Fluorescence emission spectra of compound **5** (10.0 μM) and the reaction mixture of fluorescent probe **1** (10.0 μM) with cysteine (100 μM) in 0.01 M PBS buffer ($\text{CH}_3\text{CN}/\text{water} = 6:4$, V/V, $\text{pH}=7.40$). The solid line (—) and medium dash (— — —) represent compound **5** and the reaction mixture of fluorescent probe **1** with cysteine, respectively. $\lambda_{\text{ex}} = 418$ nm.




Figure S3. The ESI-MS data of reaction mixture of probe 1 with Cys for 10 min




Figure S4. Time course of fluorescence intensity ratio ($I_{453\text{nm}}/I_{470\text{nm}}$) of probe **1** ($10\ \mu\text{M}$) in the absence (filled circles) and presence of $100\ \mu\text{M}$ cysteine (clear circles). Time points represent 0, 5, 10, 15, 20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75 and 80 min. The inset show the visual fluorescence color of probe **1** ($10\ \mu\text{M}$) before (left) and after (right) incubation with cysteine for 60 min (UV lamp, 365 nm).

NMR and MS data for compounds




Figure S5. ^1H NMR spectrum of compound **2** in CD_3OD .




Figure S6. ^1H NMR spectrum of compound 3 in $\text{DMSO-}d_6$.




Figure S7. ¹H NMR spectrum of compound 4 in CDCl₃.




Figure S8. ¹H NMR spectrum of compound 5 in CDCl₃.




Figure S9. ^1H NMR spectrum of compound 6 in $\text{DMSO-}d_6$.




Figure S10. ESI-MS spectrum of compound **6**.




Figure S11. ¹H NMR spectrum of compound **1** in CDCl₃.




Figure S12. ¹³C NMR spectrum of compound **1** in CDCl₃.




Figure S13. DEPT135 spectrum of compound **1** in CDCl₃.




Figure S14. ESI-MS spectrum of compound **1**.