

Supplementary materials

Fluorescence sensing of tyrosinase activity based on amines riched carbon dots through direct interaction in homogeneous system: detection mechanism and application

Jianzhong Zhang^{a,#}, Yuyuan Chen^{b,#}, Zongfu Zheng^c, Zhenzhen Wang^b, Yanjie Zheng^b, Xinhua Lin^b, Shaohuang Weng^{b,*}

1. Department of Gastric Surgery, Fujian Medical University Union Hospital, Fuzhou 350001, China

2. Department of Pharmaceutical Analysis, The Higher Educational Key Laboratory for Nano Biomedical Technology of Fujian Province, School of Pharmacy, Fujian Medical University, Fuzhou 350122, China

3. 476 Hospital of PLA, Fuzhou 350002, China

* Correspondence: shweng@fjmu.edu.cn

Both of these authors contributed equally to this work.




Figure S1. TEM image of N-CDs. The insets are the size distribution (upper) and high resolution TEM (below) of N-CDs.




Figure S2. The lifetime decay of three kinds of carbon dots in the absence (black) and presence (red) of the mixture of TYR and DA. A: N-CDs, B: C-2, C: C-1.

Table S1. The comparison of lifetime of different carbon dots quenching or not quenching by the mixture of DA and TYR.

Samples	N-CDs	C-2	C-1
CDs	8.82 ns	8.86 ns	12.79 ns
CDs+DA+TYR	8.67 ns	8.01 ns	12.37 ns




Figure S3. High-resolution XPS spectra of N-CDs (A) and C-2 (B) for N1s.




Figure S4. (A) Fluorescence intensity changes of N-CDs of the sensing platform with 10 U/mL TYR and 0.2 mM DA in different pH PBS, (B) Fluorescence intensity changes of N-CDs in the presence of 0.2 mM DA in different pH PBS.




Figure S5. (A) Fluorescence intensity changes of N-CDs with 10 U/mL TYR and 0.2 mM DA in different reacting time and temperature, (B) Fluorescence intensity changes of N-CDs with 10U/mL TYR and 0.2 mM DA in 25, 30 and 37°C after reacting 120 min.




Figure S6. (A) Fluorescence emission spectra of the sensing system to different concentrations of DA with 10 U/mL TYR. (B) The relationship between the fluorescence quenching efficiency and the DA concentrations ranging from 0.01 to 2.0 mM with 10 U/mL TYR.

Table S2. Comparison of analytical performance of several fluorescence assays for TYR

Applied Materials	Linear range (U/mL)	Detection limit (U/mL)	References
Dopa-Au/Ag NCs	0.045-0.3195	0.0135	[S1]
Dopa-CQDs	0.0444-0.7111	0.0177	[S2]
Dopa-CQDs	0.0232-0.7935	0.007	[S3]
Dopa-Mn:ZnS QDs	0.018-0.36	0.0018	[S4]
PMNPs	0.4-7.0	0.1	[S5]
SiNPs	0.6-6.0	0.14	[S6]
N-CDs	0.05-6.0	0.039	This Work

References

- [S1] H. Ao, Z. S. Qian, Y. Y. Zhu, M. Z. Zhao, C. Tang, Y. Y. Huang, H. Feng and A. J. Wang, *Biosens. Bioelectron.*, 2016, 86, 542.
- [S2] J. J. Hu, X. L. Bai, Y. M. Liu and X. Liao, *Anal. Chim. Acta*, 2017, 995, 99.
- [S3] L. J. Chai, J. Zhou, H. Feng, C. Tang, Y. Y. Huang and Z. S. Qian, *ACS Appl. Mater. Interfaces*, 2015, 7, 23564.
- [S4] F. M. Zhang, L. B. Xu, Q. N. Zhao, Y. Sun, X. H. Wang, P. Y. Ma and D. Q. Song, *Sensors and Actuators B: Chemical*, 2018, 256, 1069.

[S5] G. Y. Liu, J. H. Zhao, S. S. Lu, S. Wang, J. Sun and X. R. Yang, ACS Sens., 2018, 3, 1855.

[S6] Y. Z. Ding, W. F. Wang, T. Chai, Y. Qiang, Y. P. Shi and J. L. Yang, Talanta, 2019, 197, 113.