Supplementary Information

Ultrasensitive Nucleic Acid Detection Using Confocal Laser Scanning Microscope with High Crystalline Silver Dendrites Prepared by Our New Developed Technique as a Single Particle Indicator

Xuan Yang, Xuping Sun, Zhaozi Lv, Weiwei Guo, Yan Du and Erkang Wang*

^a State Key Laboratory of Electroanalytical Chemistry, Changchun Institute of Applied Chemistry, Chinese Academy of Sciences, Changchun, Jilin, 130022, P. R. China. Graduate School of the Chinese Academy of Sciences, Changchun 130022, Jilin, P. R. China

* Corresponding author: Prof. Erkang Wang, Fax: +86-431-85689711, Tel: +86-431-85262003, Email: <u>ekwang@ciac.jl.cn</u>

Chemicals and materials

All DNA strands were synthesized by Shanghai Sangon Biotechnology Co. Ltd (Shanghai, China). Hydroxylamine (50 wt.% solution in water) was purchased from Sigma-Aldrich (St. Louis, MO, USA). All the other chemicals were of analytical reagent grade and were used as received without further purification. All the solutions were prepared with double-distilled water purified by a Milli-Q system (Millipore, Bedford, MA, USA) and stored at 4 °C. The concentration of DNA was determined using the 260 nm UV absorbance and the corresponding extinction coefficient.

Apparatus

UV-visible absorbance spectra was recorded on a Cary 500 scan UV-vis-NIR spectrophotometer (Varian, Harbor City, CA) at room temperature. Scanning electron microscopy (SEM) measurements were made on a S4800 scanning electron microscope (Hitachi, Janpan) at an accelerating voltage of 10.0 kV. Transmission electron microscopy (TEM) measurements were made on a FEI TECNAI G² transmission electron microscope (Netherlands) operated at an accelerating voltage of 120 kV. X-ray diffraction (XRD) measurements were made on a PW1700 X-ray diffractometer. Confocal laser scanning microscope (Germany).

Preparation: Sample 1 of silver micro-dendrites were prepared as follows: in a typical preparation, 112 μ L of 1.63 M NH₂OH was diluted into 1.5 mL water, followed by adding 100 μ L 0.45 M AgNO₃ aqueous solution into above solution. Sample 2 of silver

micro-dendrties used in the separation assays were prepared as follows: in a typical preparation, 896 μ L of 1.63 M NH₂OH was diluted to 652 μ L water, followed by adding 800 μ L of 0.45 M AgNO₃ aqueous solution into above solution, leading to a large amount of precipitate in several seconds.



Figure S1. Typical SEM images of silver particles of sample 1 a), b) shows the TEM image of single silver dendrite and the corresponding SAED pattern (inset).



Figure S2. XRD pattern of silver micro-dendrites of sample 1.



Figure S3. SEM images of silver particles obtained under otherwise identical conditions

used for preparing sample 2 but with sonication involved.



Figure S4. SEM images of the silver particles obtained at 50 , under otherwise identical conditions used for preparing sample 2.



Figure S5. SEM images of silver particles obtained with molar ratio (A, B) 8:1 and (C, D)

2:1 of NH₂OH:AgNO₃, under otherwise identical conditions used for preparing sample 2.



Figure S6. The SEM image of sample 2 (A) and the size distribution of sample 2 (B), 20 particles are measured to get the size distribution.

Probe 1 DNA	5' -SH-(CH2)6TCCTCCAGG-3'
Probe 2 DNA	5'-CGCAAATACTCAAA-FAM-3'
Target DNA	5'- GAGTATTTGCGCATGGCCTGGAGGA-3'
One nucleotide mismatched DNA	5' - GAGTA <u>G</u> TTGCGCATGGCCTGGAGGA-3'
Two nucleotide mismatched DNA	5'- GAGTA <u>G</u> TTGCGCATGGCCTG <u>A</u> AGGA-3'
Four nucleotide mismatched DNA	5'- GAGTA <u>GG</u> TGCGCATGGCCT <u>AA</u> AGGA-3'
The absolutely mismatched DNA	5'- CTTCTTGAATTCGTATAGCGCCTCT-3'

Table 1. The sequences of the DNA used in our experiments.