

Supporting information


Index

1. The monomeric unit of polymer Xanthan	S1
2. Materials and Methods	S2
3. SERF pulse sequence	S3
4. Experimental Details	S4
5. Indication of mesophase formation	S5
6. Proton NMR spectrum of (R/S)-β-butyrolactone	S6
7. Temperature dependence of quadrupole split doublet	S7
8. Transition from isotropic to anisotropic phase	S8
9. Variation of quadrupole split doublet as a function of quantity of xanthan	S9
10. Sustainability of anisotropic phase even at 370 K	S10
11. Prochiral discrimination of Glycine-D₂	S11
12. The temperature reversibility of the mesophase	S12
13. The ²H NMR spectrum of a scalemic mixture	S13

Materials and Methods:

Xanthan polymer purchased from Sigma Aldrich [catalog No. G1253, molecular weight not known]. The liquid crystalline phase of xanthan was prepared according to the known procedure [T. Sato, T. Norisuye, H. Fujita, *Polymer J.*, **1984**, 16, 341-350.], by initially dissolving 700 mg of xanthan obtained from Sigma-Aldrich, in 35 ml of D₂O and exposing this aqueous solution to 20 kHz sonic irradiation for about six hours (sonication done with 20 sec MW pulse and 10 sec interval). Sonication was done under ice bath to avoid sample heating. The sonicated xanthan solution was centrifuged at 8000 rpm for 30 minutes and the supernatant was precipitated using large quantity of acetone and subsequently vacuum dried for about 1-2 hours. The liquid crystalline phase was prepared by adding 120 mg of this aqueous xanthan solution and 380 µl of 0.1M NaCl to NMR tube and centrifuging back and forth to get visibly homogeneous phase. The alignment of the phase was confirmed by detecting the quadrupole split doublet of the solvent D₂O.

SElective ReFocusing Experimental pulse Sequence



All the pulses are selective. Phases are, $\Phi_1=X$, $\Phi_2=\Phi_R=X, -X, -X, X$

Experimental Details

1. ^2H NMR of solvent D_2O :

The solution was prepared with 70 mg of xanthan in 300 μl of 0.1 M NaCl solution. The one dimensional spectrum was recorded using Bruker DRX-500 NMR spectrometer. The 16 transients were accumulated with a relaxation delay of 1 s. The 90° pulse width is 9.8 μs . 2 k time domain data were collected and Fourier transformed without using any window function.

2. ^1H NMR of (*R/S*)- β -butyrolactone

The solution was prepared with 120 mg of xanthan, 10 μl of (*R/S*)- β -butyrolactone in 380 μl of 0.1 M NaCl solution. The one dimensional spectrum was recorded using Bruker DRX-500 NMR spectrometer. The single transient was collected. The 90° pulse width is 8.25 μs . 3 k time domain data were collected and Fourier transformed without using any window function.

3. SERF Experiment on (*R/S*)- β -butyrolactone


The solution was prepared with 120 mg of xanthan, 10 μl of (*R/S*)- β -butyrolactone in 380 μl of 0.1 M NaCl solution. The 2D SERF spectra were recorded using Bruker DRX-500 NMR spectrometer. 90° pulse width is 8.25 μs . Two transients were accumulated with 256 t1 increments. The size of the 2D data matrix is 300 X 1024. The data was Fourier transformed using unshifted sine bell window function in both the dimensions. 3 k time domain data were collected and Fourier transformed without using any window function.

4. HSQC experiment on a tri peptide in xanthan medium


The solution of tripeptide Gly-Glu-Cys was prepared taking 5 mg of the solute in 50 mg of xanthan and 380 μl of NaCl solution. The size of the 2D data matrix was 512 X 1536. Sixteen free induction decays were accumulated and Fourier transformed using sinebell window functions in both the dimensions.

Indication of mesophase formation:

The ^2H quadrupole doublet of solvent D_2O , indicating the alignment of guest molecule. 120 mg of xanthan and 380 μl of 0.1M NaCl solution was taken. The one dimensional ^2H spectrum was recorded in Bruke DRX 500 NMR spectrometer. The doublet separation corresponds to 140 Hz indicating the formation of the mesophase.




The ^1H NMR spectrum of (*R/S*)- β -butyrolactone




Temperature dependence of quadrupole split doublet

The temperature dependence of quadrupole split doublet separation of solvent D₂O aligned in the mesophase of xanthan.




Transition from isotropic to anisotropic phase

The ^2H -NMR spectrum of solvent D_2O as a function of concentration of Xanthan in 380 μl of 0.1 M NaCl indicating the minimum quantity of xanthan required to obtain mesophase. The spectra are recorded at 61.4 MHz on AVIII 400- Bruker NMR spectrometer.




Variation of quadrupole split doublet separation as a function of quantity of xanthan

The variation of deuterium quadrupole splitting of solvent D₂O as a function of the quantity of xanthan, in 380 μl of 0.1M NaCl.




Sustainability of anisotropic phase even at 370 K

The ^2H NMR spectra of solvent D_2O at different temperatures between 280- 370 K, varied in steps of 10 degrees. Spectra recorded at 76.73 MHz frequency on Bruker DRX-500 NMR spectrometer. It is apparent that the anisotropic medium is sustained over a wide range of temperature, even when the boiling point of water is approached. Sample prepared using 70 mg of Xanthan in 300ul of 0.1 M NaCl. The measured quadrupolar splitting varied from 138 (at 280 k) to 89 Hz (at 370 K).




Prochiral discrimination of Glycine-D₂

Additional example of prochiral discrimination in Glycine-D₂. The 10 mg of glycine-D₂ was dissolved in 70 mg of xanthan in 300 μ l of 0.1 M NaCl solution. The ²H NMR spectra were recorded at different temperatures. The pro-chiral discrimination was observed from 300 -370 k. The peaks marked * are the doublets from natural abundant ²H of solvent water.



The temperature reversibility of the mesophase:

The quadrupole doublet separation was monitored for the solvent D₂O from 280-370 K and 370-280 K, in steps of five degrees. The process of measurement was repeated twice and each temperature sufficient time was given to attain temperature stability. The plots confirm not only the sustainability of the mesophase over a wide temperature range, but also reversibility of the mesophase.



The ^2H NMR spectrum of a scalemic mixture of D/L-alanine- d_3

The scalemic mixture of D/L-alanine- d_3 with 33.3 % enrichment of L enantiomer was prepared to measure the enantiomeric excess. The ^2H NMR spectrum was recorded at ambient temperature at 76.73 MHz. The ratiometric analysis of the intensities provided the enantiomeric excess of 32.78 %.

