

Supporting Information:




Fig. S1. Determining Lutrol CMT as a function of concentration.

Table S1. DSC parameters of Lutrol micellization.

Lutrol Concentration pH 7, IS 0.02 mM	T_{onset}, °C	Lutrol CMT T_m, °C
1.683	22.84	26.8
0.842	24.7	28.3
0.421	25.0	29.7
0.168	26.0	31.2
0.0842	27.6	32.5
0.0421	28.1	33.6
0.0084	30.4	37.0
0.0042	32.9	39.0