

1 **Supporting Information**

2

3 **LC-MS in Combination with DMBA Derivatization for Sialic Acid Speciation**
4 **and Distribution Analysis in Fish Tissues**

5

6 Junjie Du ^a, Qiwei Zhang ^{*a}, Jianjun Li ^b, Qi Zheng ^a

7

8 ^aKey Laboratory of Optoelectronic Chemical Materials and Devices of Ministry of Education,
9 Institute for Interdisciplinary Research, Institute of Environment and Health, Jiangnan University,
10 Wuhan 430056, China

11 ^bNational Research Council Canada, 100 Sussex Drive, Ottawa, ON, Canada K1A 0R6

12

13 *Corresponding authors:

14 Dr. Qiwei Zhang, School of Chemical and Environmental Engineering, Jiangnan University,
15 Wuhan, 430056, China, E-mail: cvcon@msn.cn


16

17 Dr. Qi Zheng, School of Chemical and Environmental Engineering, Jiangnan University, Wuhan,
18 430056, China, E-mail: zq_1101@sina.com

19

Table of Contents


21		
22	Fig. S1. The schematic diagram of DMBA derivatization of sialic acid.....	S3
23	Fig. S2. The effect of LC gradients on the separation of different sialic acids.....	S4
24	Fig. S3. The effect of column temperature on the separation of sialic acids.....	S5
25	Fig. S4. LC-MS/MS analysis of DMBA-labeled mono- <i>O</i> -acetylated <i>N</i> -glycolylneuraminic acids	
26	(<i>m/z</i> 468.20) from crucian carp blood.....	S6
27	Fig. S5. LC-MS analysis of DMB-labeled sialic acids from crucian carp blood.....	S7
28	Fig. S6. LC-MS/MS analysis of DMBA-labeled di- <i>O</i> -acetylated <i>N</i> -acetylneuraminic acids (<i>m/z</i>	
29	494.21) from crucian carp blood	S8
30	Fig. S7. LC-MS/MS analysis of DMBA-labeled tri- <i>O</i> -acetylated <i>N</i> -acetylneuraminic acids (<i>m/z</i>	
31	536.22) from crucian carp blood.....	S9
32	Table S1. Detected sialic acids in crucian carp blood with their relative retention times and	
33	characteristic ions.....	S10
34		
35		



36

37 **Fig. S1.** The schematic diagram of DMBA derivatization of sialic acid.


38



39

40 **Fig. S2.** The effect of LC gradients on the separation of different sialic acids. XICs are the sum of
 41 ions at m/z 452.20, 468.20, 494.21 and 536.22. The mobile phase B was increased from 10 % at 3
 42 min to 20 % at different time. (A) 10 min; (B) 15 min; (C) 20 min; (D) 25 min; (E) 30 min.

43




44

45 **Fig. S3.** The effect of column temperature on the separation of sialic acids. The XICs are the sum
 46 of ion at m/z 452.20, 468.20, 494.21 and m/z 536.22. (A) 35 °C; (B) 40 °C; (C) 45 °C; (D) 50 °C;
 47 (E) 55 °C.

48


49



50

51 **Fig. S4.** LC-MS/MS analysis of DMBA-labeled mono-*O*-acetylated *N*-glycolylneuraminic acids
 52 from crucian carp blood. (A) MS/MS spectrum extracted at Peak 1; (B) MS/MS spectrum extracted
 53 at Peak 2; (C) MS/MS spectrum extracted at Peak 3; (D) MS/MS spectrum extracted at Peak 4.


54




55

56 **Fig. S5.** LC-MS analysis of DMB-labeled sialic acids from crucian carp blood. (A) XIC at m/z
 57 426.15 (Neu5Ac); (B) XIC at m/z 442.15 (Neu5Gc); (C) XIC at m/z 468.16 (NeuAc₂); (D) XIC at
 58 m/z 510.17 (NeuAc₃); (E) XIC at m/z 552.18 (NeuAc₄).

59



61 **Fig. S6.** LC-MS/MS analysis of DMBA-labeled di-*O*-acetylated *N*-acetylneuraminic acids (*m/z*
 62 494.21) from crucian carp blood. (A) MS/MS spectrum extracted at Peak 11; (B) MS/MS
 63 spectrum extracted at Peak 12; (C) MS/MS spectrum extracted at Peak 13; (D) MS/MS spectrum
 64 extracted at Peak 14; (E) MS/MS spectrum extracted at Peak 15; (F) MS/MS spectrum extracted
 65 at Peak 16.



66
 67 **Fig. S7.** LC-MS/MS analysis of DMBA-labeled tri-*O*-acetylated *N*-acetylneuraminic acids (*m/z*
 68 536.22) from crucian carp blood. (A) MS/MS spectrum extracted at Peak 17; (B) MS/MS spectrum
 69 extracted at Peak 18; (C) MS/MS spectrum extracted at Peak 19.

70

71 **Table S1.** Detected sialic acids in crucian carp blood with their relative retention times and
 72 characteristic ions.

Peak No.	Sialic acid	Abbreviation	Retention time (min)	RRT ₁	RRT ₂	m/z [M+H] ⁺	Characteristic ions
1	7- <i>O</i> -Acetyl- <i>N</i> -glycolylneuraminic acid	Neu5Gc7Ac	7.49	1.01	0.98	468.20	213.10, 267.11, 279.13, 297.12, 450.18
2	8- <i>O</i> -Acetyl- <i>N</i> -glycolylneuraminic acid	Neu5Gc8Ac	8.27	1.11	1.09	468.20	213.10, 243.11, 279.11, 297.12, 450.18
3	9- <i>O</i> -Acetyl- <i>N</i> -glycolylneuraminic acid	Neu5Gc9Ac	8.42	1.13	1.21	468.20	213.10, 243.11, 279.11, 297.12, 450.18
4	4- <i>O</i> -Acetyl- <i>N</i> -glycolylneuraminic acid	Neu4Ac5Gc	8.79	1.18	1.23	468.20	213.10, 243.11, 267.11, 279.11, 297.12, 408.17
5	<i>N</i> -Acetylneuraminic acid	Neu5Ac	7.44	1.00	1.00	410.19	213.10, 243.11, 267.11, 279.11, 285.12, 297.12, 392.18
6	<i>N</i> -Glycolylneuraminic acid	Neu5Gc	6.60	0.89	0.84	426.19	213.10, 243.11, 267.11, 279.11, 285.12, 297.12, 408.17
7	7- <i>O</i> -Acetyl- <i>N</i> -acetylneuraminic acid	Neu5,7Ac ₂	8.05	1.08	1.12	452.20	201.10, 213.10, 267.11, 279.11, 285.12, 297.12, 434.19
8	8- <i>O</i> -Acetyl- <i>N</i> -acetylneuraminic acid	Neu5,8Ac ₂	9.32	1.25	1.40	452.20	213.10, 243.11, 279.11, 297.12, 434.19
9	9- <i>O</i> -Acetyl- <i>N</i> -acetylneuraminic acid	Neu5,9Ac ₂	9.64	1.30	1.59	452.20	213.10, 243.11, 279.11, 297.12, 434.19

10	4- <i>O</i> -Acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5Ac ₂	10.31	1.39	1.70	452.20	201.10, 213.10, 243.11, 267.11, 279.11, 285.12, 297.12, 392.18
11	7,8-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu5,7,8Ac ₃	10.56	1.42	1.77	494.21	213.10, 279.11, 297.12, 476.20
12	7,9-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu5,7,9Ac ₃	10.96	1.47	2.03	494.21	213.10, 279.11, 297.12, 476.20
13	4,7-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5,7Ac ₃	11.32	1.52	2.51	494.21	213.10, 279.11, 285.12, 297.12, 434.19, 476.20
14	8,9-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu5,8,9Ac ₃	13.02	1.75	2.73	494.21	213.10, 243.11, 279.11, 297.12, 434.19
15	4,8-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5,8Ac ₃	13.65	1.83	3.28	494.21	213.10, 243.11, 279.11, 297.12
16	4,9-Di- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5,9Ac ₃	15.72	2.11	4.75	494.21	213.10, 243.11, 279.11, 297.12, 476.20
17	7,8,9-Tri- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu5,7,8,9Ac ₄	14.88	2.00	3.81	536.22	213.10, 279.11, 297.12
18	4,7,8(9)-Tri- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5,7,8(9)Ac ₄	16.62	2.23	4.99	536.22	213.10, 279.11, 297.12
19	4,8,9-Tri- <i>O</i> -acetyl- <i>N</i> -acetylneuraminic acid	Neu4,5,8,9Ac ₄	20.75	2.79	9.51	536.22	213.10, 243.11, 279.11, 297.12

73 Relative retention time (RRT) is the ratio of the retention time of analyte relative to that of Neu5Ac.

74 RRT_{1S} are the values from DMBA derivatization and obtained in this study. RRT_{2S} are the values

75 from DMB derivatization that obtained from reference #28.

76