

Supplementary Information

Gas Chromatography-Tandem Mass Spectrometry based Detection of Half Nitrogen Mustards in Plasma as a New Biomarker of Nitrogen Mustards exposure

Buddhadeb Chandra^a, Kanchan Sinha Roy^{ab}, MahabulShaik^a, Chandrakant Waghmare^a, Ramarao Golime^a, and Meehir Palit^{a*}

^a*VERTOX Laboratory, Defence Research and Development Establishment,
Jhansi Road, Gwalior-474002, Madhya Pradesh, India.*

^b*Present Address: Department of Chemistry, University of Waterloo, Waterloo, ON, Canada,*

N2L 3G1

*Corresponding Author:

Fax: (+) 0751-2341148


E-mail: meehirpalit@rediffmail.com

Contents


Sl. No.	Caption	Figure/table no.
1.	Different approaches for the synthesis of halfNMs	Scheme S1
2.	Hydrolysis kinetics of NM salts: Log (C ₀ /C) vs time (min) coordinates for the hydrolysis of NMs salts in 2:1 (v/v), acetone-water	Figure S1
3.	Observed half-life time of salts of NM in 2:1 (v/v), acetone-water solution.	Table S1
4.	Screening of the temperature for the heptafluorobutyrylation of half NMs.	Figure S2
5.	Fine optimization of the temperature for the heptafluorobutyrylation of halfNMs.	Figure S3
6.	Study of the effect of time on the heptafluorobutyrylation of halfNMs.	Figure S4

7.	PCI-MS spectra of a) halfHN2HFB, b) halfHN1HFB, c) 1halfHN3HFB and d) 2halfHN3HFB respectively.	Figure S5
8.	Collision energy optimization during SRM based method development for half NMs.	Figure S6a, S6b, S6c and S6d
9.	Parameters for Multiple Reaction Monitoring (MRM)	Table S2
10.	Calibration curve for half NMs taking plasma as a matrix.	Figure S7
11.	SRM chromatograms, with transition 348 → 241 (quantifier ion) and 312 → 241 (qualifier ion), obtained during <i>in-vitro</i> analysis of halfHN1 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.	Figure S8
12.	SRM chromatograms, with transition 382 → 241 (quantifier ion) and 346 → 241 (qualifier ion), obtained during <i>in-vitro</i> analysis of 1halfHN3 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.	Figure S9
13.	SRM chromatograms, with transition 560 → 241 (quantifier ion) and 346 → 241 (qualifier ion), obtained during <i>in-vitro</i> analysis of 2halfHN3 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.	Figure S10
14.	Response of halfHN2 with respect to time after exposure of Wister rat with 1X LD ₅₀ of HN2 (<i>in vivo</i> analysis).	Figure S11

1st approach:



2nd approach:



Scheme S1: Different approaches for the synthesis of halfNMs.




Figure S1: Log (C_0/C) vs time (min) co-ordinates for the hydrolysis of (a) HN1 salt (b) HN2 salts and (c) HN3 salts during hydrolysis of each salt in 2:1 (v/v), acetone-water solution. Here C_0 is the concentration of NM at time = 0 min. and C is the concentration of NM at time = t min.

Table S1: Observed half-life time of salts of NM in 2:1 (v/v), acetone-water solution.

SL No.	Salts of NM	Slope (m)	$K = m \times 2.303$	Half-life time ($t_{1/2}$) = $0.693/K$
1	HN1	0.00111	0.00255633	271.09 min
2	HN2	0.00116	0.00267148	259.41 min
3	HN3	0.00385	0.00886655	78.16 min




Figure S2: Screening of the temperature for the heptafluorobutyrylation of half NMs.




Figure S3: Fine optimization of the temperature for the heptafluorobutyrylation of halfNMs.




Figure S4: Study of the effect of time on the heptafluorobutyrylation of halfNMs.




Figure S5: PCI-MS spectra of a) halfHN2HFB, b) halfHN1HFB, c) 1halfHN3HFB and d) 2halfHN3HFB respectively




Figure S6a: Optimization of collision energy for the SRM transitions halfHN1HFB.




Figure S6b: Optimization of collision energy for the SRM transitions halfHN2HFB.




Figure S6c: Optimization of collision energy for the SRM transitions 1halfHN3HFB.




Figure S6d: Optimization of collision energy for the SRM transitions 2halfHN3HFB.

Table S2: Parameters for Multiple Reaction Monitoring (MRM)

Analytes	Nature of ions	Ionic transition	Fragmentation energy (eV)
halfHN1HFB	Q	348 > 241	30
	q	312 > 241	30
halfHN2HFB	Q	334 > 241	25
	q	298 > 241	25
1halfHN3HFB	Q	382 > 241	30
	q	346 > 241	30
2halfHN3HFB	Q	560 > 241	40
	q	346 > 241	25
HCB (IS)	Q	284.9 > 249.8	20
	q	284.9 > 248	20




Figure S7: Calibration curve for (a) halfHN1, (b) halfHN2, (c) 1halfHN1 and (d) 2halfHN3 taking plasma as a matrix.




Figure S8: SRM chromatograms, with transition 348 → 241 (quantifier ion) and 312 → 241 (qualifier ion), obtained during *in-vitro* analysis of halfHN1 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.




Figure S9: SRM chromatograms, with transition 382 → 241 (quantifier ion) and 346 → 241 (qualifier ion), obtained during the *in-vitro* analysis of 1halfHN3 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.




Figure S10: SRM chromatograms, with transition 560 → 241 (quantifier ion) and 346 → 241 (qualifier ion), obtained during the *in-vitro* analysis of 2halfHN3 in human plasma; (a) blank sample, (b) HN1 spiked sample, and (c) Standard reference sample.




Figure S11: Response of halfHN2 with respect to time after exposure of Wister rat with 1X LD₅₀ of HN2 (in *vivo* analysis).