

SUPPLEMENTARY INFORMATION

A simple three-dimensional gut model constructed in a restricted ductal microspace induces intestinal epithelial cell integrity and facilitates absorption assays

Tadaaki Nakajima^a, Katsunori Sasaki^b, Akihiro Yamamori^b, Kengo Sakurai^b, Kaori Miyata^b, Tomoyuki Watanabe^{b, †}, Yukiko T. Matsunaga^{a,*}.

^a Institute of Industrial Science, The University of Tokyo, 4-6-1 Komaba, Meguro-ku, Tokyo 153-8505, Japan

^b Sumitomo Chemical, Co., Ltd., 27-1, Shinkawa 2-chome, Chuo-ku, Tokyo 104-8260, Japan

[†]Current address: Matsuura Yakugyo Co., Ltd., 24-21, Enjo-cho, Showa-ku, Nagoya 466-0054, Japan

* Correspondence to:

Yukiko T. Matsunaga, Ph.D.

Institute of Industrial Science, The University of Tokyo 4-6-1 Komaba, Meguro-ku, Tokyo 153-8505, Japan

Tel.: +81-3-5452-6470; Fax: +81-3-5452-6471

E-mail: mat@iis.u-tokyo.ac.jp

Table of Contents

1. Supplementary Figures

Figure S1. Confocal images for EdU assay and immunofluorescence for differentiation marker proteins in 2D or 3D gut model cultured for 16 days.

Figure S2. Fluorescent images for leak check at absorption assay in 2D and 3D gut models.

2. Supplementary Table

Figure S1. Primer sets for qRT-PCR.

1. Supplementary Figures




Figure S1. Confocal images for EdU assay and immunofluorescence for differentiation marker proteins in 2D or 3D gut model cultured for 16 days. (A) Confocal microscopies of EdU assay. EdU-positive nuclei (red) and nuclei (blue) were visualized. (B) Binarized immunofluorescent images for villin (red) and lysozyme (green). Scale bar: 50 μm .




Figure S2. Fluorescent images for leak check at absorption assay in 2D and 3D gut models.

(A) Fluorescent and bright-field (BF) images at 20 min after FITC-4 kDa dextran treatment in the 3D gut model cultured for 4 or 16 days with or without 10 mM NaN_3 treatment. Scale bar: 500 μm . (B) Fluorescent images at 20 min after FITC-4 kDa dextran treatment in the bottom well of the 2D gut model cultured for 16 days with or without 10 mM NaN_3 treatment to check for absence of leakage. Background: without FITC-dextran treatment. In the region on the filter, fluorescent intensity was saturated in all samples. Scale bar: 200 μm .

2. Supplementary Table

	Forward	Reverse
<i>ACTB</i>	GGGCATGGGTCAGAAGGATT	AGGTCTCAAACATGATCTGGGT
<i>LGR5</i>	AGCCTTCAATCCCTGCGTCT	TAACGCATTGTCATCCAGCCA
<i>LYZ</i>	ATCAGCCTAGCAAAGTGGATG	TGACAGGCATTAAGTCTCCT
<i>MUC2</i>	ACGTGGCTGTTTCAGGACTAC	TCCGTCCTCCCATGAAGAT
<i>VIL1</i>	GATATGGAGGATCGAGGCCAT	GTCCTGGCCAATCCAGTAGT
