

Supporting Information

Auto-affitech: An Automatic Binding Ligand Affinity Evaluation Platform Using Digital Microfluidics with Bidirectional Magnetic Separation Method

Jingjing Guo^a, Li Lin^a, Kaifeng Zhao^a, Yanling Song^b, Mengjiao Huang^a, Zhi Zhu^a, Leiji Zhou^{*a}, Chaoyong Yang^{*ab}

a MOE Key Laboratory of Spectrochemical Analysis & Instrumentation, Key Laboratory of Chemical Biology of Fujian Province, State Key Laboratory of Physical Chemistry of Solid Surfaces, Collaborative Innovation Center of Chemistry for Energy Materials, College of Chemistry and Chemical Engineering, Xiamen University, Xiamen, 361005, China

b Institute of Molecular Medicine, Renji Hospital, School of Medicine, Shanghai Jiao Tong University, Shanghai, 200127, China

E-mail: cyyang@xmu.edu.cn, ljzhou@xmu.edu.cn

1. Buffers

For the EpCAM-aptamer system, the buffers used were as follows: binding buffer (0.55 mM MgCl_2 in PBS, pH=7.2 - 7.4), washing buffer (0.55 mM MgCl_2 , 0.1% BSA in PBS, pH=7.2 - 7.4). For the H5N1 antibody-antigen system, the buffers used were as follows: washing buffer (TBS, 20 mM Tris, 150 mM NaCl, pH=7.2 - 7.4), binding buffer (0.5% BSA in washing buffer, pH=7.2 - 7.4). All the buffers used in the DMF system also included 0.1% F127.

2. The calibration curve of HRP-beads for bead retention efficiency characterization




Fig. S1 The calibration curve of the absorbance response to different concentrations of HRP-beads for characterizing bead retention efficiency. Error bars indicate the standard deviations of three samples.

3. The feasibility of SA-HRP with substrate




Fig. S2 The PMT response to different concentrations of SA-HRP. Error bars indicate the standard deviations of three samples.

4. Dissociation constant (K_d) measured by flow cytometry




Fig. S3 (A) K_d fitting curve of SYL3C aptamer against EpCAM measured by flow cytometry. (B) K_d fitting curve of H5N1 antibody against H5N1 antigen measured by flow cytometry. Error bars indicate the standard deviations of three samples.

5. Comparison of Auto-affitech and flow cytometry in processing time for one detection cycle

Methods	Processing time					Total time
	Incubating with targets	Washing	Incubating with SA-PE /SA-HRP	Washing	Substrate	
Flow cytometry	30 min	5 min	30 min	5 min	/	70 min
Auto-affitech	10 min	1 min	10 min	2 min	30 s	23.5 min

Table S1 Comparison of Auto-affitech and flow cytometry in processing time for one detection cycle.

6. Comparison of Auto-affitech and flow cytometry in reagent consumption for one detection cycle

Methods	Reagent consumption					Total volume
	Ligand reagent	Washing	Incubating with SA-PE /SA-HRP	Washing	Substrate	
Flow cytometry	100 μ L	300 μ L	100 μ L	300 μ L	/	800 μ L
Auto-affitech	1.7 μ L	6.8 μ L	1.7 μ L	10.2 μ L	3.4 μ L	23.8 μ L

Table S2 Comparison of Auto-affitech and flow cytometry in reagent consumption for one detection cycle.