

## Supporting Information

# 4-Hydrazinoquinazoline acting as reactive matrix for rapid and sensitive analysis of neutral and sialylated glycans using MALDI MS

Ling Ling, Shaoning Yu, Chuanfan Ding\*

*Key Laboratory of Advanced Mass Spectrometry and Molecular Analysis of Zhejiang Province, Institute of Mass Spectrometry, School of Material Science and Chemical Engineering, Ningbo University, Ningbo 315211, China*

\*Corresponding Author:


Chuanfan Ding

Key Laboratory of Advanced Mass Spectrometry and Molecular Analysis of Zhejiang Province,


Institute of Mass Spectrometry,

School of Material Science and Chemical Engineering, Ningbo University,  
Ningbo 315211, China


Email: dingchuanfan@nbu.edu.cn




**Fig. S1.** MALDI mass spectrum of 4-HQ.




**Fig. S2.** UV-vis absorption spectra of 4-HQ and DHB (each at concentrations of  $10^{-4}$  M in water).




**Fig. S3.** S/N ratio variations of DP7 (5 pmol) with increasing storage time: a) 4-HQ matrix solution was stored at 4°C; b) sample was placed on MALDI target at room temperature.



**Fig. S4.** MALDI mass spectra of different amount of DP7 using **a)** 4-HQ as a reactive matrix and **b)** DHB as the matrix.



**Fig. S5.** MALDI mass spectra of A3 glycan using **a)** 4-HQ as a reactive matrix and **b)** DHB as the matrix.

















**Fig. S6.** MALDI mass spectra of N-glycans released from human serum using 4-HQ as a reactive matrix in **a)** positive and **b)** negative ion mode.

**Table S1.** List of identified N-glycans released from ovalbumin using 4-HQ or DHB as matrix.

The m/z values are corresponding to the ions of  $[M+Na]^+$ .

Peak No. <sup>a</sup>	Theoretical m/z of native glycans	Observed m/z of native glycans	Theoretical m/z of 4-HQ-labeled glycans	Observed m/z of 4-HQ-labeled glycans	Composition
1	933.3	933.3	1075.4	1075.4	
2	1095.4	-	1237.4	1237.4	
3	1136.4	1136.4	1278.5	1278.4	
4	1257.4	1257.4	1399.5	1399.5	
5	1298.4	1298.4	1440.5	1440.5	
6	1339.5	1339.5	1481.5	1481.5	
7	1419.5	1419.5	1561.5	1561.5	














8	1460.5	-	1602.5	1602.5	
9	1501.5	1501.5	1643.6	1643.6	
10	1542.6	1542.6	1684.6	1684.6	
11	1581.5	-	1723.6	1723.6	
12	1663.6	1663.6	1805.6	1805.6	
13	1704.6	1704.6	1846.7	1846.7	
14	1745.6	1745.6	1887.7	1887.7	
15	1866.7	-	2008.7	2008.7	
16	1907.7	1907.7	2049.7	2049.7	
17	1948.7	1948.7	2090.8	2090.7	
18	2028.7	-	2170.8	2170.7	
19	2069.7	-	2211.8	2211.8	
20	2110.8	2110.7	2252.8	2252.8	
21	2151.8	2151.7	2293.8	2293.8	

22	2272.8	-	2414.9	2414.9	
23	2313.8	2313.8	2455.9	2455.9	
24	2434.9	-	2576.9	2576.9	
25	2475.9	2475.8	2617.9	2617.9	

<sup>a</sup>The numbers are corresponding to the peak numbers labeled in Fig. 5a, b.

**Table S2.** List of identified N-linked glycans released from human serum using 4-HQ as a reactive matrix in the positive ion mode. The m/z values are corresponding to the ions of  $[M+(n+1)Na-nH]^+$  (n represents the amount of sialic acid residues in glycan).

Peak No. <sup>a</sup>	Theoretical m/z of native glycans	Theoretical m/z of 4-HQ-labeled glycans	Observed m/z of 4-HQ-labeled glycans	Structure
1	1136.4	1278.4	1278.4	
2	1257.4	1399.5	1399.5	
3	1298.4	1440.5	1440.5	
4	1339.5	1481.5	1481.5	
5	1419.5	1561.5	1561.5	
6	1460.5	1602.6	1602.6	
7	1485.5	1627.6	1627.6	
8	1501.5	1643.6	1643.6	
9	1542.6	1684.6	1684.6	
10	1581.5	1723.6	1723.6	
11	1622.6	1764.6	1764.6	

12	1647.6	1789.6	1789.6	
13	1663.6	1805.6	1805.6	
14	1688.6	1830.7	1830.7	
15	1704.6	1846.7	1846.7	
16	1743.6	1885.6	1885.6	
17	1809.6	1951.7	1951.7	
18	1814.6	1956.7	1956.7	
19	1850.7	1992.7	1992.7	
20	1866.7	2008.7	2008.7	
21	1905.6	2047.7	2047.7	
22	1976.7	2118.7	2118.7	
23	2012.7	2154.8	2154.8	
24	2017.7	2159.8	2159.8	
25	2122.7	2264.8	2264.8	
26	2179.8	2321.8	2321.8	
27	2289.8	2431.8	2431.8	
28	2325.8	2467.9	2467.9	
29	2435.8	2577.9	2577.9	
30	2638.9	2780.9	2780.9	

<sup>a</sup>The numbers are corresponding to the peak numbers labeled in Fig. S3a.

**Table S3.** List of identified N-linked glycans released from human serum using 4-HQ as a reactive matrix in the negative ion mode. The m/z values are corresponding to the ions of

$[M+(n-1)Na-nH]^-$  (n represents the amount of sialic acid residues in glycan).

Peak No. <sup>a</sup>	Theoretical m/z of native glycans	Theoretical m/z of 4-HQ-labeled glycans	Observed m/z of 4-HQ-labeled glycans	Structure
1	1565.6	1707.6	1707.6	
2	1727.6	1869.7	1869.7	
3	1768.6	1910.7	1910.7	
4	1889.7	2031.7	2031.7	
5	1914.7	2056.8	2056.8	
6	1930.7	2072.8	2072.8	
7	1971.7	2113.8	2113.8	
8	2076.7	2218.8	2218.8	
9	2117.8	2259.8	2259.8	
10	2133.8	2275.8	2275.8	
11	2186.8	2328.8	2328.8	
12	2243.8	2363.8 2385.8	2363.8 2385.8	
13	2279.8	2421.9	2421.9	
14	2295.8	2437.9	2437.9	
15	2389.8	2531.9	2531.9	
16	2592.9	2735.0	2735.0	
17	2922.0	3064.0	3064.0	
18	3068.1	3210.1	3210.1	

<sup>a</sup>The numbers are corresponding to the peak numbers labeled in Fig. S3b.