

An integrated strategy using LC-MS/MS combined with *in vivo* microdialysis for simultaneous determination of lignans of *Schisandra chinensis* (Turcz.) Bail.

Fructus and endogenous neurotransmitters: Application to pharmacokinetic and pharmacodynamic study

Jun Fu^{a,b}, Hongxu Zhang^{a,b}, Shu Liu^b, Jiajie Wu^{a,b}, Yuying Zhang^{a,b}, Yang Gao^{a,b}, Fengrui Song^b, Yuhua Qin^c, Xiuli Hu^a, Zhongying Liu^{a*}

Affiliation

^a *School of Pharmaceutical Sciences, Jilin University, Changchun 130021, China.*

^b *National Center of Mass Spectrometry in Changchun & Jilin Provincial Key Laboratory of Chinese Medicine Chemistry and Mass Spectrometry, Changchun Institute of Applied Chemistry, Chinese Academy of Sciences, Changchun, 130022, China.*

^c *Hainan Tropical Marine University, Sanya, 572022*

* Corresponding author: Zhongying Liu, School of Pharmaceutical Sciences, Jilin University, Changchun 130021, China.

E-mail address: liuzy@jlu.edu.cn

Table of contents

1. Table S1 in page 3.

Table S1 The *in vitro* recoveries for all analytes determined by gain method (n=3)

Comp	Con (ng/mL)	Blood prob					Brain prob				
		Flow rate ($\mu\text{L}/\text{min}$)					Flow rate ($\mu\text{L}/\text{min}$)				
		1.0	1.5	2.0	2.5	3.0	1.0	1.5	2.0	2.5	3.0
Gomisin A	10	31.52 \pm 2.20	30.82 \pm 8.81	27.22 \pm 3.21	23.60 \pm 3.84	19.75 \pm 6.77	30.35 \pm 3.78	27.79 \pm 1.33	25.49 \pm 4.01	18.83 \pm 2.07	12.69 \pm 5.38
	200	32.58 \pm 5.37	30.17 \pm 3.03	28.73 \pm 3.67	24.35 \pm 7.79	20.31 \pm 5.71	30.90 \pm 7.58	28.06 \pm 5.05	25.78 \pm 5.69	18.73 \pm 5.35	12.03 \pm 5.96
	400	32.21 \pm 5.23	30.41 \pm 5.69	27.47 \pm 6.93	23.73 \pm 8.86	20.15 \pm 8.80	30.82 \pm 7.83	28.17 \pm 3.84	24.55 \pm 5.75	19.10 \pm 5.82	12.42 \pm 5.01
Gomisin D	2	31.21 \pm 3.81	29.76 \pm 8.64	27.59 \pm 6.65	24.79 \pm 4.56	20.23 \pm 5.02	29.88 \pm 9.54	25.55 \pm 1.07	23.93 \pm 3.77	16.41 \pm 3.07	10.69 \pm 5.57
	200	32.33 \pm 5.65	31.82 \pm 1.61	27.66 \pm 4.48	25.26 \pm 5.87	21.76 \pm 3.84	30.18 \pm 7.68	25.19 \pm 2.27	23.44 \pm 5.65	15.26 \pm 1.47	11.76 \pm 1.20
	400	31.81 \pm 2.61	30.43 \pm 5.44	26.08 \pm 3.29	25.08 \pm 4.92	20.62 \pm 5.71	29.75 \pm 4.79	25.97 \pm 3.50	24.45 \pm 5.49	16.85 \pm 5.46	11.71 \pm 3.66
Gomisin E	2	34.55 \pm 2.04	30.26 \pm 2.61	25.84 \pm 5.10	24.18 \pm 8.04	20.13 \pm 5.42	31.79 \pm 7.90	27.57 \pm 3.88	24.75 \pm 4.07	17.71 \pm 3.92	11.69 \pm 1.99
	200	34.70 \pm 5.05	30.24 \pm 7.53	26.69 \pm 5.17	24.73 \pm 3.89	20.70 \pm 6.73	31.32 \pm 6.41	27.19 \pm 4.95	24.56 \pm 3.78	16.75 \pm 1.57	11.65 \pm 2.03
	400	33.21 \pm 6.50	30.33 \pm 1.74	26.35 \pm 4.47	23.96 \pm 3.61	19.95 \pm 7.12	31.24 \pm 3.84	28.01 \pm 5.17	24.40 \pm 4.78	18.15 \pm 4.38	12.16 \pm 5.52
Gomisin G	1	32.83 \pm 7.05	31.09 \pm 8.18	28.29 \pm 3.37	24.88 \pm 4.42	21.09 \pm 4.73	30.48 \pm 8.45	28.35 \pm 4.36	25.69 \pm 4.96	15.94 \pm 2.92	11.55 \pm 3.82
	200	33.28 \pm 2.99	31.56 \pm 5.22	28.25 \pm 4.77	25.11 \pm 8.93	21.87 \pm 5.32	30.17 \pm 5.39	28.51 \pm 5.87	25.29 \pm 5.96	15.83 \pm 4.46	11.51 \pm 1.31
	400	33.70 \pm 6.61	30.28 \pm 9.64	28.57 \pm 5.46	25.26 \pm 4.81	21.32 \pm 8.57	31.01 \pm 7.17	27.58 \pm 2.98	25.38 \pm 4.46	16.67 \pm 2.30	12.29 \pm 2.16
Gomisin K1	1	33.74 \pm 2.90	30.15 \pm 3.39	27.05 \pm 6.85	24.27 \pm 6.60	21.09 \pm 4.49	32.89 \pm 3.68	28.16 \pm 2.30	23.86 \pm 2.01	19.34 \pm 2.40	11.83 \pm 3.99
	200	34.51 \pm 8.01	30.07 \pm 9.10	27.14 \pm 3.31	23.79 \pm 6.56	21.51 \pm 7.50	33.58 \pm 4.79	28.45 \pm 1.55	24.43 \pm 4.38	19.90 \pm 4.03	10.56 \pm 3.01
	400	34.20 \pm 5.99	30.21 \pm 3.46	27.43 \pm 3.81	23.63 \pm 5.58	20.82 \pm 6.44	33.18 \pm 7.90	27.81 \pm 2.08	23.47 \pm 4.37	18.84 \pm 3.14	12.02 \pm 3.21
Benzoyl-gomisin-H	0.5	34.35 \pm 5.45	30.55 \pm 4.33	27.66 \pm 4.63	22.46 \pm 7.80	20.42 \pm 8.31	32.95 \pm 6.41	25.11 \pm 2.88	21.19 \pm 4.60	15.86 \pm 3.65	11.73 \pm 5.21
	200	34.21 \pm 3.37	30.42 \pm 5.00	27.43 \pm 5.32	23.61 \pm 6.78	21.77 \pm 4.11	33.09 \pm 3.84	25.23 \pm 2.67	21.15 \pm 5.47	15.21 \pm 1.70	11.06 \pm 1.84
	400	35.09 \pm 7.85	31.38 \pm 6.49	27.54 \pm 5.18	23.59 \pm 5.36	21.10 \pm 4.68	33.25 \pm 8.45	24.75 \pm 3.12	20.53 \pm 1.53	15.70 \pm 5.26	10.50 \pm 1.77
Angeloylgomisin Q	2	33.32 \pm 3.31	30.74 \pm 1.42	26.78 \pm 4.75	24.13 \pm 7.70	20.98 \pm 3.27	29.92 \pm 5.39	27.54 \pm 1.51	24.12 \pm 3.61	18.88 \pm 2.10	12.48 \pm 4.81
	200	33.75 \pm 6.45	30.97 \pm 2.44	26.96 \pm 8.77	24.21 \pm 5.78	21.18 \pm 5.72	30.57 \pm 4.79	28.78 \pm 1.15	23.25 \pm 3.28	19.31 \pm 2.23	11.71 \pm 5.47
	400	32.64 \pm 8.60	30.61 \pm 3.86	27.19 \pm 3.57	23.88 \pm 6.66	21.52 \pm 6.65	30.20 \pm 4.38	28.86 \pm 3.91	23.35 \pm 3.24	19.42 \pm 2.18	12.22 \pm 4.83

Schizandrin	2	33.97±7.17	30.21±7.74	26.63±4.71	23.62±4.23	21.39±5.24	31.13±6.41	26.26±4.57	21.69±1.67	18.91±3.43	11.37±4.61
	200	35.12±2.17	30.32±3.06	27.93±4.04	24.49±4.65	21.30±3.70	31.19±3.84	25.75±1.21	20.98±3.76	17.63±2.19	11.79±1.89
	400	34.31±2.56	31.16±1.77	27.43±4.14	24.12±6.71	21.01±7.92	31.55±3.68	27.18±2.03	21.51±5.98	18.75±3.42	11.42±1.65
Schisandrin A	1	33.41±9.16	30.45±3.02	28.85±3.56	23.45±7.46	20.65±7.89	31.39±4.07	25.37±5.42	21.86±3.41	17.63±4.58	10.74±5.88
	200	34.36±4.96	30.70±2.84	28.38±6.53	23.22±7.48	21.27±4.62	31.71±6.20	25.47±5.12	22.26±4.75	17.65±5.92	10.77±2.97
	400	34.18±4.03	30.09±8.35	27.96±6.92	23.87±4.83	21.15±3.39	31.82±6.75	25.14±2.50	22.29±5.92	18.51±3.86	11.26±2.89
Schisandrin B	1	32.34±5.99	29.73±3.18	27.19±3.42	23.19±6.57	20.71±4.42	29.92±4.24	25.66±3.58	22.13±3.87	19.57±5.43	11.68±4.50
	200	32.49±4.14	30.15±4.14	27.02±3.25	23.27±5.72	21.43±7.42	30.20±3.92	25.85±3.26	21.52±5.97	19.55±3.56	11.47±4.51
	400	32.14±6.41	30.25±6.96	26.86±5.89	22.64±6.74	21.39±5.37	30.31±7.08	25.63±4.49	21.99±3.37	19.78±2.94	11.91±2.71
NE	50	34.18±4.75	30.59±3.64	27.48±3.90	25.36±7.29	20.74±4.72	32.64±3.22	27.51±1.62	22.53±1.01	18.25±2.84	12.14±1.69
	400	34.27±2.33	31.13±3.68	27.90±6.38	24.74±7.94	21.65±5.74	32.79±5.47	27.45±1.53	22.97±1.73	17.62±3.40	12.05±5.82
	800	33.91±2.73	30.83±3.93	26.92±4.49	25.18±6.83	21.17±6.24	31.96±3.83	27.15±1.28	23.22±5.80	18.19±3.77	11.61±3.73
5-HT	2	35.79±5.38	31.97±5.57	29.37±4.37	25.81±4.38	19.54±4.35	32.77±5.47	27.73±2.66	22.79±1.80	14.66±3.51	11.03±1.92
	200	35.93±4.95	32.07±3.06	29.30±4.35	25.64±5.04	19.04±4.42	33.12±5.28	27.91±1.78	22.43±4.36	15.23±4.36	10.41±4.60
	400	35.57±1.80	32.12±9.43	28.94±8.73	25.59±5.05	19.97±5.31	33.19±5.96	28.07±1.98	21.87±4.43	14.71±2.13	10.84±1.26
Ach	20	34.92±7.90	30.32±5.94	25.95±7.82	22.24±4.96	19.32±5.68	31.68±7.35	27.65±2.29	21.48±2.38	18.54±4.48	10.78±2.96
	400	35.33±1.48	30.64±6.64	25.68±5.92	22.49±5.64	18.52±7.72	32.13±3.92	28.12±4.15	21.55±2.14	18.32±3.72	10.51±2.83
	800	35.21±3.14	30.61±1.58	26.17±7.89	21.65±6.07	19.28±5.92	32.26±3.66	28.24±3.07	21.16±4.58	19.22±5.73	11.22±5.24
GABA	5	34.95±5.04	31.02±5.45	28.39±8.48	23.57±7.54	19.58±5.57	30.54±5.30	27.07±4.45	23.95±2.18	17.59±1.29	11.77±5.32
	200	35.06±3.70	30.41±1.30	28.61±7.69	23.42±4.20	19.59±7.30	30.72±4.22	26.63±2.54	23.96±3.63	17.57±3.56	12.52±3.07
	400	34.73±6.52	30.92±1.54	28.41±7.14	23.75±4.99	19.84±4.76	29.62±2.80	26.96±2.13	23.77±2.19	18.22±3.18	12.17±1.01
Asp	20	33.66±7.82	30.45±3.92	25.95±5.81	22.74±3.84	20.31±8.96	31.53±6.65	27.46±1.85	20.48±3.79	16.06±4.17	12.96±4.88
	400	33.81±6.19	30.04±2.62	26.38±8.61	23.27±3.38	20.84±6.47	31.81±5.92	27.09±1.42	20.46±2.04	16.13±5.61	12.91±3.82
	800	33.72±7.71	30.56±4.21	26.29±3.47	23.26±8.41	21.24±6.35	31.69±5.20	27.18±1.75	20.93±2.47	16.29±3.32	12.32±1.88
GLU	20	35.68±5.64	30.81±4.94	27.82±3.27	23.85±4.81	20.72±8.84	33.63±4.49	28.16±3.21	24.94±2.82	19.69±1.27	12.56±3.12

	400	35.40±9.45	31.19±7.79	27.56±8.82	23.61±6.26	20.35±5.04	34.82±8.03	28.02±4.77	25.09±4.66	18.94±3.70	12.41±2.82
	800	35.97±6.31	31.52±4.43	27.53±4.45	23.94±7.05	20.58±8.52	34.58±4.01	27.91±2.11	25.26±2.42	19.76±3.68	12.88±2.33
GLY	50	34.94±5.55	31.57±9.46	27.65±4.03	24.64±8.88	19.72±3.57	32.94±3.35	28.16±3.41	23.58±5.45	15.91±1.09	12.69±3.01
	400	34.71±5.68	31.96±6.63	27.71±6.22	24.50±5.11	19.90±8.46	32.97±3.98	28.57±3.38	23.97±2.55	16.24±3.95	12.62±4.94
	800	34.87±6.08	31.79±4.71	26.82±4.28	24.71±3.99	19.61±6.77	32.53±3.51	27.91±5.26	23.83±3.18	16.22±3.66	12.55±1.33
