

Supplementary Information

An amino-substituted 2-(2'-hydroxyphenyl)benzimidazole for the fluorescent detection of phosgene based on ESIPT mechanism

Zi-Jie Li,^a Wen-Jie Zhang,^a Wen-Zhu Bi,^{*a} Qiu-Juan Ma,^{*a} Su-Xiang Feng,^{*abc} Xiao-Lan Chen,^d and Ling-Bo Qu^d

^a School of Pharmacy, Henan University of Chinese Medicine, Zhengzhou, 450046, China.

^b Collaborative Innovation Center for Respiratory Disease Diagnosis and Treatment & Chinese Medicine Development of Henan Province, Zhengzhou, 450046, China

^c Zhengzhou Key Laboratory of Chinese Medicine Quality Control and Evaluation, Zhengzhou, 450046, China

^d College of Chemistry, Zhengzhou University, Zhengzhou, 450052, China.

Corresponding E-mail: biwenzhu2018@hactcm.edu.cn

Contents

1. Evaluation of HBI as the fluorescent probe for the detection of phosgene.....	S2
2. Measurement of the fluorescence quantum yields.....	S3
3. UV-Vis spectra of P1 upon addition of phosgene.....	S4
4. Measurement of the LoD for P1.....	S4
5. The selectivity of P1 for the detection of phosgene.....	S5
6. Fluorescence spectra of P1 in the presence of interfering compounds.....	S6
7. Table S1 Determination of phosgene in the presence of interfering compounds..	S6
8. HPLC-MS chromatogram and spectrum of the reaction mixture.....	S6
9. ¹ H NMR, ¹³ C NMR and HRMS copies of P1 and P1-CO.....	S7
10. Table S2 Comparison of recently reported probes for phosgene detection.....	S11

1. Evaluation of HBI as the fluorescence probe for the detection of phosgene

HBI (2-(2'-hydroxyphenyl)benzimidazole) was synthesized according to literature methods with minor modifications.¹ The method for the evaluation of HBI as the fluorescent

probe for the detection of phosgene was same as that of **P1**.




Figure S1 (A) Time-dependent fluorescence intensity at 351 nm of HBI (10 μ M) treated with triphosgene (10 μ M) in the presence (black line) or absence (red line) of TEA (20 μ L, 0.1 vol%), λ_{ex} = 309 nm, λ_{em} = 351 nm, slit width = 5/5 nm. (B) Fluorescence spectra of 10 μ M HBI solutions containing TEA (20 μ L, 0.1 vol%) upon addition of triphosgene (0-5 μ M), λ_{ex} = 309 nm, slit width = 5/5 nm. (C) Fluorescence spectra of HBI solutions before and after addition of phosgene or various analytes (0.2 mM) for 10 min: (1) blank, (2) triphosgene (5 μ M)/TEA (20 μ L, 0.1 vol%), (3) (COCl)₂, (4) CH₃COCl, (5) SOCl₂, (6) TsCl, (7) DCP, (8) HOAc, (9) POCl₃, (10) SO₂Cl₂, (11) triphosgene. λ_{ex} = 309 nm. (D) Images of HBI solution before (left) and after (right) addition of phosgene under 365 nm UV-light.

Reference

1. K. Akutsu, S. Mori, K. Shinmei, H. Iwase, Y. Nakano and Y. Fujii, *Talanta*, 2016, **146**, 575.

2. Measurement of the fluorescence quantum yields




Figure S2 Measurement of the fluorescence quantum yields (Φ_f) of **P1** and **P1-CO**. **P1** and **P1-CO** were determined in CH_2Cl_2 with solvent refractive index correction. Quinine sulfate in 1.0 M H_2SO_4 was used as the reference ($\Phi_f = 54\%$) at an excitation wavelength of 346 nm. The fluorescence quantum yield was calculated by the following equation: $\Phi_x = \Phi_s (F_x/F_s)(A_s/A_x)(n_x/n_s)^2$. Where x and s indicate the determined and reference, respectively, F is the area of the fluorescence peak, A is the optical density at the excitation wavelength, and n is the refractive index of the solvent.

3. UV-Vis spectra of P1 upon addition of phosgene




Figure S3. UV-Vis spectra of 10 μM **P1** solutions containing TEA (20 μL , 0.1 vol%) upon addition of triphosgene (0-6 μM).

4. Measurement of the LoD for P1




Figure S4 Measurement of the LoD for **P1** to triphosgene. (A) The ratio of the emission intensities at 358 and 540 nm vs. the triphosgene concentration. (B) Ten times of the blank experiment to evaluate the standard deviation σ . The triphosgene detection limit was determined to be 5.3 nM ($\text{LoD} = 3\sigma/k$, where σ is the standard deviation of the blank experiment, and k is the slope of the relationship between the emission-intensity ratio and the phosgene concentration).

5. The selectivity of P1 for the detection of phosgene.




Figure S5. Fluorescence spectra of **P1** (10 μM) containing TEA (20 μL , 0.1 vol%) in 2 mL CH_2Cl_2 upon addition of phosgene or various analytes (50 μM) for 10 min. λ_{ex} = 305 nm.

6. Fluorescence spectra of **P1** in the presence of interfering compounds




Figure S6 Fluorescence spectra of 10 μM **P1** containing TEA (20 μL , 0.1 vol%) and interferents (A, 5.0 μM ; B, 10 μM) in the presence of 3.0 μM of triphosgene in 2 mL CH_2Cl_2 . λ_{ex} = 305 nm.

7. Table S1 Determination of phosgene in the presence of interfering compounds

Interfering compounds (5.0 μ M)	Phosgene added (μ M)	Phosgene found (μ M)	Recovery (%)	Interfering compounds (10 μ M)	Phosgene added (μ M)	Phosgene found (μ M)	Recovery (%)
(COCl) ₂	9.0	8.7	96.7	(COCl) ₂	9.0	5.15	57.1
CH ₃ COCl	9.0	9.8	108.9	CH ₃ COCl	9.0	4.7	52.0
SOCl ₂	9.0	8.1	90.0	SOCl ₂	9.0	4.55	50.8
TsCl	9.0	8.15	90.5	TsCl	9.0	4.75	52.8
DCP	9.0	9.45	105.2	DCP	9.0	5.6	62.4
HOAc	9.0	9.35	104.1	HOAc	9.0	2.45	27.5
POCl ₃	9.0	7.9	87.8	POCl ₃	9.0	4.2	46.6
SO ₂ Cl ₂	9.0	8.15	90.7	SO ₂ Cl ₂	9.0	1.85	20.6

8. HPLC-MS chromatogram and spectrum of the reaction mixture




Figure S7 HPLC chromatogram (up) of the reaction mixture and MS spectrum (down) of the peak at 21.63 min. Chromatographic conditions: chromatographic column: Accucore™ C18 (150 × 4.6 mm, 2.6 μ m); mobile phase: 0.1% formic acid (A), 0.1% formic acid water

acetonitrile (B); gradient elution: 0 ~ 1 min, 95% A; 1~25 min, 95%~0% A; 25~29 min, 0% A; 29~29.1 min, 0%~95% A, 29.1~34 min, 95% A; flow rate: 0.2 mL/min; column temperature: 30 °C. Mass spectrometry conditions: ion source: heated electrospray ionization source (H-ESI); auxiliary gas pressure: 1.0 Mpa; capillary temperature: 300 °C, auxiliary heating temperature: 275 °C.

9. ¹H NMR, ¹³C NMR and HRMS copies of P1 and P1-CO




Figure S8 ¹H NMR of P1




Figure S9 ¹³C NMR of P1




Figure S10 HRMS of P1




Figure S11 ^1H NMR of P1-CO




Figure S12 ^{13}C NMR of P1-CO


























Figure S13 HRMS of P1-CO







10. Table S2 Comparison of recently reported probes for phosgene detection.







Structure	Mechanism	Detection type	Reaction process	$\lambda_{ex}/\lambda_{em}$ (nm)	Detection limit	Response time	Reference
	FRET	Ratiometric	intermolecular reaction of two fluorophores	$\lambda_{ex} = 343$ nm $\lambda_{em} = 464$ nm	50 μ M	--	<i>Chem. Commun.</i> , 2007, 12 , 1238.
	--	Turn-on	ring opening reaction of benzimidazole-fused rhodamine dye	$\lambda_{ex} = 560$ nm $\lambda_{em} = 590$ nm	50 nM (Triphosgene)	20-30 s	<i>Chem Commun.</i> , 2012, 48 , 1895.
	--	Turn-on	intramolecular reaction of cinnamic acids	$\lambda_{ex} = 330$ nm $\lambda_{em} = 382$ nm	1 nM	--	<i>Anal. Chem.</i> , 2012, 84 , 4594
	PET	Turn-on	twice carbamylation reactions	$\lambda_{ex} = 270$ nm $\lambda_{em} = 308$ nm	0.7 ppb (Triphosgene)	2 min	<i>ACS Appl. Mater. Interfaces.</i> , 2016, 8 , 22246.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{ex} = 580$ nm $\lambda_{em} = 593$ nm	20 nM (Triphosgene)	2 min	<i>Angew. Chem. Int. Edit.</i> , 2016, 55 , 4729.




	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 410 \text{ nm}$ $\lambda_{\text{em}} = 511/442 \text{ nm}$	1.3 nM (Triphosgene)	20 min	<i>Chem. Commun.</i> , 2017, 53 , 1530.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 368 \text{ nm}$ $\lambda_{\text{em}} = 446 \text{ nm}$	3 nM (Triphosgene)	0.5 min	<i>ACS Sens.</i> , 2017, 2 , 178.
	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 390/465 \text{ nm}$ $\lambda_{\text{em}} = 445/512 \text{ nm}$	0.12 nM (Phosgene)	1.5 s	<i>ACS Appl. Mater. Interfaces</i> , 2017, 9 , 13920.
	--	Turn-on	Dehydration of oxime to nitrile	$\lambda_{\text{ex}} = 530 \text{ nm}$ $\lambda_{\text{em}} = 570 \text{ nm}$	0.09 ppb (Triphosgene)	10 s	<i>Anal. Chem.</i> , 2017, 89 , 12837.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 450 \text{ nm}$ $\lambda_{\text{em}} = 530 \text{ nm}$	2.7 nM (Triphosgene)	15 s	<i>Anal. Chem.</i> , 2017, 89 , 4192.

	ESIPT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 375 \text{ nm}$ $\lambda_{\text{em}} = 445/495 \text{ nm}$	0.14 ppm (Phosgene)	5 min	<i>Anal. Chem.</i> , 2017, 89 , 12596.
	ICT	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 380 \text{ nm}$ $\lambda_{\text{em}} = 508 \text{ nm}$	20 nM (Phosgene)	20 min	<i>Org. Chem. Front.</i> , 2017, 4 , 1719.
	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 440 \text{ nm}$ $\lambda_{\text{em}} = 482/550 \text{ nm}$	27 nM (Phosgene)	2 min	<i>Anal. Chim. Acta.</i> , 2018, 1029 , 97.
	ICT	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 400 \text{ nm}$ $\lambda_{\text{em}} = 468 \text{ nm}$	0.2 nM (Triphosgene)	30 s	<i>Chem. Eur. J.</i> , 2018, 24 , 5652.
	--	Turn-on	conversion of amide to nitrile	$\lambda_{\text{ex}} = 653 \text{ nm}$ $\lambda_{\text{em}} = 679 \text{ nm}$	8.9 nM (Triphosgene)	4 min	<i>J. Mater. Chem. C.</i> , 2018, 6 , 10472.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 460 \text{ nm}$ $\lambda_{\text{em}} = 511 \text{ nm}$	179 nM (Triphosgene)	10 s	<i>Chem. Eur. J.</i> , 2018, 24 , 3136.

	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{ex} = 434/502 \text{ nm}$ $\lambda_{em} = 482/615 \text{ nm}$	2.3 nM (Phosgene)	5 min	<i>Anal. Chem.</i> , 2018, 90 , 8686.
	--	Turn-on	Spirocyclic ring -open reaction	$\lambda_{ex} = 530 \text{ nm}$ $\lambda_{em} = 578 \text{ nm}$	3.2 ppb (Triphosgene)	2 min	<i>Anal. Chem.</i> , 2018, 90 , 3382.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{ex} = 480 \text{ nm}$ $\lambda_{em} = 516 \text{ nm}$	24 pM (Phosgene)	within 3 s	<i>Sens. Actuator B-Chem.</i> , 2019, 283 , 458.
	ESIPT	Turn-on	conversion of oxime to nitrile	$\lambda_{ex} = 438 \text{ nm}$ $\lambda_{em} = 474 \text{ nm}$	0.48 nM (Phosgene)	20 min	<i>Dyes Pigment.</i> , 2019, 163 , 483.
	ESIPT	Ratiometric	conversion of oxime to isoxazole	$\lambda_{ex} = 382 \text{ nm}$ $\lambda_{em} = 495/577 \text{ nm}$	0.087 ppm (Phosgene)	1.43 s	<i>J. Mater. Chem. A.</i> , 2019, 7 , 1756.
	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{ex} = 470 \text{ nm}$ $\lambda_{em} = 520/610 \text{ nm}$	0.09 nM (Phosgene)	Within 20 s	<i>Dyes Pigment.</i> , 2019, 163 , 489.

	ICT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 400 \text{ nm}$ $\lambda_{\text{em}} = 488/548 \text{ nm}$	0.3 nM (Phosgene)	60 s	<i>J. Mater. Chem. C.</i> , 2019, 7 , 1510.
	ESIPT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 335 \text{ nm}$ $\lambda_{\text{em}} = 393/469 \text{ nm}$	0.14 ppm (Phosgene)	30 s	<i>Talanta</i> , 2019, 200 , 78.
	ICT	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 390 \text{ nm}$ $\lambda_{\text{em}} = 422/526 \text{ nm}$	3.2 nM (Phosgene)	within 10 s	<i>Anal Chem.</i> , 2019, 91 , 5690.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 460 \text{ nm}$ $\lambda_{\text{em}} = 525 \text{ nm}$	1.2 nM (Triphosgene)	within 20 s	<i>Anal. Methods</i> , 2019, 11 , 4600.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 350 \text{ nm}$ $\lambda_{\text{em}} = 430 \text{ nm}$	0.4 μM (Triphosgene)	20 s	<i>New J. Chem.</i> , 2019, 43 , 11743.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 440 \text{ nm}$ $\lambda_{\text{em}} = 500 \text{ nm}$	72 nM (Phosgene)	2 min	<i>Anal Chem.</i> , 2019, 91 , 12070.

	PET	Ratiometric	chloroformylation reaction of N atom	$\lambda_{\text{ex}} = 400 \text{ nm}$ $\lambda_{\text{em}} = 422/445/477 \text{ nm}$	1.54 nM (Phosgene)	within 50 s	<i>New J. Chem.</i> , 2019, 43 , 14991.
	--	Turn-on	Conversion of amide to nitrile	$\lambda_{\text{ex}} = 342 \text{ nm}$ $\lambda_{\text{em}} = 440 \text{ nm}$	5.56 nM (Phosgene)	1.4 min	<i>ChemistrySelect</i> , 2019, 4 , 2968.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 432 \text{ nm}$ $\lambda_{\text{em}} = 484 \text{ nm}$	4.6 nM (Phosgene)	15 s	<i>Chem. Commun.</i> , 2019, 55 , 13753.
	PET	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 467 \text{ nm}$ $\lambda_{\text{em}} = 530 \text{ nm}$	0.15 nM (Phosgene)	1.5 s	<i>Org. Lett.</i> , 2019, 21 , 9497.
	ICT	Turn-on	Beckmann rearrangement of ketoxime	$\lambda_{\text{ex}} = 367 \text{ nm}$ $\lambda_{\text{em}} = 448 \text{ nm}$	6.3 nM (Phosgene)	15 min	<i>Dyes Pigment.</i> , 2020, 173 , 107854.
	--	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 440 \text{ nm}$ $\lambda_{\text{em}} = 492 \text{ nm}$	84.2 nM (Phosgene)	2 min	<i>Dyes Pigment.</i> , 2020, 173 , 107933.

	--	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 370 \text{ nm}$ $\lambda_{\text{em}} = 416 \text{ nm}$	1.27 nM (Phosgene)	1 min	<i>Talanta</i> , 2021, 221 , 121477.
	ICT	Turn-on	twice carbamylation reactions	$\lambda_{\text{ex}} = 390 \text{ nm}$ $\lambda_{\text{em}} = 442/483/517 \text{ nm}$	1.65 nM (Phosgene)	3 min	<i>Sens. Actuator B-Chem.</i> , 2021, 326 , 128837.
 <p>This work</p>	ESIPT	Ratiometric	twice carbamylation reactions	$\lambda_{\text{ex}} = 305 \text{ nm}$ $\lambda_{\text{em}} = 358/540 \text{ nm}$	5.3 nM (Trihosgene)	Within 50 s	