

Supporting Information

Enhancing the Performance of Photoelectrochemical Glucose Sensor via the Electron Cloud Bridge of Au in SrTiO₃/PDA Electrodes

Yadong Wang^a, Jinxin Ma^a, Nan Zhang^a, Delun Chen^a, Jinchun Tu^a, Yang Cao^{a, c},
Qiang Wu^{b, *}, Xiaolin Zhang^{c, *}, Wanjun Hao^{a, *}

^a *Key Laboratory of Advanced Materials of Tropical Island Resources, State Key Laboratory of Marine Resource Utilization in South China Sea, Hainan University, Haikou 570228, P. R. China;*

^b *School of Tropical Medicine and Laboratory Medicine, Key Laboratory of Emergency and Trauma of Ministry of Education, Research Unit of Island Emergency Medicine, Chinese Academy of Medical Sciences (No. 2019RU013), Hainan Medical University, Haikou, 571199, China;*

^c *Key Laboratory of Child Cognition & Behavior Development of Hainan Province, Qiongtai Normal University, Haikou 571127, P. R. China*

Glucose oxidase (GOx) is specific catalysis for glucose, so the electrode can only work in the presence of glucose. Then, the performance toward selectivity of photoelectrochemistry (PEC) electrode is well-defined. The catalytic process of GOx and the photocurrent response can be described as follows [1]:







Fig. S1. EDX spectrum of the 3D hollowed out SrTiO₃/Au/PDA nanoarray.




Fig. S2. HERTEM image of the 3D hollowed out SrTiO₃/PDA nanoarray.




Fig. S3. XRD pattern of the precursor 3D hollowed out TiO₂.




Fig. S4. The high-resolution spectra of C1s: (a) SrTiO₃/Au; (b) SrTiO₃/Au/PDA.




Fig. S5. Photocurrent response.




Fig. S6. The 3D hollowed out SrTiO₃/PDA/GOx PEC sensor. (a) Photocurrent responses at different glucose concentrations (0–20 mM) in 0.1 M PBS (pH = 7.4). (b) Calibration curve of glucose concentration and photocurrent density.




Fig. S7. (a) TEM of SrTiO₃/Au/PDA. (b) TEM of SrTiO₃/Au/PDA (excess PDA).




Fig. S8. The 3D hollowed out SrTiO₃/Au/GOx PEC sensor. (a) Photocurrent responses at different glucose concentrations (0–20 mM) in 0.1 M PBS (pH = 7.4). (b) Calibration curve of glucose concentration and photocurrent density. The 3D hollowed out SrTiO₃/Au/PDA (excess PDA)/GOx PEC sensor. (c) Photocurrent responses at different glucose concentrations (0–20 mM) in 0.1 M PBS (pH = 7.4). (d) Calibration curve of glucose concentration and photocurrent density.




Fig. S9. Long-term (10 days) stability test of the hollowed out SrTiO₃/Au/PDA/GOx PEC sensor in 0.1 M PBS (pH = 7.4) (presence of 20 mM glucose).

Table S1. The surface chemical states of C element and corresponding percentages in the samples.

Peak	C-O or C-N (%)
SrTiO ₃ /Au	28.57
SrTiO ₃ /Au/PDA	62.60

Table S2. Determination of glucose concentration in human serum samples (n = 5).

Sample	Known/mM	Discovered/mM	Added glucose/mM	Total discovered/mM	RSD/%	Recovery/%
1	0.840	0.851	0.500	1.351	3.200	100.831
2	6.080	6.133	0.500	6.633	3.100	100.805
3	7.410	7.396	0.500	7.896	2.900	99.823
4	11.370	11.190	0.500	11.690	3.100	98.487
5	16.940	16.679	0.500	17.179	3.200	98.507

References

- [1] B.D. Yan, X.R. Zhao, D.L. Chen, Y. Cao, C.Z. Lv, J.C. Tu, X.H. Wang, Q. Wu, Enhanced photoelectrochemical biosensing performance for Au nanoparticle-polyaniline-TiO₂ heterojunction composites, *RSC Adv.* 10(72) (2020) 43985-43993.