

SUPPORTING INFORMATION

Rapid diagnosis of egg allergy by targeting ovalbumin specific IgE and IgG4 in serum at a disposable electrochemical immunoplatform

José M. Gordón Pidal^{δ,1,2}, Alejandro Valverde^{δ,1}, Sara Benedé^{3,*}, Elena Molina³, María Moreno-Guzmán⁴, Miguel Ángel López^{2,5}, José M. Pingarrón¹, Alberto Escarpa^{2,5,*}, Susana Campuzano^{1,*}

^δThese authors contributed equally.

¹ Department of Analytical Chemistry, Faculty of Chemistry, Complutense University of Madrid, 28040 Madrid, Spain

² Department of Analytical Chemistry, Physical Chemistry and Chemical Engineering, University of Alcalá, Alcalá de Henares E-28871, Madrid, Spain

³ Institute of Food Science Research (CIAL), C/ Nicolás Cabrera 9, 28049, Madrid, Spain

⁴ Department of Chemistry in Pharmaceutical Sciences, Analytical Chemistry, Faculty of Pharmacy, Complutense University of Madrid, Av. Complutense s/n, 28040 Madrid, Spain

⁵ Chemical Research Institute "Andrés M. del Río", University of Alcalá, Alcalá de Henares E-28871, Madrid, Spain

*To whom correspondence should be addressed: s.benede@csic.es; alberto.escarpa@uah.es;
susanacr@quim.ucm.es

CONTENTS	PAGE
Fig. S1	S3
Fig. S2	S4
Fig. S3	S5
Table S1	S5
Table S2	S6
Fig. S4	S7
Fig. S5	S7
Fig. S6	S8




Fig. S1 Optimization of key experimental variables involved in the preparation and functioning of the bioplatform prepared for the amperometric determination of OVA-specific IgE. Dependence of the amperometric responses measured in the absence (empty bars, B) and in the presence of 0.30 kU L^{-1} (full bars, S) IgE serum standard solutions, and the corresponding S/B ratios (in red), with OVA concentration and incubation time (a, b); number of steps involved in the preparation of the Av-HRP-b-anti-IgE-IgE-OVA-MBs (c); b-anti-IgE dilution (d); incubation time with IgE serum standard and b-anti-IgE mixture solutions (e); and Av-HRP concentration and incubation time (f, g).




Fig. S2 Optimization of key experimental variables involved in the preparation and functioning of the bioplatform prepared for the amperometric determination of OVA-specific IgG4. Dependence of the amperometric responses measured in the absence (empty bars, B) and in the presence of $0.036 \mu\text{g L}^{-1}$ (full bars, S) IgG4 serum standard solutions, and the corresponding S/B ratios (in red), with OVA concentration and incubation time (a, b); number of steps involved in the preparation of Av-HRP-b-anti-IgG4-IgG4-OVA-MBs (c); b-anti-IgG4 and Av-HRP concentrations (d, e), and incubation time with IgG4 serum standard, b-anti-IgG4 and Av-HRP mixture solution (f).




Fig. S3 Amperometric responses (and corresponding S/B ratios) measured with the prepared bioplatforms for the single determination of IgE (a) or IgG4 (b) in the absence (empty bars, B) and in the presence of 0.30 kU L^{-1} IgE or $0.036 \mu\text{g L}^{-1}$ IgG4 (full bars, S) serum standard solutions using enzyme labelling of the magnetic bioconjugates with Strep-HRP or Av-HRP.

Table S1 Optimization of key experimental variables involved in the single amperometric determination of OVA-specific IgE and IgG4 with the developed bioplatforms

Experimental variable	Selected value	
	IgE	IgG4
[OVA], $\mu\text{g mL}^{-1}$	50	50
Incubation time with OVA, min	30	30
Number of incubation steps	2	1
b-anti-IgE, dil.	1/1,000	--
Incubation time with (IgE serum standard + b-anti-IgE) mixture solution, min	60	--
Av-HRP, dil.	1/5,000	--
Incubation time with Av-HRP, min	30	--
b-anti-IgG4, dil.	--	1/1,000
Av-HRP, dil.	--	1/5,000
Incubation time with (IgG4 standard + b-anti-IgG4 + Av-HRP) mixture solution, min	--	60

Table S2 Comparison of the available methodologies to determine IgE and/or IgG4 in serum

Basis	Detection type	Target Igs	LOD	Assay time	Ref.*
ELISA	Absorbance	IgG4	--	4 h	[22]
LICA based on nanomicrospheres	Chemiluminescence	Children's milk protein-specific-IgE	--	3 h	[23]
ImmunoCAP	Fluorescence	Hen's egg white-specific IgE	0.1 kU L ⁻¹ IgE (human serum); 0.0007 µg L ⁻¹ IgG4 (plasma)	1 h 40 min	[24,25]
Microarray-based method	Fluorescence	Egg-white-specific IgE and IgG4	0.25 kU L ⁻¹ IgE; 100 µg L ⁻¹ IgG4	21 h 30 min	[42]
ELISA methodologies involving the same immunoreagents as the bioplatforms	Absorbance	OVA specific IgE and IgG4	0.03 kU L ⁻¹ IgE; 0.006 µg L ⁻¹ IgG4	4 h	This work
Immunoassays performed on the surface of OVA-MBs	Electrochemical	OVA specific IgE and IgG4	LOD: 0.003 kU L ⁻¹ IgE; 0.0002 µg L ⁻¹ IgG4	1 h 30 min	This work

LICA: light-initiated chemiluminescent assay; MBs: magnetic beads; OVA: ovalbumin

* These references correspond to those given in the manuscript




Fig. S4 Amperometric responses provided with the bioplatforms prepared from the stored OVA-MBs (at 4 °C in filtered PBS pH 7.4) in the absence (empty circles) and in the presence of 0.30 kU L⁻¹ IgE (a) and 0.036 µg L⁻¹ IgG4 (b) serum standard solutions (full circles). Control limits (red dashed lines) were calculated as $\pm 3 \times$ mean values of three amperometric responses obtained the day of OVA-MBs preparation (day 0).




Fig. S5 Amperometric responses (and corresponding S/B ratios) provided by the developed bioplatforms for the determination of OVA-specific IgE (a) and IgG4 (b) for 0 (white bars, B), and 0.30 kU L⁻¹ IgE or 0.036 µg L⁻¹ IgG4 (full bars, S) serum standard solutions prepared in the absence and in the presence of 50 mg mL⁻¹ HSA and 1 mg mL⁻¹ HIgG undiluted or 50 and 500 times diluted.




Fig. S6 Correlation plots for the concentrations of OVA-specific IgE (a) and IgG4 (b) in the analysed serum samples (replicates included in the plot) obtained with the developed electrochemical bioplatforms and the ELISA method. Red lines: regression lines with r = Pearson regression coefficient.