

Supporting information

Gold-based paper for antigen detection of monkeypox virus

Liya Ye^{a,b}, Xianlu Lei^{a,b}, Xinxin Xu^{a,b*}, Liguang Xu^{a,b}, Hua Kuang^{a,b}, Chuanlai Xu^{a,b}

^aState Key Laboratory of Food Science and Technology, Jiangnan University, Wuxi, People's Republic of China;

^bInternational Joint Research Laboratory for Biointerface and Biodetection, and School of Food Science and Technology, Jiangnan University, Wuxi, People's Republic of China;

*Corresponding author. Email: xxx89826@163.com;

Contents

Fig. S1. Soluble expression of recombinant A29 protein. (A): M: The molecular weight standard of protein. 1: Supernatant induced at 37°C. 2: Precipitation induced at 37°C. 3: Ni²⁺ affinity chromatography purification. (B): Gel filtration chromatography purification.

Table S1. Biological properties of MPXV-A29-specific mAbs generated in this study.

Table S2. The screening of the matched mAbs against A29 protein by the checkerboard method.

Fig. S2. Optimization of colloidal gold and antibody binding conditions. (A) The state of colloidal gold with different amount of antibody. (B) The state of colloidal gold at different pH value.




Fig. S1. Soluble expression of recombinant A29 protein. (A): M: The molecular weight standard of protein. 1: Supernatant induced at 37°C. 2: Precipitation induced at 37°C. 3: Ni²⁺ affinity chromatography purification. (B): gel filtration chromatography purification.

Table S1. Biological properties of MPXV-A29-specific mAbs generated in this study.

mAbs	Isoytppe	ELISA titer
1D6	IgG1	10 ⁻⁵
2D2	IgG1	10 ⁻⁶
3D6	IgG1	10 ⁻⁴
4C7	IgG1	10 ⁻⁴
5D8	IgG1	10 ⁻⁵
6D12	IgG2a	10 ⁻⁴
7C5	IgG1	10 ⁻⁶
8F8	IgG1	10 ⁻⁶
9D8	IgG1	10 ⁻⁵
10F8	IgG1	10 ⁻⁴
11F5	IgG1	10 ⁻⁴
12G7	IgG2b	10 ⁻⁶
13G4	IgG2a	10 ⁻⁵
14F12	IgG1	10 ⁻⁵
15F11	IgG1	10 ⁻⁴
16C4	IgG1	10 ⁻⁵
17C9	IgG1	10 ⁻⁶
18A8	IgG1	10 ⁻⁶
19D11	IgG2a	10 ⁻⁵
20B8	IgG2a	10 ⁻⁴
21B4	IgG1	10 ⁻⁵
22D3	IgG1	10 ⁻⁵
23F3	IgG2b	10 ⁻⁶
24H2	IgG1	10 ⁻⁵
25C1	IgG1	10 ⁻⁴

Table S2. The screening of the matched mAbs against A29 protein by the checkerboard method.

Gold labelled antibodies	Coated antibodies							
	3D6	5D8	7C5	8F8	14F12	15F11	20B8	25C1
3D6	-	-	++	-	+	-	-	-
5D8	-	-	+++	-	+	+	+	-
7C5	+	-	+	+	-	-	-	+
8F8	+	-	-	-	-	-	+	-
14F12	-	-	-	-	+	-	-	+
15F11	+	-	-	-	-	-	-	+
20B8	-	-	-	-	-	-	-	+
25C1	++	++	-	-	-	+++	-	+




Fig. S2. Optimization of colloidal gold and antibody binding conditions. (A) The state of colloidal gold with different amount of antibody. (B) The state of colloidal gold at different pH value.