

Supplementary Information

Batch analysis of microplastics in water using multi-angle static light scattering and chemometric methods

Mehrdad Lotfi Choobbari^a, Leonardo Ciaccheri^b, Tatevik Chalyan^c, Barbara Adinolfi^b, Hugo Thienpont^c, Wendy Meulebroeck^c, and Heidi Ottevaere^{c*}

^a Vrije Universiteit Brussel, Department of Applied Physics and Photonics, Brussels Photonics, Pleinlaan 2, 1050 Brussels, Belgium

^b CNR-Istituto di Fisica Applicata "Nello Carrara", Via Madonna del Piano 10 - 50019, Sesto Fiorentino (FI) -Italy

^c Vrije Universiteit Brussel and Flanders Make, Department of Applied Physics and Photonics, Brussels Photonics, Pleinlaan 2, 1050 Brussels, Belgium

*Corresponding author

Heidi Ottevaere

Phone/fax: +32 472 38 67 12

Email: Heidi.Ottevaere@vub.be

Table of content

Fig. S1. Scattering patterns of all monodisperse samples and their logarithmic derivative....	S2
Fig. S2. Loading plot of PCA analysis.....	S2
Fig. S3. Scree plot of PCA analysis.....	S3
Fig. S4. Scattering patterns of polydisperse samples and their logarithmic derivative.....	S3
Fig. S5. Projection of PCA maps of polydisperse samples onto the PCA map of monodisperse samples.....	S3
Table S1. Classification results.....	S4




Fig. S1. Logarithmic scattering intensity patterns of different concentrations of PS MPs (a), and their corresponding 1st derivative (b).




Fig. S2. Loading plots of 4 PCs that were used to discriminate between the size of the particles in monodisperse samples. (a) shows the loading of PC1 and PC2, and (b) shows the loading of PC3 and PC4.




Fig. S3. The contribution of each PC in explaining the total variance of data in the calibration set.




Fig. S4. Logarithmic scattering intensity patterns of polydisperse samples (a), and their corresponding 1st derivative (b).




Fig. S5. Projection of the PCA maps of the small (a) and large (b) PS polydisperse samples onto the PCA maps of their constituent particles. The pink diamonds and the pale blue stars show the PCA maps of PS polydisperse samples with equal and unequal concentrations of each size, respectively.

Table S1. The classification of PS, PE and PMMA polydisperse samples using the PCA-LDA model that was trained by different numbers of PCs. S, M, and L stand for “Small”, “Medium” and “Large”, respectively.

Type of MPs	Mixture	Assigned class		
		2PCs*-LDA	3PCs-LDA	4PCs-LDA
PS	0.50 - 0.75 - 1 μm	S	S	S
	0.50 - 0.75 μm	S	S	S
	0.50 - 1 μm	S	S	S
	0.75 - 1 μm	S	S	S
	2 - 3 - 4.5 μm	M	M	M
	2 - 3 μm	M	M	M
	2 - 4.5 μm	M	M	M
	3 - 4.5 μm	M	M	M
	3 - 10 μm	L	L	L
	3 - 6 - 10 μm	L	L	L
	3 - 6 μm	L	L	L
	6 - 10 μm	L	L	L
	0.50 - 0.75 - 1 μm (221)**	S	S	S
	0.50 - 0.75 - 1 μm (211)	S	S	S
	0.50 - 0.75 - 1 μm (122)	S	S	S
	0.50 - 0.75 - 1 μm (112)	S	S	S
	2 - 3 - 4.5 μm (221)	M	M	M
	2 - 3 - 4.5 μm (211)	M	M	M
	2 - 3 - 4.5 μm (122)	M	M	M
	2 - 3 - 4.5 μm (112)	M	M	M
	3 - 6 - 10 μm (221)	L	L	L
	3 - 6 - 10 μm (211)	L	L	L
	3 - 6 - 10 μm (122)	L	L	L
3 - 6 - 10 μm (112)	L	L	L	
PE	1 – 4 μm	M	M	M
	3 – 16 μm	L	L	L
	20 – 27 μm	L	L	L
PMMA	3 -10 μm	L	L	L

* Is the number of PCs that is used to train the LDA classifier.

** Demonstrates the concentration ratio of each corresponding size in a polydisperse sample.