## Hierarchically patterned protein scaffolds with nano-fibrillar and micro-lamellar structures modulate neural stem cell homing and promote neuronal differentiation

Jiaqi Hu<sup>a</sup>, Chenlin Li<sup>a</sup>, Zhangze Yang<sup>a</sup>, Qi Wu<sup>a</sup>, Jie Wang<sup>a, b</sup>, Zongpu Xu<sup>a, b</sup>, Yuyin

Chen<sup>a, b</sup>, Quan Wan<sup>a, b</sup>, Yajun Shuai<sup>\*, a, b</sup>, Shuxu Yang<sup>\*, c</sup>, Mingying Yang<sup>\*, a, b</sup>

a Institute of Applied Bioresource Research, College of Animal Science, Zhejiang
University, Hangzhou, Zhejiang 310058, China

key Laboratory of Silkworm and Bee Resource Utilization and Innovation of
Zhejiang Province, Hangzhou, Zhejiang, 310058, China

c Department of Neurosurgery, Sir Run Run Shaw Hospital, School of Medicine, Zhejiang University, Hangzhou, Hangzhou, 310016, China

\* Corresponding author.

## Supplementary materials

Experimental methods

Supplementary Figures 1-12

## **Experimental methods**

*Preparation of SF-GO film:* SF film was prepared by extension method with 2% solution. After ultrasonic dispersion for 30 min, 2% GO solution was well-dispersed. Then SF film was immersed into the GO solution, slowly oscillating at room temperature to form the nano structure on the film surface. After processing for two hours, the SF-GO film was washed twice and air dried.

*The viability of NSCs on the SF-GO film:* After co-cultured with the material for 24 h, the LDH cytotoxicity of NSCs was detected by LDH Cytotoxicity Assay Kit (Beyotime, China) according the manufacturer's instructions. The cell cultured in serum-free medium on the plate without any treatment was treated as the control group. To examine the viability of NSCs on the SF-GO film, the Live/Dead assay was conducted using 1  $\mu$ mol/L calcein AM and 4  $\mu$ mol/L ethidium homodimer-1 (Life Technologies kit) for 30 min at 37 °C.

*RNA extraction, Library preparation, and Sequencing:* Total RNA was extracted from the NSCs growing on the scaffolds using TRIzol® Reagent according the manufacturer's instructions. Then RNA quality was determined by 5300 Bioanalyser (Agilent) and quantified using the ND-2000 (NanoDrop Technologies). Only high-quality RNA sample (OD260/280=1.8~2.2, OD260/230≥2.0, RIN≥6.5, 28S:18S≥1.0, >1µg) was used to construct sequencing library. mRNA libraries were generated by Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus. They were then sequenced using an Illumina NovaSeq 6000.



**Supplementary Figure 1.** The XRD patterns of GO, FSF and FSF<sub>99</sub>GO<sub>1</sub> scaffolds. The peak value of GO, FSF and FSF<sub>99</sub>GO<sub>1</sub> appeared at 11°, 21.4° and 20.6°, respectively.



**Supplementary Figure 2.** The Raman spectra of GO, CSF,  $CSF_{99}GO_1$ , FSF,  $FSF_{199}GO_1$ ,  $FSF_{99}GO_1$  and  $FSF_{49}GO_1$  scaffolds. FSF and CSF scaffolds below showed random coil structure. When GO was added, the G band and D band appeared.



**Supplementary Figure 3.** Fluorescence staining images of NSCs proliferating on SF-GO film for 3 d. (Blue: nucleus; Green: cytoskeleton) Scale bars represent 20  $\mu$ m. Dotted box: cells showing spindle spread. Arrow: cell pseudopodia.



**Supplementary Figure 4.** LDH cytotoxicity of NSCs co-cultured with the material for 24 h. Sample pairs (CSF vs  $CSF_{99}GO_1$ , FSF vs  $FSF_{99}GO_1$ , CSF vs FSF,  $CSF_{99}GO_1$  vs  $FSF_{99}GO_1$ ) were analyzed with the Student's t-test. One-way ANOVA analysis was applied among FSF-GO groups. \**P*<0.05, (Mean ± SD, n=3)



**Supplementary Figure 5.** Immunofluorescence-stained neural cells cultured on SF-GO films for 14 d. (A) Neurons are stained green (beta III tubulin: green), astrocytes are stained red (GFAP: red) and nuclei are stained blue (DAPI: blue). (B) Percentage of tubulin versus DAPI-stained area. (C) Percentage of GFAP versus DAPI-stained area. \*P<0.05, \*\*P<0.01. (Mean ± SD, n=3)



**Supplementary Figure 6.** Immunofluorescence-stained neural cells cultured for 14 d. Astrocyte stained red (GFAP: red) and nuclei stained blue (DAPI: blue).