Electronic Supplementary Material (ESI) for RSC Chemical Biology. This journal is © The Royal Society of Chemistry 2023

Supporting Information Gut et. al.

## **Electronic Supporting Information (ESI)**

Designed Ankyrin Repeat Proteins for detecting prostate specific antigen expression in vivo

Melanie Gut, <sup>1</sup> Birgit Dreier, <sup>2</sup> Sven Furler, <sup>2</sup> Jens Sobek, <sup>3</sup> Andreas Plückthun, <sup>2</sup> and Jason P. Holland <sup>1\*</sup>

<sup>1</sup> University of Zurich, Department of Chemistry, Winterthurerstrasse 190, 8057 Zurich, Switzerland

<sup>2</sup> University of Zurich, Department of Biochemistry, Winterthurerstrasse 190, 8057 Zurich, Switzerland

<sup>3</sup> Functional Genomics Center Zurich, Eidgenössische Technische Hochschule (ETH) Zurich and University of Zurich, Zurich, Switzerland

## \* Corresponding Author:

Prof. Dr Jason P. Holland

Tel: +41.44.63.53.990

E-mail: jason.holland@chem.uzh.ch

Website: www.hollandlab.org

Twitter: @HollandLab\_

## First author:

Melanie Gut

E-mail: melanie.gut@uzh.ch

## **Table of Contents**

General information	3
Chemicals	
Mass spectrometry	3
Enzyme-linked immunosorbent assay (ELISA)	3
Size-exclusion chromatography	4
Spectroscopy	
Size-exclusion chromatography coupled to multi-angle light scattering (SEC-MALS) detector	4
Sodium dodecyl sulfate - polyacrylamide gel electrophoresis (SDS-PAGE)	4
Native-PAGE analysis of proteins	5
Radiochemical Methods	5
Gallium-68	5
Radio-TLC	5
Radioactive size-exclusion chromatography (radio-SEC)	<i>t</i>
Quantification of radioactivity	<i>t</i>
Cell culture	6
General methods in cell culture	<i>t</i>
LNCaP cells	<i>t</i>
Animals and xenograft models	
Ex vivo biodistribution studies	
Effective half-life measurement	
Statistical analysis	
DARPin expression and selection	7
Biotinylation of PSA	
Ribosome display and FRET analysis	8
ELISA of 27 DARPin constructs	9
Other exclusion criteria including SEC-MALS	
Cloning of PSA-DARPins with a Cys C-terminus	11
DARPin expression	11
DARPin purification	
Size-exclusion chromatography coupled to multi-angle light scattering (SEC-MALS) detector	13
Bioconjugation of DARPin-Cys proteins with a NODAGA-maleimide	15
NODAGA-maleimide bioconjugation	
SDS-PAGE gel with DARPin-Cys constructs	16
Synthesis of natGa-NODAGA-DARPin complexes	16
Qualitative assessment of PSA binding of NODAGA-DARPin constructs	22
Evaluation of PSA binding of DARPin constructs after NODAGA modification by automated	SEC chromatography
Native-PAGE gel analysis of PSA-DARPin-Cys complexes	22
Surface plasmon resonance measurements for quantitative assessment of PSA binding of DARP	in-Cys and NODAGA
DARPins	23
Radiolabelling of NODAGA-DARPins	27
<sup>68</sup> Ga-radiolabelling of NODAGA-DARPins	27
Molar activities of <sup>68</sup> Ga-radiollabelling of [ <sup>68</sup> Ga]GaNODAGA-DARPins	27
Stability of [68Ga]GaNODAGA-DARPins	29
Incubation in PBS	29
Incubation with excess EDTA	29
Incubation in human serum	29
Radioactive PSA-binding assay	
Tumor uptake studies and ex vivo biodistribution studies of [68Ga]GaNODAGA-DARPin radiotraces	
Effective half-life measurement	32
Ex vivo biodistribution studies	32
D. C.	2.0

## **General information**

#### Chemicals

Unless otherwise stated, all other chemicals were of reagent grade and purchased from SigmaAldrich (St. Louis, MO), Merck (Darmstadt, Germany), Tokyo Chemical Industry (Eschborn, Germany), abcr (Karlsruhe, Germany) or CheMatech (Dijon, France). Water (>18.2 M $\Omega$ ·cm at 25 °C, Merck Millipore Milli-Q® Reference A+, with an additional 0.22  $\mu$ m Millipack Express40 filter) was used without further purification for chemical and biological experiments, and additionally purified by aid of a Chelex® resin (100 sodium form, 100 – 200 mesh particle size, Sigma Aldrich) for radiochemical experiments.

## Mass spectrometry

Liquid chromatography high resolution electrospray ionisation mass spectra (LC-HR-ESI-MS) were measured by the mass spectrometry service at the Department of Chemistry, University of Zurich.

Acquity UPLC (Waters, Milford, USA) connected to an Acquity e $\lambda$  diode array detector and a Synapt G2 HR-ESI-QTOF-MS (Waters, Milford, USA); injection of 1 – 5 μL sample (concentration ca. 10-100 μg mL<sup>-1</sup> in the indicated solvent); UV spectra were recorded from 200 – 600 nm at 1.2 nm resolution and at a sampling rate of 20 points s<sup>-1</sup>; ESI: positive ionisation mode, capillary voltage 3.0 kV, sampling cone 40 V, extraction cone 4 V, N<sub>2</sub> cone gas 4 L h<sup>-1</sup>, N<sub>2</sub> desolvation gas 800 L min<sup>-1</sup>, source temperature 120 °C; mass analyser in resolution mode: mass range 100 – 2,000 m/z with a scan rate of 1 Hz; mass calibration to <2 ppm within 50–2500 m/z with a 5 mM aqueous solution of NaHCO<sub>2</sub>, lock-masses: m/z 195.0882 (caffein, 0.7 ng mL<sup>-1</sup>) and 556.2771 (leucine-enkephalin, 2 ng mL<sup>-1</sup>).

Acquity BEH C18 HPLC column (1.7  $\mu$ m particle size, 2.1  $\times$  50 mm, Waters) kept at 30 °C; elution at a flow rate of 400  $\mu$ L min<sup>-1</sup> with solvent A: H<sub>2</sub>O + 0.02% TFA and solvent B: CH<sub>3</sub>CN + 0.02% TFA, linear gradient from 10 – 95% B within 3 min., then isocratic for 2 min. at 95% B.

Acquity BEH C8 HPLC column (1.7  $\mu$ m particle size, 2.1  $\times$  100 mm, Waters) kept at 30 °C; elution at a flow rate of 400  $\mu$ L min<sup>-1</sup> with solvent A: H<sub>2</sub>O + 0.02% TFA and solvent B: CH<sub>3</sub>CN + 0.02% TFA, linear gradient from 10 – 70% B within 10 min.

#### Enzyme-linked immunosorbent assay (ELISA)

Biotinylated PSA was immobilised as follows: 384-well high binding microplates (MTP384 Greiner Bio-One GmbH, Frickenhausen, Germany) were coated with 20 μL of streptavidin (66 nM in PBS) overnight at 4 °C and blocked with 0.2% BSA (*w/v*) (45 μL PBS, supplemented with 0.1 [v/v] Tween20 and 0.2% [w/v] BSA [PBS-TB]) for 1 h at 4 °C with orbital shaking (Titramax 1000, at 900 rpm). After washing, biotinylated PSA (20 μL of 50 nM in PBS-TB supplemented with 1 mM DTT) was immobilised for 2 h, 4 °C, 900 rpm. IMAC-purified DARPin constructs (expressed and selected as described below; 20 μL, 100 nM in PBS supplemented with 0.1 [v/v] Tween20 [PBS-T] and 1 mM DTT) were applied to wells with or without immobilised PSA for 1.5 h, 4 °C, 900 rpm. After extensive washing, binding was detected by incubation with anti-FLAG antibody (FLAG-tagged DARPins [F-DARPins]; 20 μL, 1:5'000 diluted in PBS-TB; mouse anti-FLAG IgG, F3165, Sigma-Aldrich, Buchs, Switzerland) or anti-RGS-His antibody (DARPin-Cys; 20 μL, 1:5'000 diluted in PBS-TB; mouse anti

RGA (His)<sub>4</sub> IgG1, No. 34650, Qiagen, GmbH, Hilden, Germany, which specifically discriminates between the *N*-terminal RGS-His6 tags of the DARPins and the *N*-terminal His<sub>6</sub> tag of the antigen) for 1 h, 4 °C, 900 rpm. After washing, anti-mouse IgG-alkaline phosphatase conjugate (20 μL, 1:10,000 diluted in PBS-TB; Pierce, Thermo Fisher Scientific) was applied to wells and the plate was incubated for 1 h, 4 °C, 900 rpm. After extensive washing *p*-nitrophenylphosphate (20 μL of 50 μM solution containing 1 M MgCl<sub>2</sub> and 1 M NaHCO<sub>3</sub>; Sigma-Aldrich, Buchs, Switzerland) was added and the plate was incubated for 30 min., 37 °C, 900 rpm, and the OD<sub>405-540nm</sub> was measured (Micro Plate Reader, BIO-TEK Synergy). All washing steps were performed with 110 μL PBS-T per well on an ELx405 microplate washer (BioTek, Winooski, USA). Dispensing steps were performed on a MicroFlo Select liquid dispenser (BioTek). DARPin constructs were applied with a Liquidator96 system (Mettler-Toledo GmbH, Greifensee, Switzerland).

### Size-exclusion chromatography

Automated size-exclusion chromatography was performed using an Agilent AdvanceBio SEC 130 Å column (Agilent Technologies, AdvanceBio SEC 130 Å, 2.7 μm 4.6 × 300 mm) connected to a Hitachi Chromaster Ultra Rs system equipped with a UV/visible detector (electronic absorption measured at 220 and 280 nm) as well as a radioactivity detector (FlowStar² LB 514, Berthold Technologies, Zug, Switzerland). Note that the UV/Vis detector and radioactivity detector were arranged serially with an offset time of approximately 0.10 min. The identity of the radiolabelled <sup>68</sup>Ga-compound was confirmed by co-injection with an authenticated sample of <sup>nat</sup>Ga-compound. Dependent on the system, different elution methods were used:

*Method A*: Isocratic elution with phosphate buffered saline (PBS, pH 7.4) and 20 mM arginine, flow rate: 0.35 mL min<sup>-1</sup>. Electronic absorption was measured at 220 nm.

*Method B*: Isocratic elution with phosphate buffered saline (PBS, pH 7.4), 20 mM arginine and 10 mM EDTA, flow rate: 0.35 mL min<sup>-1</sup>. Electronic absorption was measured at 280 nm.

## Spectroscopy

Electronic absorption spectra were recorded using a Nanodrop<sup>TM</sup> One<sup>C</sup> Microvolume UV-Vis Spectrophotometer (ThermoFisher Scientific, supplied by Witec AG, Sursee, Switzerland). The same system was used to evaluate protein concentrations in accordance with the manufacturer's protocol.

#### Size-exclusion chromatography coupled to multi-angle light scattering (SEC-MALS) detector

The absolute mass of protein samples was determined by using a liquid chromatography system (Agilent LC1100 Agilent Technologies, Santa Clara, USA) coupled to an Optilab rEX refractometer (Wyatt Technology, Santa Barbara, USA) and a miniDAWN three-angle light-scattering detector (Wyatt Technology). For protein separation, a 24 mL Superdex 200 10/30 column (GE Healthcare Biosciences) was run in Tris-buffered saline (TBS150, pH 7.5) at 0.5 mL min<sup>-1</sup>. Analysis was performed by using the ASTRA software (version 5.2.3.15; Wyatt Technology).

## Sodium dodecyl sulfate - polyacrylamide gel electrophoresis (SDS-PAGE)

Samples (15 – 20  $\mu$ L) were prepared for SDS-PAGE by mixing single proteins or reaction mixtures (1 – 3  $\mu$ g of protein, or 20 – 30 mg of total protein for cell lysates), sample buffer (4 × lithium dodecyl sulfate (LDS)-Buffer

for NuPage, supplemented with 0.05 M DTT unless otherwise stated) and phosphate buffered saline (PBS, pH 7.4). The resulting solution was heated at 80 °C for 10 min. (Labnet Accu Therm<sup>TM</sup>). Samples  $(10-15~\mu\text{L})$  including a protein ladder (Spectra<sup>TM</sup> Multicolor Broad Range Protein Ladder, ThermoFisher Scientific) were applied to a SDS-polyacrylamide gel (4-12% Bis-Tris, unless otherwise stated). Unless otherwise stated, gels were run in NuPage MES-SDS buffer (ThermoFisher Scientific), at 150 V for 60 min. Gels were washed in deionised water before staining or transfer procedures.

### Native-PAGE analysis of proteins

Samples  $(15-20 \,\mu\text{L})$  were prepared for native-PAGE by mixing single proteins or reaction mixtures with sample buffer (Native-TrisGly Sample Buffer 2×, NOVEX Life Technologies) and phosphate buffered saline (PBS, pH 7.4). The resulting solutions were mixed and samples  $(10-15 \,\mu\text{L})$  including a protein ladder (Spectra<sup>TM</sup> Multicolor Broad Range Protein Ladder, ThermoFisher Scientific) were applied to a Bis-Tris Gel (Native-PAGE 4–16%, ThermoFisher Scientific). The gel was run in MES running buffer (0.05 M MES, pH 6) at 150 V for 60 min. The gel was washed with deionised water before the staining.

## **Radiochemical Methods**

All instruments for measuring radioactivity were calibrated and maintained in accordance with previously reported routine quality control procedures.<sup>1</sup> For analysis appropriate background and decay corrections were applied as necessary.

## Gallium-68

[68Ga][Ga(H<sub>2</sub>O)<sub>6</sub>]Cl<sub>3</sub> (aq.) was obtained from <sup>68</sup>Ge/<sup>68</sup>Ga-generators (Eckert&Ziegler, Model IGG100 Gallium-68 Generator), and eluted with 0.1 M HCl (aq.). The eluted <sup>68</sup>Ga activity was trapped and purified by using a strong cation exchange column (Strata-XC, [SCX], Eckert&Ziegler). [<sup>68</sup>Ga][Ga(H<sub>2</sub>O)<sub>6</sub>]Cl<sub>3</sub> (aq.) was eluted from the SCX cartridge by using a solution containing 0.13 M HCl (aq.) and approx. 5 M NaCl (aq.) (SCX eluent). For radiolabelling experiments, the <sup>68</sup>Ga stock solution (in ~0.1 M HCl, ca. 5 – 20 MBq) was typically added as the limiting reagent to an aqueous solution of unconjugated compounds (ca. 4 nmol for DARPins) buffered with NaOAc (0.1 M, pH 4.4). The reaction was monitored by using radio-iTLC and characterised by analytical size-exclusion chromatography.

## Radio-TLC

Radioactive reactions were monitored by using instant thin-layer chromatography (iTLC). Glass-fibre iTLC plates impregnated with silica-gel (iTLC-SG, Agilent Technologies) were developed in citrate buffer (0.2 M, pH 4.5, >18.2 M $\Omega$ ·cm H<sub>2</sub>O) and analysed on a radio-TLC detector (SCAN-RAM, LabLogic Systems Ltd, Sheffield, United Kingdom). Radiochemical conversion (RCC) was determined by integrating the data obtained by the radio-TLC plate reader and determining both the percentage of radiolabelled product (retention factor  $R_{\rm f} = 0.0 - 0.1$ ) and 'free' <sup>68</sup>Ga ( $R_{\rm f} = 0.8 - 1.0$ ). Integration and data analysis were performed by using the software Laura version 5.0.4.29 (LabLogic).

## Automated radioactive size-exclusion chromatography (radio-SEC)

Samples were analysed by using two methods of SEC. First, radiochemical purities (RCPs) of labelled protein samples were determined by using automated size-exclusion chromatography (Agilent AdvanceBio SEC 130 Å column, method B) as outlined above. Radiochemical purity (RCP) was determined by integrating the data obtained by the automated radio-SEC and determining both the percentage of radiolabelled product and 'free'  $^{68}$ Ga after appropriate background and decay corrections. Integration and data analysis were performed by using the software RadioStar version 5.0. As a second method, manual SEC procedures were used involving PD-10 (Sephadex G-25 resin, 85-260 µm, 14.5 mm ID × 50 mm, >30 kDa, GE Healthcare) or NAP<sup>TM</sup>-25 columns (Cytiva) desalting columns. For analytical procedures, PD-10 and NAP-25 columns were eluted with sterile PBS. A total of 20 × 200 µL (PD-10), or 30 × 300 µL (NAP-25) fractions were collected up to a final elution volume of 6.5 mL (PD-10) or 11.5 mL NAP-25). The loading/dead-volume of the PD-10 and NAP-25 columns is precisely 2.50 mL which was discarded prior to aliquot collection. PD-10 and NAP-25 SEC columns were also used for preparative purification and reformulation of radiolabelled products by collecting a fraction of the eluate corresponding to the high molecular weight protein (>30 kDa fraction) eluted in ranges indicated for each experiment.

## Quantification of radioactivity

For quantification of radioactivity, each fraction was measured on a gamma counter (HIDEX Automatic Gamma Counter, Hidex AMG, Turku, Finland) by using a counting time of 30 s and an energy window between 480 – 558 keV for <sup>68</sup>Ga (511 keV emission), and 15 – 2047 keV for <sup>111</sup>In. Appropriate background and decay corrections were applied throughout.

## Cell culture

#### General methods in cell culture

All cells were cultured at 37 °C in a humidified 5%  $CO_2$  atmosphere. All cell media was supplemented with fetal bovine serum (FBS, 10% [v/v], ThermoFisher Scientific) and penicillin/streptomycin (P/S, 1% [v/v] of penicillin 10,000 U mL<sup>-1</sup> and streptomycin 10 mg mL<sup>-1</sup>). Cells were pelleted (100 g, 5 min.) and resuspended in media after trypsinisation.

## LNCaP cells

The human prostate cancer cell line LNCaP was obtained from the American Type Culture Collection (ATCC-CRL-1740<sup>TM</sup>, Manassas, VA). Cells were cultured at 37 °C in a humidified 5% CO<sub>2</sub> atmosphere in RPMI-1640 (without phenol-red) containing [+]-L-glutamine (2.5 mM). Cells were grown by serial passage and were harvested by using trypsin (0.1%).

## Animals and xenograft models

All experiments involving mice were conducted in accordance with an animal experimentation license approved by the Zurich Canton Veterinary Office, Switzerland (Jason P. Holland). Experimental procedures also complied with guidelines issued in the *Guide for the Care and Use of Laboratory Animals*. Female athymic nude mice  $(Crl:NU(NCr)-FoxnI^{nu}, 17-30 \text{ g}, 6-8 \text{ weeks old})$  were obtained from Charles River Laboratories Inc. (Freiburg

International). Tumours developed after a period of between 3-8 weeks. Tumour volume ( $V/mm^3$ ) was estimated by external Vernier caliper measurements of the longest axis, b/mm. The tumours were assumed to be spheroidal and the volume was calculated in accordance with Equation S1. All mice injected with cancer cells developed tumours at the University of Zurich Laboratory Animal Services Center vivarium for at least 1 week prior to implanting tumour cells. Mice were provided with food and water ad libitum. Tumours were induced on the right shoulder or flank by subcutaneous (s.c.) injection of LNCaP ( $5.5 \times 10^6$ ) cells. The cells were injected in a 150-200 μL suspension of a 1:1 v/v mixture of cell medium (RPMI-1640) and reconstituted basement membrane (Corning® Matrigel® Basement Membrane Matrix, obtained from VWR International). Tumours developed after a period of between 3-8 weeks. Tumour volume ( $V/mm^3$ ) was estimated by external Vernier caliper measurements of the longest axis, a/mm, and the axis perpendicular to the longest axis, b/mm. The tumours were assumed to be spheroidal and the volume was calculated in accordance with Equation S1. All mice injected with cancer cells developed tumours and the average volume of the LNCaP tumours was  $624.56 \pm 714.01 \text{ mm}^3$  ( $151.25 \pm 78.46 \text{ mm}^3$ , n = 35 mice).

Equation S1. Approximation of tumour volume.

$$V = \frac{4\pi}{3} \cdot \left(\frac{a}{2}\right)^2 \cdot \left(\frac{b}{2}\right)$$

#### Ex vivo biodistribution studies

Biodistribution studies were conducted at the specified time points after completing the final measurements for estimating the effective half-live  $t_{1/2}$ (eff). Animals (n = 3 - 4 mice per group) were anaesthetised individually by isoflurane and euthanised by isoflurane asphyxiation followed by terminal exsanguination. A total of 15 tissues (including the tumour) were removed, rinsed in water (except the fat, skin, muscle, and bone), dried in air for approx. 1 min., weighed and counted on a calibrated gamma counter (HIDEX Automated AMG Gamma Counter) for accumulation of activity. The mass of radiotracer formulation injected into each animal was measured and used to determine the total number of counts per minute (cpm) injected into each mouse by comparison to a standard syringe of known activity and mass. Count data were background- and decay-corrected, and the tissue uptake for each sample (determined in units of percentage injected dose per gram [%ID g<sup>-1</sup>]) was calculated by normalisation to the total amount of activity injected for each individual animal.

## Effective half-life measurement

The effective half-live  $t_{1/2}$ (eff) was calculated from the measurement of total internal radioactivity in the mouse models over time by using a dose calibrator (ISOMED 2010, NuviaTech Healthcare).

#### Statistical analysis

Data were analysed using GraphPad Prism and Excel. Unpaired, two-tailed Student's t-test were used with confidence levels of p <0.05, p <0.01, and p <0.001 suggesting low, medium, and high statistical significance differences.

## **DARPin** expression and selection

## **Biotinylation of PSA**

PSA (250  $\mu$ L, 1.5 mg mL<sup>-1</sup>, obtained from BBI<sup>TM</sup> Solutions) was re-buffered (Amicon 0.5, 10 kDa MWCO, 3,000 g, 15 min.) in PBS buffer (pH 7.4), and the concentration was measured to yield purified PSA (156.5  $\mu$ L, 2.33 mg mL<sup>-1</sup>). Biotin-NHS (6.2  $\mu$ L, 62 nM 5 eq., 10 mM in DMSO) was added to purified PSA (156.5  $\mu$ L, 368  $\mu$ g, 12.58 nmol). The reaction was mixed by gentle pipetting and incubated at room temperature (30 min., 300 rpm). The

reaction mixture was purified, to remove excess biotin, by spin filtration (Amicon 0.5, 10 kDa MWCO, 3,000 g, 15 min.,  $3 \times 0.5$  mL) into PBS + 0.1% w/v NaN<sub>3</sub>. The purified PSA-biotin conjugate was then stored at 4 °C until further use.

#### Ribosome display and FRET analysis

To generate PSA-specific DARPin binders, biotinylated PSA was immobilized on either MyOne T1 streptavidin-coated beads (Pierce) or Sera-Mag neutravidin-coated beads (GE), depending on the selection round, and these beads were alternated. Ribosome display selections were performed similar to those described earlier,<sup>4</sup> using a semi-automatic KingFisher Flex MTP96 well platform.


The library includes N3C-DARPins with the original randomization strategy as reported<sup>5</sup> but includes a stabilized C-cap.<sup>6–8</sup> Additionally, the library is a mixture of DARPins with randomized and non-randomized N- and C-terminal caps, respectively<sup>7,9</sup> and successively enriched pools were cloned as intermediates in a ribosome display-specific vector.<sup>9</sup> Selections were performed over four rounds with decreasing target concentration and increasing washing steps to enrich for binders with high affinities.

The final enriched pool was cloned as fusion construct into a bacterial pQE30 derivative vector with an *N*-terminal MRGS(H)<sub>8</sub> tag and *C*-terminal FLAG tag *via* unique BamHI and HindIII sites containing *laclq* for expression control. After transformation of *E. coli* XL1-blue, 380 single DARPin clones selected to bind to PSA were expressed in 96-well format and lysed by addition of B-Per Direct detergent plus Lysozyme and Nuclease (Pierce). These bacterial crude extracts of single DARPin clones were subsequently used in a Homogeneous Time Resolved Fluorescence (HTRF)-based screen to identify potential binders. Binding of the FLAG-tagged DARPins (abbreviated as F-DARPin) to streptavidin-immobilized biotinylated PSA was measured using FRET (donor: Streptavidin-Tb cryptate (610SATLB, Cisbio), acceptor: mAb anti FLAG M2-d2 (61FG2DLB, Cisbio). Further HTRF measurement against 'No Target' allowed for discrimination of PSA-specific hits. Experiments were performed at room temperature in white 384-well Optiplate plates (PerkinElmer) using the Taglite assay buffer (Cisbio) at a final volume of 20 µl per well. FRET signals were recorded after an incubation time of 30 minutes using a Varioskan LUX Multimode Microplate (Thermo Scientific). HTRF ratios were obtained by dividing the acceptor signal (665 nm) by the donor signal (620 nm) and multiplying this value by 10,000 to derive the 665/620 ratio. The background signal was determined by using reagents in the absence of DARPins.

From the initial hits of the 380 analysed F-DARPin clones, 32 were chosen for further analysis and their sequence determined. 27 were identified as true single clones. They were expressed in a 96-well format and IMAC purified. The purified F-DARPins were used in a hit validation to bind immobilized PSA by ELISA using an anti-FLAG antibody and SEC (see sections below).

## **ELISA of 27 DARPin constructs**

ELISA binding of the 27 F-DARPins was evaluated by using the standard protocol (using anti-FLAG antibody; see above).



**Figure S1**. Bar chart showing the results from ELISA measurements. The specific binding of the 27 constructs towards PSA was tested. The 12 selected F-DARPins are highlighted in dark cyan. Data presented show the relative absorbance for PSA-binding F-DARPins ( $A - \S$ ) and the non-binding control (F-DARPin 13, E3\_5). Red bars show the non-specific binding control where no PSA was coated on the plate. As a positive control DARPin Off7 was added on a Maltose-Binding Protein (MBP) coated well. Error bars represent one standard deviation around the mean (n = 2 measurements per sample).

## Other exclusion criteria including SEC-MALS

SEC-MALS measurements of the 27 single clones were conducted in accordance with the standard procedure described above.

**Table S1.** Table of the 27 constructs showing a summary of the selection criteria: ELISA binding, SEC profile and sequence alignment. The 12 selected F-DARPins are highlighted in green, the exclusion reasons are highlighted in red.

Initial DARPins	ELISA binding	SEC	Selected F- DARPins
A	high	monomeric	F-DARPin 1
В	high	not monomeric	
С	low	broad	
D	high	monomeric	F-DAPRin 2
Е	low	monomeric	
F	low	monomeric	
G	high	monomeric	F-DARPin 3
Н	low	monomeric	F-DARPin 4
I	high	not monomeric	
J	low	monomeric	F-DARPin 5
K	medium	not monomeric	
L	medium	monomeric	F-DARPin 6
M	low	monomeric	
N	high	not monomeric	
О	low	monomeric	F-DARPin 7
P	low	monomeric	
Q	low	monomeric	
R	medium	monomeric	F-DARPin 8
S	medium	not monomeric	
T	low	monomeric	
U	high	monomeric	F-DARPin 9
V	medium	not monomeric	
W	high	monomeric	F-DARPin 10
X	medium	monomeric	F-DARPin 11
Y	high	monomeric	F-DARPin 12
Z	low	monomeric	
§	high	not monomeric	
E3_5	negative control		F-DARPin 13

Based on ELISA and SEC-MALS results 12 PSA-binding F-DARPins and one negative control (E\_5, F-DARPin 13) were selected. The corresponding amino acid sequences are shown in **Figure S2**.

	N-cap		1 <sup>st</sup> repeat				
F-DARPin	M R G S H H H H H H H H G S D L G K K L L E A A X X G Q D D E V R	ILMANGADVNA	XDXXGYTPLHLAAXXGHLEI	VEVLLKXGADVNA			
1	M R G S H H H H H H H H G S D L G K K L L E A A D N G Q D D E V R	ILMANGADVNA					
2	M R G S H H H H H H H H G S D L G K K L L E A A E T G Q D D E V R	ILMANGADVNA					
3	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA					
4	M R G S H H H H H H H H G S D L G K K L L E A A W F G H L D E V R	ILMANGADVNA					
5	M R G S H H H H H H H H G S D L G K K L L E A A W A G Q L D E V R	ILMANGADVNA					
6	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA	SDYVGYTPLHLAADYGHLEI	VEVLLKTGADVNA			
7	M R G S H H H H H H H H G S D L G K K L L E A A T M G H L D E V R	ILMANGADVNA					
8	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA					
9	M R G S H H H H H H H H G S D L G K K L L E A A M R G H L D E V R	ILMANGADVNA					
10	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA					
11	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA					
12	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA					
13	M R G S H H H H H H H H G S D L G K K L L E A A R A G Q D D E V R	ILMANGADVNA	TDNDGYTPLHLAASNGHLEI	VEVLLKNGADVNA			
	2 <sup>nd</sup> repeat		3 <sup>rd</sup> repeat				
F-DARPin	X D X X G X T P L H L A A X X G H L E I V E V L L K X G A D V N A	XDXXGXTPLHL	AAXXGHLEIVLLKHGADVNA				
1	KDHIGSTPLHLAAIEGHLEIVEVLLKAGADVNA	TDWFGWTPLHL	AADWGHLEIVLLKHGADVNA	<b>、</b>			
2	WDVMGQTPLHLAALEGHLEIVEVLLKTGADVNA	V D A F G W T P L H L	AAVDGHLEIVLLKHGADVNA	<b>、</b>			
3	RDRLGQTPLHLAAIDGHLEIVEVLLKAGADVNA	TDQFGFTPLHL	AAQEGHLEIVLLKHGADVNA	•			
4	D D L L G Q T P L H L A A M H G H L E I V E V L L K A G A D V N A	RDYFGYTPLHL	AAIWGHLEIVLLKHGADVNA	•			
5	WDFIGQTPLHLAAYRGHLEIVEVLLKTGADVNA	QDYIGWTPLHL	AAMQGHLEIVLLKHGADVNA	•			
6	LDRIGQTPLHLAAWQGHLEIVEVLLKTGADVNA	IDYFGWTPLHL	AASEGHLEIVLLKHGADVNA	<b>、</b>			
7	Q D M L G R T P L H L A A H N G H L E I V E V L L K T G A D V N A	SDWFGWTPLHL	AAFDGHLEIVLLKHGADVNA	•			
8	H D Y I G Q T P L H L A A S Y G H L E I V E V L L K A G A D V N A	QDFMGLTPLHL	A A EWGHLEIVLLKHGADVNA	4			
9			AAVYGHLEIVLLKAGADVNA				
10	H D R L G Q T P L H L A A M L G H L E I V E V L L K T G A D V N A						
11	WDIIGNTPLHLAAMQGHLEIVEVLLKTGADVNA						
12	WDWLGQTPLHLAALEGHLEIVEVLLKTGADVNA						
13	SDLTGITPLHLAAATGHLEIVEVLLKHGADVNA		AAKYGHLEIVLLKHGADVNA	4			
	<del> </del>	FLAG-tag					
F-DARPin	Q D K F G K T X F D L A I D N G N E D I A E V L Q K A A K L N D Y						
1	Q D K F G K T A F D L A I D N G N E D I A E V L Q K A A K L N D Y						
2	Q D V W G E T A F D L A A W W G N E D I A E V L Q K A A K L N D Y						
3	Q D K F G K T A F D L A I D N G N E D I A E V L Q K A A K L N D Y						
4	Q D K F G K T A F D L A I D N G N E D I A E V L Q K A A K L N D Y						
5	Q D K F G K T A F D L A I D N G N E D I A E V L Q K A A K L N D Y						
6	Q D K F G K T A F D L A I D N G N E D I A E V L Q K A A K L N D Y						
7	Q D K F G K T A F D I S I D N G N E D I A E V L Q K A A K L N D Y						
8	Q D K F G K T P F D L A I D N G N E D I A E V L Q K A A K L N D Y						
9	Q DWY G E T P F D L A AWE G N E D I A E V L Q K A A K L N D Y						
10	Q D K F G K T A F D I S I D N G N E D I A E V L Q K A A K L N D Y						
11	Q D K F G K T P F D L A I D N G N E D I A E V L Q K A A K L N D Y						
12	Q D K F G K T P F D L A I D N G N Q D I A E V L Q K A A K L N D Y						
13	Q D K F G K T A F D I S I D N G N E D L A E I L Q K L N D Y	KUUDDK					

Figure S2. Amino acid sequence (single letter representation) of FLAG-tagged DARPins (F-DARPins) 1-13 showing the different repeats. The randomized regions are highlighted in blue. Note: The N-cap sequence of the F-DARPins features a His8-tag.

#### Cloning of PSA-DARPins with a Cys C-terminus

For the expression of Cys-*C*-terminally modified PSA-DARPins (DARPin-Cys) the cDNA encoding for PSA-DARPin was cloned into the *E. coli* expression plasmid pQI\_E3\_5\_GGC (Containing DARPin E3\_5 from the unselected library previously described<sup>5</sup>) *via* a *BamH*I/*Hind*III under the control of a T5 lac promoter to yield the expression vectors MRGS(H)<sub>6</sub> tag\_PSA\_DARPin\_GGC (general amino acid sequence shown in **Figure S3**, with the modified regions highlighted in green).<sup>10–12</sup>

	N-cap	1 <sup>st</sup> repeat			
DARPin-Cys	MRGSHHHHHHGSDLGKKLLEAAXXGQDDE	V R I LMANGADVNA	XDXXGYTPLHLAAXXGHLEI	VEVLLKXGADVNA	
	2 <sup>nd</sup> repeat		3 <sup>rd</sup> repeat		
DARPin-Cys	X D X X G X T P L H L A A X X G H L E I V E V L L K X G A D V	NAXDXXGXTPLHL	AAXXGHLEIVLLKHGADVNA		
	C-cap	Cys-tag			
DARPin-Cys	QDKFGKTXFDXXIDNGNEDIAEVLQKAAKLN	G G C			

Figure S3. Amino acid sequence (single letter representation) of DARPin-Cys 1 - 13 showing the different repeats. The randomized regions are highlighted in blue, with the sequence being equivalent to FLAG-tagged DARPins 1 - 13 (Figure S2); the modified region is highlighted in green. Note: The *N*-cap sequence of the DARPin-Cys clones features a His6-tag.

#### **DARPin** expression

For the expression of the DARPin-Cys constructs, the *E. coli* strain XL1-Blue (Stratagene, USA) was transformed with the plasmids described above. After overnight expression single clones were picked and used for the inoculation of 2YT medium (15 or 30 mL) supplemented with ampicillin (100 µg mL<sup>-1</sup>) and glucose (1%). The

inoculated medium was incubated overnight at 200 rpm and 37 °C. These cultures were used to inoculate 200 mL or 500 mL cultures (2YT media supplemented with 100  $\mu$ g mL<sup>-1</sup> ampicillin and 1% glucose) at an OD<sub>600</sub> of 0.1. Expression was induced with isopropyl  $\beta$ -D-thiogalactopyranoside (1 mM) at OD<sub>600</sub> of 0.8 – 1. After 4 h incubation at 200 rpm and 37 °C, the cultures were harvested by centrifugation for 15 min. (5,000 g, 4 °C). The pellet was washed with ice-cold PBS (30 mL, pH 7.4), centrifuged for 15 min (4,000 g, 4 °C) and stored at -20 °C until further purification.

## **DARPin** purification

All DARPin-Cys proteins were purified by their His6-tag using Ni-NTA immobilised metal-ion affinity chromatographic (IMAC) methods. All steps were performed at 4 °C. DARPins expressed in 200 mL cultures were thawed and resuspended in ice cold Ni-Lysis buffer (3.6 mL, 20 mM Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub>, 500 mM NaCl, 1 mM DTT, pH 7.4) supplemented with 25 U mL<sup>-1</sup> Pierce<sup>TM</sup> Universal Nuclease for Cell Lysis (ThermoFischer Scientific). The cells were lysed by sonication and centrifuged for 30 min. (20,000 g, 4 °C). The lysates were supplemented with imidazole (20 mM) and glycerol (10% v/v). The proteins were purified via automated IMAC using a HisTrap FF crude (GE Healthcare) column on an ÄKTA Pure system (GE Healthcare). Columns were equilibrated with washing buffer (20 mM Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub>, 500 mM NaCl, 20 mM imidazole, 10% v/v glycerol, 1 mM DTT, pH 7.4). Lysates were applied to the columns, washed with 5 column volumes (CV) of lowsalt buffer (20 mM Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub>, 137 mM NaCl, 20 mM imidazole, 10% v/v glycerol, 1 mM DTT, pH 7.4), 5 CV of high-salt buffer (20 mM Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub>, 1 M NaCl, 20 mM imidazole, 10 % glycerol, 1 mM DTT, pH 7.4) and 5 CV of washing buffer. The proteins were eluted with 5 CV elution buffer (20 mM Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub>, 500 mM NaCl, 250 mM imidazole, 10 % v/v glycerol, 1 mM DTT, pH 7.4). The proteins were buffer-exchanged using a HiTrap 26/10 Desalting (GE Healthcare) column. Protein-containing fractions were evaluated with a NanoDrop spectrometer, pooled, aliquoted, flash-frozen in liquid nitrogen and stored at -80 °C until further use. Protein purity and monomeric behaviour were determined to be >90% by Coomassiestained SDS-PAGE and by MALS-SEC.

Bacterial pellets from 500 mL cultures were thawed and resuspended in ice-cold TBS<sub>400</sub> (15 – 20 mL, 50 mM Tris-HCl, 400 mM NaCl, pH 7.4) supplemented with 20 mM imidazole, glycerol (10% *v/v*), and 5 mM DTT. Lysozyme (1 g L<sup>-1</sup>) was added before the cells were lysed by sonication and centrifuged for 20 min. (15,000 g, 4 °C). DARPin-Cys were purified *via* manual IMAC. Ni-NTA Superflow Resin (Quiagen) was packed in 1 – 2 mL benchtop columns (PD-10). Columns were equilibrated with lysis buffer (50 mM Tris-HCl, 400 mM NaCl, 20 mM imidazole, 1 mM DTT, 10% *v/v* glycerol, pH 7.4). Lysates were applied to the columns, washed with 2 × 10 CV washing buffer (50 mM Tris-HCl, 400 mM NaCl, 20 mM imidazole, 1 mM DTT, 10 % glycerol, pH 7.4), 10 CV of low-salt buffer (50 mM Tris-HCl, 20 mM NaCl, 20 mM imidazole, 10% *v/v* glycerol, 1 mM DTT, pH 7.4), 10 CV of high-salt buffer (50 mM Tris-HCl, 1 M NaCl, 20 mM imidazole, 1 mM DTT, 10% *v/v* glycerol, pH 7.4) and 10 CV of washing buffer (50 mM Tris-HCl, 400 mM NaCl, 20 mM imidazole, 1 mM DTT, 10% *v/v* glycerol, pH 7.4). The proteins were eluted with 2.5 CV elution buffer (50 mM Tris-HCl, 400 mM NaCl, 250 mM imidazole, 1 mM DTT, 10% *v/v* glycerol, pH 7.4). The eluate was dialysed overnight against PBS (pH 7.4, supplemented with 1 mM DTT). DARPin-Cys were evaluated for the purity by SDS-PAGE (**Figure S4**), aliquoted, flash-frozen in liquid nitrogen and stored at -80 °C until further use.





Figure S4. SDS PAGE analysis of PSA-DARPin-Cys where: (lane 1) molecular weight marker, (lane 2-14) PSA-DARPin-Cys 1 - 13. Protein bands are stained with Coomassie blue.

## Size-exclusion chromatography coupled to multi-angle light scattering (SEC-MALS) detector

The absolute mass of protein samples was determined by using a liquid chromatography system (Agilent LC1100 Agilent Technologies, Santa Clara, USA) coupled to an Optilab rEX refractometer (Wyatt Technology, Santa Barbara, USA) and a miniDAWN three-angle light-scattering detector (Wyatt Technology). For protein separation, a 24 mL Superdex 200 10/30 column (GE Healthcare Biosciences) was run in TBS150 (pH 7.5) at 0.5 mL min<sup>-1</sup>. Analysis was performed by using the ASTRA software (version 5.2.3.15; Wyatt Technology). Note: For some of the smaller DARPin constructs a S75 column might have improved results of these measurements.

#### Molar Mass vs. volum e

1	₩	BSA[Samples]evalu	✓ .	DARPin1[Samples]evalu	✓ .	DARPin2[Samples]evalu	✓ .	DARPin5[Samples]evalu	<b>~</b>	DARPin4[Samples]evalu
	<b>~</b>	DARPin3[Samples]evalu	₩	DARPin6[Samples]evalu	✓	DARPin7[Samples]evalu	<b>~</b>	- DARPin8[Samples]evalu	<b>~</b>	DARPin9[Samples]evalu
	<b>~</b>	DARPin10[Samples]evalu	<b>~</b>	DARPin11[Samples]evalu	<b>~</b>	DARPin12[Samples]evalu	·	<ul> <li>DARPin13[Samples]evalu</li> </ul>		



**Figure S5.** SEC-MALS traces showing the elution of BSA (red dotted line) as a standard and the 12 selected DARPin-Cys constructs with elution times between 17 and 20 min. DARPin-Cys **2** shows higher molecular weight aggregates, and DARPin-Cys **4**, **5** and **7** show peaks corresponding to a lower molecular weight species.

**Table S2.** Summary of the SEC-MALS analysis. The determined size and calculated mass are displayed for each DARPin-Cys construct. BSA was used as a standard.

			Monom	er (Peak 1)			Dimer	(Peak 2)	
	Expected Mw (kDa)	Mw (kDa)	Uncertainty	Calculated mass (µg)	Mass fraction (%)	Mw (kDa)	Uncertainty	Calculated mass (µg)	Mass fraction (%)
BSA	66.4	64.3	0.10%	97.71	93.5	128.2	0.30%	6.74	6.5
DARPin-Cys 1	15.3	18.6	5.30%	6.38	100				
DARPin-Cys 2	15.2	28.4	8.50%	3.49	109.9	223.9	17.00%	-0.31	-9.9
DARPin-Cys 3	15.3	15.4	7.00%	9	100				
DARPin-Cys 4	15.4	6.1	4.50%	20.25	100				
DARPin-Cys 5	15.4	5.7	4.80%	17.84	103.4			-0.59	-3.4
DARPin-Cys 6	18.8	20.3	2.40%	7.85	100				
DARPin-Cys 7	15.4	7.4	3.50%	17.09	100				
DARPin-Cys 8	15.2	16.4	3.70%	7.7	100				
DARPin-Cys 9	11.8	18.8	7.10%	3.43	172			-1.43	-72
DARPin-Cys 10	15.2	19.3	2.80%	5.47	100				
DARPin-Cys 11	15.3	14.6	4.10%	7.33	100				
DARPin-Cys 12	15.3	15.7	2.60%	6.62	100				
DARPin-Cys 13	18.3	18	1.50%	8.15	100				
Average		19.2		15.59	105.6	176.1		1.1	-19.7
Standard deviation		14.3		24.22	19.4	67.7		3.79	35.5
% Standard deviation		74.6		155.3	18.4	38.5		344.38	-180.1
Minimum		5.7		3.43	93.5	128.2		-1.43	-72
Maximum		64.3		97.71	172	223.9		6.74	6.5

## Bioconjugation of DARPin-Cys proteins with a NODAGA-maleimide

## NODAGA-maleimide bioconjugation

The DARPin-Cys proteins (40 nmol) were reduced by incubation with an excess of DTT (8  $\mu$ L, 0.5 M stock in water, 4  $\mu$ mol, 100 eq. see **Table S3**A) for 30 min. (24 °C, 300 rpm). The excess DTT was removed by two subsequent spin filtrations (Zeba desalting column, 7 kDa MWCO, 0.5 mL). The DARPin-Cys proteins (21 to 39 nmol, see **Table S3**B) were then incubated with NODAGA-maleimide (2.14 to 3.92  $\mu$ L, 2 eq., 20 mM in DMSO, see **Table S3**B) in PBS at pH 7.9 for 2 h (24 °C, 300 rpm). To remove excess NODAGA-maleimide, the reactions were purified with NAP-25 columns (Cytiva) pre-equilibrated with the running buffer (PBS, pH 7.4). The dead volume (2.6 mL) was discarded, and the high molecular weight fraction (1.6 mL) was collected and rebuffered into PBS and concentrated by spin-filtration (Amicon ultra 0.5 mL, 3 kDa MWCO, 5,000 g, 3 × 10 min.). Protein concentrations of the purified NODAGA-DARPin conjugates were measured by using a NanoDrop spectrometer, to estimate the recovered protein fractions (21 – 56% yield, see **Table S3**C). Aliquots of the crude reaction mixtures were retained for further radiochemical analysis.

**Table S3.** Table of (A) the utilised reaction conditions for DTT reduction, (B) maleimide conjugation reactions for DARPin-Cys 3, 6, 9, 12, and 13; and (C) The protein concentrations and recovery yields for each NODAGA-DARPin construct.

	DARPin 3	DARPin 6	DARPin 9	DARPin 12	DARPin 13				
(A) DTT treatment									
DARPin conc. / mg mL <sup>-1</sup>	9.76	6.12	5.43	8.50	3.93				
DARPin conc. / μM	640.00	324.88	458.42	553.94	214.77				
DARPin amount / μL	62.50	123.12	87.26	72.21	186.25				
DARPin amount / nmol	40.00	40.00	40.00	40.00	40.00				
DTT (0.5 M) / μL	8.00	8.00	8.00	8.00	8.00				
DTT (0.5 M) / nmol	4000.00	4000.00	4000.00	4000.00	4000.00				
	(B) NODAGA	A-maleimide conjuga	ation after rebufferi	ng					
DARPin conc. / mg mL <sup>-1</sup>	6.05	4.54	3.18	5.91	4.95				
DARPin conc. / μM	396.36	241.16	268.48	384.85	270.69				
DARPin amount / μL	54.00	120.00	90.00	74.00	145.00				
DARPin amount / nmol	21.40	28.94	24.16	28.48	39.25				
NODAGA-maleimide (20uM in DMSO) /μL	2.14	2.89	2.42	2.85	3.92				
NODAGA-maleimide (20uM in DMSO) /nmol	42.81	57.88	48.33	56.96	78.50				
PBS / μL	68.86	2.11	32.58	48.15	0.00				
Conc. DMSO / %	1.71	2.32	1.93	2.28	2.64				
		(C) After purific	ation						
DARPin conc. / mg mL <sup>-1</sup>	2.05	3.59	1.73	2.42	3.15				
DARPin conc. / µM	134.55	190.45	145.79	157.58	172.26				
DARPin amount / μL	92.57	116.80	108.87	54.53	91.40				
DARPin amount / nmol	12.45	22.24	15.87	8.59	15.74				
yield recovered DARPin / %	31.14	55.61	39.68	21.48	39.35				

## SDS-PAGE gel with DARPin-Cys constructs

Purified DARPin-Cys constructs were analysed by SDS-PAGE using the standard protocol. The gel (Bis-Tris, 12%) was run in MES-SDS buffer (see above). DTT was omitted from the loading buffer to analyse the extent of dimerisation *via* the formation of disulfide bridge.


## Synthesis of natGa-NODAGA-DARPin complexes

NODAGA-DARPins ( $\sim$ 0.02 mg,  $\sim$ 1 nmol) were diluted in H<sub>2</sub>O (50  $\mu$ L). Ga(NO<sub>3</sub>)<sub>3</sub> (130 mM in H<sub>2</sub>O, 1  $\mu$ L, 130 eq.) was added. After incubation of the reaction for 15 min. at room temperature, the reaction was analysed by analytical HR-MS. The mass of the non-complexed conjugates was no longer be observed.


The different DARPin-Cys, NODAGA-DARPins and natGa-NODAGA-DARPin complexes were analysed by liquid chromatography high resolution electrospray ionization mass spectrometry (LC-HR-ESI-MS). The obtained results are summarised in **Table S4** and shown in **Figure S6** – **Figure S10**.

**Table S4.** High-resolution mass spectrometry analysis results with the calculated and observed m/z of the five DARPin-Cys, NODAGA-DARPin and <sup>nat</sup>Ga-NODAGA-DARPin. Showing a shift in molecular mass of 597 m/z up on NODAGA conjugation and 66 m/z up on complexation with <sup>nat</sup>Ga. Note: No peaks associated with non-conjugated DARPin-Cys were observed in the measurements for NODAGA-DARPins.


	DARP	in-Cys	NODAGA	A-DARPin	nat GaNODAGA-DARP in		
	Calculated /	Observed /	Calculated /	Observed /	Calculated /	Observed /	
	m/z	m/z	m/z	m/z	m/z	m/z	
DARPin 3	15255.13	15254.79	15752.34	15752.93	15818.24	15820.91	
DARPin 6	18842.18	18841.65	19339.39	19338.83	19405.29	19404.48	
DARPin 9	11854.27	11855.83	12351.48	12351.01	12417.38	12419.13	
DARPin 12	15353.28	15351.82	15850.49	15850.99	15916.39	15915.91	
DARPin 13	18285.34	18287.33	18782.55	18784.44	18848.45	18848.32	




**Figure S6.** Mass-spectrum trace convoluted and de-convoluted HR-MS spectra of DARPin-Cys **3**, NODAGA-DARPin **3** and  $^{nat}$ Ga-NODAGA-DARPin **3**. Showing a shift in molecular mass of 597 m/z upon NODAGA conjugation and 66 m/z upon complexation with  $^{nat}$ Ga.




**Figure S7.** Mass-spectrum trace convoluted and de-convoluted HR-MS spectra of DARPin-Cys **6**, NODAGA-DARPin **6** and  $^{\text{nat}}$ Ga-NODAGA-DARPin **6**. Showing a shift in molecular mass of 597 m/z upon NODAGA conjugation and 66 m/z upon complexation with  $^{\text{nat}}$ Ga.



**Figure S8.** Mass-spectrum trace convoluted and de-convoluted HR-MS spectra of DARPin-Cys **9**, NODAGA-DARPin **9** and  $^{nat}$ Ga-NODAGA-DARPin **9**. Showing a shift in molecular mass of 597 m/z upon NODAGA conjugation and 66 m/z upon complexation with  $^{nat}$ Ga.



**Figure S9.** Mass-spectrum trace convoluted and de-convoluted HR-MS spectra of DARPin-Cys **12**, NODAGA-DARPin **12** and <sup>nat</sup>Ga-NODAGA-DARPin **12**. Showing a shift in molecular mass of 597 *m/z* upon NODAGA conjugation and 66 *m/z* upon complexation with <sup>nat</sup>Ga.



**Figure S10.** Mass-spectrum trace convoluted and de-convoluted HR-MS spectra of DARPin-Cys **13**, NODAGA-DARPin **13** and <sup>nat</sup>Ga-NODAGA-DARPin **13**. Showing a shift in molecular mass of 597 *m/z* upon NODAGA conjugation and 66 *m/z* upon complexation with <sup>nat</sup>Ga.

## Qualitative assessment of PSA binding of NODAGA-DARPin constructs

Evaluation of PSA binding of DARPin constructs after NODAGA modification by automated SEC chromatography

The purified NODAGA-DARPin conjugates (7.5  $\mu$ L, 125  $\mu$ M, 937.5 pmol) were incubated with PSA (8  $\mu$ L, 937.5 pmol, 1 eq.) at 24 °C (30 min., 400 rpm). These reaction mixtures (10  $\mu$ L) were then analysed by using analytical SEC (method A). In the chromatograms, higher molecular weight complexes were observed at shorter retention times for all four PSA binding DARPins (3, 6, 9, and 12), but not for the negative control (DARPin 13).




Figure S11. SEC analysis of PSA binding. The NODAGA-DARPins 3, 6, 9, 12, and 13 (panels A – E, cyan) were incubated with 1 eq. of PSA (panels A – E, black) leading to the formation of a higher molecular weight complex with slightly shorter retention times (panels A – E, purple), except for the non-binding control DARPin 13, where only a mixture of both proteins was detected and no complex formation is expected (theoretical overlay of PSA and DARPin 13 trace in pink dotted line). (F) Table of retention times of the NODAGA-DARPins and the NODAGA-DARPin-PSA complexes. Note: The absorption was measured at 220 nm.

## Native-PAGE gel analysis of PSA-DARPin-Cys complexes

The purified NODAGA-DARPin conjugates (1.5  $\mu$ L, 125  $\mu$ M, 187.5 pmol) were incubated with PSA (1.6  $\mu$ L, 187.5 pmol, 1 eq.) at 24 °C (30 min., 400 rpm). These reaction mixtures were then mixed with native loading buffer (Native-TrisGly Sample Buffer 2×, NOVEX Life Technologies), and the binding of the NODAGA-DARPins to PSA was evaluated by using native PAGE, as described above.

## Surface plasmon resonance measurements for quantitative assessment of PSA binding of DARPin-Cys and NODAGA-DARPins

The measurements and analyses were performed at the Functional Genomics Center Zürich (FGCZ)

SPR measurements were performed at 20 °C in HBS buffer (10 mM HEPES, 150 mM NaCl, 0.005% Tween-20, 50  $\mu$ M EDTA, pH 7.6) using a Biacore T200/S200 instrument (GE Healthcare/Cytiva). Proteins were immobilised on a multidentate poly-nitrilotriacetic acid (NTA)-derivatised linear polycarboxylate hydrogel NiHC30M chip, according to a protocol recommended by the chip manufacturer (Xantec, Düsseldorf, Germany). Firstly, a 0.35 M EDTA solution (pH 8.5) was injected for 300 s, followed by a 120 s buffer injection, a surface activation step with 5 mM NiCl<sub>2</sub> (in HBS, 20  $\mu$ M EDTA) for 60 – 120 s, and a 40 – 60 s injection of DARPin-Cys (200 nM) or NODAGA-DARPins (125 nM), respectively, both containing a His<sub>6</sub> tag. This was followed by a 120 s injection of immobilization buffer. All immobilization steps were performed in immobilization buffer (10 mM HEPES, 150 mM NaCl, 0.005% Tween-20, and 50  $\mu$ M EDTA) at a flow rate of 5  $\mu$ L min<sup>-1</sup>. This results in a surface density of 35 – 140 RU.

In the double-referenced binding experiments, a two-fold dilution series of five concentrations of PSA, PSA-ACT and ACT were injected for 120 s, at a flow rate of 30  $\mu$ L min<sup>-1</sup> in HBS buffer. PSA concentrations were in the range of 69 – 142 nM, PSA-ACT concentrations in the range 1 – 16 nM, and ACT (from human plasma, Sigma-Aldrich) concentration in the range 1 – 16 nM, respectively. The chip surface was regenerated by injection of 0.35 M EDTA for 300 s. A concentration of 0.4 mg mL<sup>-1</sup> for the PSA-ACT complex stock solution (average of the concentration provided by the supplier: 0.3 – 0.5 mg mL<sup>-1</sup>, Prostate Specific Antigen– $\alpha$ 1-Antichymotrypsin Complex human, Sigma-Aldrich) and a molecular weight of 60 kDa were assumed.


The DARPin 13 coated flow cell 1 and 3 of the sensor chip were used as a reference. The use of DARPin-Cys 13 as a control stabilised the baseline and reduced non-specific binding of the analyte to the chip. Data were evaluated by using Biacore software version 2.0.3. Due to unsatisfactory fits with a 1:1 model, most sensorgrams were fitted by using a sum of two exponentials (Sum2Exp, in BiaEvaluate v2.03 called heterogeneous ligand) kinetic fit. All obtained values are displayed in **Table S5** and **Table S6**; all measured sensorgrams are displayed in **Figure S12** – **Figure S15**.

**Table S5.** Results of single-cycle kinetic surface plasmon resonance (SPR) measurements of the DARPins-Cys 3, 6, 9, and 12. Determined rate constant of association  $(k_{on} / M^{-1} s^{-1})$ , rate constant of dissociation  $(k_{off} / s^{-1})$  and the dissociation constant  $(K_D / M^{-1})$  are given for each construct. For all ACT measurements no binding was detected (not detected, *n.d.*)


Probe	Probe Density	Temp/°C	Buffer	Analyte	Model	$k_{on} / 10^5$ $M^{-1} s^{-1}$	k <sub>off</sub> / 10 <sup>-3</sup> s-1	K <sub>D</sub> / nM	RU(max)	Chi <sup>2</sup> / RU <sup>2</sup>					
				PSA	Sum2Exp	$1.515 \pm 0.003$	$2.440 \pm 0.001$	16.11 ± 0.06	117.1 ± 0.3	0.207					
DARPin-Cys 3	70	20	TIDG	HBS	HBS	HBS	HBS	HBS	rsA	Sunzexp	$7.137 \pm 0.065$	$12.92 \pm 0.005$	$18.10 \pm 0.18$	$55.92 \pm 0.37$	0.207
DARPIn-Cys 3	/0	20	20						PSA-ACT	1:1	$6.966 \pm 0.132$	$44.37 \pm 0.30$	63.69 ± 1.28	$8.509 \pm 0.115$	0.0259
				ACT			n.	d.							
		20 HBS		PSA	C2E	$8.646 \pm 0.038$	$13.88 \pm 0.045$	$16.06 \pm 0.09$	34.67 ± 0.11	0.0659					
DADDin Con 6	112		HBS	HDC	PSA	Sum2Exp	$1.201 \pm 0.003$	$2.498 \pm 0.004$	$20.79 \pm 0.06$	$61.15 \pm 0.09$	0.0639				
DARPin-Cys 6	112		20	10 HBS	HBS	PSA-ACT	1:1	$6.711 \pm 0.083$	$18.27 \pm 0.238$	$27.22 \pm 0.49$	$13.43 \pm 0.05$	0.0387			
				ACT	T n.d.										
				PSA	G 2F	$1.104 \pm 0.006$	$0.1460 \pm 0.0064$	$1.322 \pm 0.059$	$70.29 \pm 0.27$	0.440					
								PSA	Sum2Exp	$8.924 \pm 0.632$	$2.484 \pm 0.007$	$2.783 \pm 0.021$	105.3 ± 0.3	0.440	
DARPin-Cys 9	110	20	HBS	DC 4 CT	G 2F	$6.308 \pm 0.096$	$4.400 \pm 0.047$	$6.976 \pm 0.130$	$1.556 \pm 0.016$	0.0302					
				PSA-ACT	Sum2Exp	$10.57 \pm 0.070$	$53.00 \pm 0.26$	50.13 ± 0.41	$11.15 \pm 0.04$	0.0302					
				ACT			n.	d.							
				DC 4	G 2F	$9.484 \pm 0.063$	$3.022 \pm 0.009$	$3.186 \pm 0.023$	$41.51 \pm 0.09$	0.0600					
				PSA	Sum2Exp	$0.8821 \pm 0.0043$	$0.2977 \pm 0.0057$	$3.375 \pm 0.067$	$33.91 \pm 0.08$	0.0690					
DARPin-Cys 12	90	20		HBS PSA-ACT		$7.961 \pm 0.029$	$2.281 \pm 0.013$	$2.865 \pm 0.020$	$25.27 \pm 0.15$	0.245					
					Sum2Exp	$4.084 \pm 0.035$	12.61 ± 0.13	30.87 ± 0.41	$20.38 \pm 0.12$	0.345					
				ACT			n.	d.							

**Table S6.** Results of single-cycle kinetic surface plasmon resonance (SPR) measurements of the NODAGA-DARPins 3, 6, 9, and 12. Determined rate constant of association ( $k_{on} / M^{-1} s^{-1}$ ), rate constant of dissociation ( $k_{off} / s^{-1}$ ) and the dissociation constant ( $K_D / M^{-1}$ ) are given for each construct. For all ACT measurements no binding was detected (not detected, n.d.)


Probe	Probe Density	Temp/°C	Buffer	Analyte	Model	k <sub>on</sub> / 10 <sup>5</sup> M <sup>-1</sup> s <sup>-1</sup>	k <sub>off</sub> / 10 <sup>-3</sup> s-1	K <sub>D</sub> / nM	RU(max)	Chi <sup>2</sup> / RU <sup>2</sup>	
				PSA	C2E	8.824 ± 0.417	$0.8728 \pm 0.0073$	$0.9891 \pm 0.009$	$53.73 \pm 0.35$	0.715	
				PSA	Sum2Exp	$2.243 \pm 0.034$	4.046 ± 0.043	18.04 ± 0.33	$25.98 \pm 0.30$	0./15	
NODAGA- DARPin 3	49	20	HBS	PSA-ACT	Sum2Exp	$0.2488 \pm 0.0089$	$0.01391 \pm 0.00079$	$0.5691 \pm 0.0037$	61.21 ± 1.93	2.53	
D. 114 11.0				PSA-ACT	SumzExp	$22.50 \pm 0.18$	$2.923 \pm 0.029$	$1.299 \pm 0.014$	$108.8 \pm 0.1$	2.53	
				ACT			d.				
				PSA	C2E	$5.399 \pm 0.027$	$0.9626 \pm 0.0062$	$1.783 \pm 0.015$	$67.94 \pm 0.30$	0.618	
				IJA	Sum2Exp	1.432 ± 0.041	$6.184 \pm 0.130$	43.18 ± 1.54	$16.40 \pm 0.20$	0.618	
NODAGA- DARPin 6	54	20	HBS	DC 4 ACT	C2E	$10.26 \pm 0.076$	$0.9818 \pm 0.0039$	$0.9566 \pm 0.0080$	$90.69 \pm 0.2378$	1.05	
<i>D.</i> <b>III. II. V</b>				PSA-ACT	Sum2Exp	$1.033 \pm 0.037$	$6.542 \pm 0.089$	63.33 ± 2.40	$46.38 \pm 0.95$	1.05	
				ACT			n.	d.			
				PSA	Sum2Exp	$46.95 \pm 0.39$	$16.04 \pm 0.098$	$3.416 \pm 0.036$	$31.94 \pm 0.07$	0.295	
Non. G.				PSA	SumzExp	$2.428 \pm 0.010$	$2.536 \pm 0.011$	$10.44 \pm 0.06$	$29.05 \pm 0.07$	0.295	
NODAGA- DARPin 9	44	20	HBS	PSA-ACT	C2E	$2.028 \pm 0.006$	$0.07143 \pm 0.00173$	$0.3523 \pm 0.0086$	$133.1 \pm 0.1$	0.0515	
				PSA-ACT	Sum2Exp	$20.21 \pm 0.15$	$2.049 \pm 0.005$	$1.014 \pm 0.008$	$60.64 \pm 0.10$	0.0515	
				ACT			n.	d.			
				PSA	Sum2Exp	$16.05\pm0.08$	$9.007 \pm 0.030$	$5.613 \pm 0.039$	$28.05 \pm 0.09$	0.214	
				ISA	sA Sum2Exp	$1.683 \pm 0.006$	$1.438 \pm 0.005$	$8.543 \pm 0.045$	$34.69 \pm 0.09$	0.214	
NODAGA- DARPin 12	62	20	HBS	HBS	DCA ACT	NGA AGT G 25	$0.8394 \pm 0.0080$	$0.6479 \pm 0.0018$	$7.706 \pm 0.077$	$129.2 \pm 0.8$	0.0580
				PSA-ACT	Sum2Exp	$5.804 \pm 0.174$	$12.35 \pm 0.11$	$21.29 \pm 0.66$	$50.83 \pm 0.29$	0.0380	
				ACT			n.	d.			



**Figure S12.** Surface plasmon resonance (SPR) sensorgrams of DARPin-Cys **3**, **6**, **9**, **12** indicating binding to the analyte PSA. Experimental data are displayed in cyan; the black dashed line shows Sum2Exp fit.



**Figure S13.** Surface plasmon resonance (SPR) sensorgrams of DARPin-Cys **3**, **6**, **9**, **12** indicating binding to the analyte PSA-ACT. Experimental data are displayed in cyan; the black dashed line shows Sum2Exp fit.  $^{\Delta}$  Data were fitted with a 1:1 model.



**Figure S14.** Surface plasmon resonance (SPR) sensorgrams of NODAGA-DARPin 3, 6, 9, 12 indicating binding to the analyte PSA. Experimental data are displayed in cyan; the black dashed line shows Sum2Exp fit.




Figure S15. Surface plasmon resonance (SPR) sensorgrams of NODAGA-DARPin 3, 6, 9, 12 indicating a binding to the analyte PSA-ACT. Experimental data are displayed in cyan; the black dashed line shows Sum2Exp fit.

## Radiolabelling of NODAGA-DARPins

## <sup>68</sup>Ga-radiolabelling of NODAGA-DARPins


An excess of NODAGA-DARPins (3.0  $\mu$ L, ~6  $\mu$ g, ~0.40 nmol) was buffered in NaOAc (50  $\mu$ L, 0.1 M, pH 4.4). Radiolabelling was achieved by addition of an aliquot of [ $^{68}$ Ga][Ga(H<sub>2</sub>O)<sub>6</sub>]Cl<sub>3</sub> (8.5  $\mu$ L, ~4 MBq) and adjustment of the pH to 5.3 – 5.5 with Na<sub>2</sub>CO<sub>3</sub> (1.5  $\mu$ L, 1 M). After incubation of the reaction for 10 min. at room temperature, the [ $^{68}$ Ga]GaNODAGA-DARPin complexes were analysed by radio-TLC, manual NAP-25 and automated radio-SEC (**Figure S16**).

## Molar activities of <sup>68</sup>Ga-radiollabelling of [<sup>68</sup>Ga]GaNODAGA-DARPins

Varying amounts of NODAGA-DARPins (0.005-0.500 nmol) were radiolabelled with the same amounts of radioactivity ([ $^{68}$ Ga][Ga( $H_2O$ ) $_6$ ]Cl $_3$  [ $8.5~\mu$ L,  $\sim$ 4 MBq]) by using the standard procedure described above. After incubation of the reaction mixtures for 10 min. at room temperature, the reaction was analysed by automated radio-SEC (method B) and the RCPs were determined by integration of the peaks. The maximum molar activities were determined to be 30.72 to 54.37 MBq nmol $^{-1}$  of protein (**Table S7**).

**Table S7.** Summary of molar activity values in MBq nmol<sup>-1</sup> obtained by titrations of different amounts of the five NODAGA-DARPins with a constant amount of  $[^{68}Ga][Ga(H_2O)_6]Cl_3$ .

[68Ga]GaNODAGA-DARPin	3	6	9	12	13
Specific activity / MBa nmol <sup>-1</sup>	36.21	37.16	30.72	54.37	39.82



**Figure S16.** (Panels A, D, G, J, and M) Radio-TLC analysis data showing characterisation of [<sup>68</sup>Ga]GaNODAGA-DARPins **3**, **6**, **9**, **12**, and **13** after purification (cyan), as a crude (green) in citrate buffer. The purple and black traces show controls with [<sup>68</sup>Ga]Ga-NODAGA-maleimide, and [<sup>68</sup>Ga]Ga-citrate. (Panels B, E, H, K, and N) Manual radio-NAP-25 SEC analysis of the purified (cyan) and the crude (green) constructs, and [<sup>68</sup>Ga]Ga-NODAGA-maleimide (purple). (Panels C, F, I, L, and O) Automated radio-SEC analysis showing co-elution of the NODAGA-DARPins **3**, **6**, **9**, **12**, and **13** UV trace (light blue) with the [<sup>68</sup>Ga]GaNODAGA-DARPins **3**, **6**, **9**, **12**, and **13** radiotrace (cyan). The crude (green) construct, [<sup>68</sup>Ga]Ga-NODAGA-maleimide (purple), and [<sup>68</sup>Ga]Ga-EDTA control complex (black).

## Stability of [68Ga]GaNODAGA-DARPins

The stability of the five [<sup>68</sup>Ga]GaNODAGA-DARPins was investigated with respect to change in radiochemical purity (RCP). The results are summarised below (**Figure S17**).

#### Incubation in PBS

Radiolabelled [ $^{68}$ Ga]GaNODAGA-PSA-DARPins (30  $\mu$ L,  $\sim$ 3  $\mu$ g,  $\sim$ 0.20 nmol,  $\sim$ 2 MBq) were incubated with PBS (pH 7.4, 60  $\mu$ L) at 37 °C. Aliquots of the reactions were analysed by automated radio-SEC (method B) at time intervals of 0, 30, 60, 90, and 120 min. The RCP are summarised in **Table S8**.

**Table S8.** Table showing the RCP of [<sup>68</sup>Ga]GaNODAGA-DARPins over time after incubation in PBS at 37 °C as determined by automated radio-SEC.

Time /	[ <sup>68</sup> Ga]GaNODAGA-				
min.	DARPin 3 RCP / %	DARPin 6 RCP / %	DARPin 9 RCP / %	DARPin 12 RCP / %	DARPin 13 RCP / %
0	100	100	100	100	100
30	100	100	100	100	100
60	99	100	100	98	100
90	99	99	100	97	100
120	98	99	100	99	100

#### Incubation with excess EDTA

Radiolabelled [ $^{68}$ Ga]GaNODAGA-PSA-DARPins (30  $\mu$ L,  $\sim$ 3  $\mu$ g,  $\sim$ 0.20 nmol,  $\sim$ 2 MBq) were incubated with EDTA (pH 7.1, 30  $\mu$ L, 50 mM, 1.5  $\mu$ mol, 3,700 fold excess) at room temperature, and the RCP was monitored by radio-iTLC at 0, 30, 60, 90, and 120 min. post EDTA addition. The RCPs are summarised in **Table S9**.

Table S9. Table showing the RCP of [68Ga]GaNODAGA-DARPins over time after EDTA addition (evaluated by radio-TLC).


Time / min.	[ <sup>68</sup> Ga]GaNODAGA-				
mın.	DARPin 3 RCP / %	DARPin 6 RCP / %	DARPin 9 RCP / %	DARPin 12 RCP / %	DARPin 13 RCP / %
0	100	100	100	100	100
30	100	98	97	99	98
60	97	98	95	97	99
90	95	98	92	97	100
120	95	96	92	96	99

#### Incubation in human serum

Radiolabelled [ $^{68}$ Ga]GaNODAGA-PSA-DARPins (40  $\mu$ L,  $\sim$ 4  $\mu$ g,  $\sim$ 0.27 nmol,  $\sim$ 2.7 MBq) were incubated with human serum (80  $\mu$ L) at 37 °C, pH 7.4. Aliquots of the reactions were analysed by automated radio-SEC (method B) at time intervals of 0, 30, 60, 90, and 120 min. The RCP are summarised in **Table S10**.

**Table S10.** Table showing the RCP of [<sup>68</sup>Ga]GaNODAGA-DARPins over time after incubation with human serum at 37 °C, (evaluated by automated radio-SEC).

Time /	[ <sup>68</sup> Ga]GaNODAGA-				
min.	DARPin 3 RCP / %	DARPin 6 RCP / %	DARPin 9 RCP / %	DARPin 12 RCP / %	DARPin 13 RCP / %
0	100	100	100	100	100
30	96	98	97	99	100
60	97	98	97	99	94
90	99	99	97	100	99
120	99	99	97	100	99



**Figure S17**. (Panels A, D, G, J and M) Automated radio-SEC chromatograms following the change in RCP of [68Ga]GaNODAGA-DARPins over time after PBS (pH 7.4) addition at 37 °C. (Panels B, E, H, K and N) Radio-TLC chromatograms developed with citrate eluent (0.2 M, pH 4.5) following the change in RCP of [68Ga]GaNODAGA-DARPins over time after EDTA addition at 24 °C, pH 7.1. (Panels C, F, I, L and O) Automated radio-SEC following the change in RCP of [68Ga]GaNODAGA-DARPins over time after human serum addition at 37 °C, pH 7.4.

## Radioactive PSA-binding assay

Since PSA is a secreted protein, the immunoreactivity could not be determined by using a cellular-based procedure adapted from Lindmo and co-workers.<sup>13</sup> Instead, a radioactive immune capture-type assay was developed to confirm binding.

[68Ga]GaNODAGA-PSA-DARPins were synthesised by using the standard procedure stated above (reaction volume 60 μL). The reaction was diluted (1 mL, PBS + 0.3% w/v BSA), and the diluted radiotracer [68Ga]GaNODAGA-PSA-DARPin (6 – 9 μL, 2 pmol, ~26 kBq) was pre-incubated with biotinylated-PSA (85 μg mL<sup>-1</sup>, 2.9 μM, 6.8 μL, 20 pmol, 10 eq.) at 4 °C (300 rpm, 30 min.). For blocking experiments, non-radiolabelled DARPin-Cys (250 eq. 1 – 2 μL, 500 pmol) were added to biotin-PSA 15 min before the addition of the radiotracer. As a second control, biotin-PSA was omitted from the assay and substituted with PBS. Three samples of the same specified quantity of radiotracer were also prepared to serve as standards for total activity added to each sample. Then streptavidin-coated magnetic beads (5 μL, washed 2 × with PBS + w/v 0.3% BSA prior to use) were added and the reactions were incubated at 4 °C (10 min., 800 rpm). The magnetic beads were then pulled down and washed 3 times (0.5 mL, PBS + w/v 0.3% BSA), with re-pulling down and the supernatant was removed between each wash. Bead-associated radioactivity of the washed pull down was determined by using a gamma counter and the relative bead-bound fraction was determined relative to the loading standard. All reactions were performed in triplicate and the mean was determined. Results are summarised in **Table S11**.

**Table S11.** Relative associated activities (RAA) in % of the standard loading dose for all the evaluated [68Ga]GaNODAGA-DARPins.

	[ <sup>68</sup> Ga]GaNODAGA-				
	DARPin 3	DARPin 6	DARPin 9	DARPin 12	DARPin 13
	RAA/ %	RAA/ %	RAA/%	RAA/ %	RAA/ %
binding	$47.61 \pm 1.27$	$37.03 \pm 0.52$	$38.60 \pm 1.05$	$39.64 \pm 1.35$	$4.75 \pm 0.99$
no PSA	$7.10 \pm 2.26$	$4.72\pm0.88$	$11.81 \pm 3.42$	$13.23 \pm 0.51$	$3.33 \pm 1.77$
block	$5.56\pm0.88$	$3.73 \pm 1.28$	$11.67 \pm 1.36$	$13.82 \pm 0.87$	$2.51 \pm 0.52$

# Tumor uptake studies and *ex vivo* biodistribution studies of [<sup>68</sup>Ga]GaNODAGA-DARPin radiotracers in LNCaP xenografts

NODAGA-DARPins (~1.30 nmol) were buffered in NaOAc (100  $\mu$ L, 0.1 M, pH 4.4). Radiolabelling was achieved by addition of an aliquot of [ $^{68}$ Ga][Ga(H<sub>2</sub>O)<sub>6</sub>]Cl<sub>3</sub> (50  $\mu$ L) and adjustment of the pH to >5.3 with Na<sub>2</sub>CO<sub>3</sub> (10  $\mu$ L, 1M). After incubation of the reaction for 10 min. at room temperature, the [ $^{68}$ Ga]GaNODAGA-DARPin complexes were analysed by radio-TLC and automated radio-SEC. The crude reactions were purified *via* NAP-5 SEC chromatography to obtain the purified radiotracer in >99% RCP in PBS (pH 7.4) with decay-corrected RCYs of 57 – 78% (n = 2). Tumour uptake studies were conducted in mouse models with LNCaP xenografts. Tumourbearing mice were randomised before the study with n = 3 – 4, animals per group. Doses of [ $^{69}$ Ga]GaNODAGA-DARPin were administered *via* tail-vein injection (150  $\mu$ L in sterile PBS, amounts **Table S12**). A competitive blocking study was performed with a lower molar activity [ $^{69}$ Ga]GaNODAGA-DARPin, produced by supplementing the normal formulation with non-labelled DARPin (blocking group, 150  $\mu$ L in sterile PBS amounts given in **Table S12**).

**Table S12.** Summary of the decay corrected radiochemical yields (d.c. RCY / %) for the isolated, formulated radiotracers after NAP-5 purification (n = 2) and the amounts ( $\mu$ g, nmol, and MBq), and the molar activity (MBq nmol<sup>-1</sup>) of the administered radiotracers. Values are displayed as geometric mean  $\pm$  one standard deviation.

	[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>3</b>		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>6</b>		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>9</b>		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>12</b>		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>13</b>	
D.c. RCY / % (n = 2)	77.86 ± 1.10		$71.03 \pm 8.32$		$56.67 \pm 3.90$		$63.89 \pm 2.85$		$70.60 \pm 1.13$	
	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 3)$	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 4)$	Normal group $(n = 4)$	-
Amount / μg	3.00	598.96	3.84	601.14	2.37	322.44	3.03	392.17	3.74	-
Amount / nmol	0.21	31.79	0.20	31.90	0.20	27.84	0.20	20.81	0.20	-
Activity / MBq	$2.88 \pm 0.29$	$2.66\pm0.27$	$2.40\pm0.17$	$2.70\pm0.17$	1.56 ± 0.18	1.875 ± 0.158	$1.70\pm0.18$	$2.40 \pm 0.25$	$2.18\pm0.26$	-
Molar. act. / MBq nmoΓ <sup>1</sup>	14.02 ± 1.40	$0.084 \pm 0.009$	$11.79 \pm 0.84$	$0.085 \pm 0.005$	$7.60 \pm 0.86$	0.067 ± 0.006	$8.36 \pm 0.88$	$0.116 \pm 0.012$	$10.64 \pm 1.25$	-

#### Effective half-life measurement

The total internal radioactivity over time was measured using a dose calibrator. With this measurement, the effective half-lives ( $t_{1/2}(eff)$ ) of [ $^{68}$ Ga]GaNODAGA-DARPins were calculated by using a one-phase decay model in GraphPad Prism.

**Table S13.** Total internal radioactivity relative to internal activity at t = 5, 30, 60 and 120 min. of [ $^{68}$ Ga]GaNODAGA-DARPins in athymic nude mice bearing s.c. LNCaP xenografts and the calculated  $t_{1/2}$  (eff) values obtained from a one-phase decay model. Note: Data were analysed in GraphPad Prism.

	[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>3</b> Relative internal activity		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>6</b> Relative internal activity		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>9</b> Relative internal activity		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>12</b> Relative internal activity		[ <sup>68</sup> Ga]GaNODAGA- DARPin <b>13</b> Relative internal activity	
Time / min.	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 3)$	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 4)$	Normal group $(n = 4)$	-
5	1.00	1.00	1.00	1	1	1	1	1	1	-
30	$0.71 \pm 0.03$	$0.74 \pm 0.13$	$0.69 \pm 0.02$	$0.73\pm0.09$	$0.74 \pm 0.02$	$0.65 \pm 0.06$	$0.70 \pm 0.06$	$0.74 \pm 0.11$	$0.83 \pm 0.05$	-
60	$0.57 \pm 0.04$	$0.49 \pm 0.06$	$0.48 \pm 0.04$	$0.51 \pm 0.07$	$0.51 \pm 0.01$	$0.50 \pm 0.03$	$0.42 \pm 0.02$	$0.51 \pm 0.05$	$0.59 \pm 0.02$	-
120	$0.28 \pm 0.02$	$0.27 \pm 0.05$	$0.27 \pm 0.04$	$0.29 \pm 0.05$	$0.27 \pm 0.01$	$0.26 \pm 0.01$	$0.22 \pm 0.01$	$0.24 \pm 0.02$	$0.31 \pm 0.02$	-
t 1/2 (eff) / min.	61.26 ± 12.79	55.92 ± 22.02	39.21 ± 4.75	51.98 ± 16.95	52.48 ± 3.61	$38.65 \pm 7.38$	$40.04 \pm 5.55$	64.55 ± 20.24	112.59 ± 32.29	-

## Ex vivo biodistribution studies

Animals were anesthetised by inhalation of 2-3% isoflurane (Baxter Healthcare, Deerfield, IL) / oxygen gas mixture ( $\sim$ 5 L min<sup>-1</sup>) and were euthanised by exsanguination under anaesthesia. Tissue samples were removed, rinsed in water, dried in air for  $\sim$ 1 min., weighed and counted on a gamma-counter for accumulation of <sup>69</sup>Garadioactivity. Tumour uptake values (% ID g<sup>-1</sup>) are given in **Table S14**. Tumour-to-tissue ratios are given in **Table S15**. The statistical analysis of the comparisons of the different groups is represented in **Table S16** – **Table S18** (P-values of unpaired, two-tailed Student's t-test).

**Table S14.** Biodistribution data showing the accumulation of [ $^{68}$ Ga]GaNODAGA-DARPins **3**, **6**, **9**, **12**, and **13** in different tissues at 2 h post-administration in the normal and blocking groups of athymic nude mice bearing s.c. LNCaP xenografts mice. The tissue associated uptake is represented in %ID  $g^{-1}$ . Note: S.D. = 1 standard deviation. Grey shading indicates that n = 3, as one value had to be excluded from the experiment.

	[ <sup>68</sup> Ga]GaNODA	GA-DARPin 3	[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>6</b>		[ <sup>68</sup> Ga]GaNOD	[ <sup>68</sup> Ga]GaNODAGA-DARPin 9		[ <sup>68</sup> Ga]GaNODAGA-DARPin 12		[ <sup>68</sup> Ga]GaNODAGA-DARPin 13	
	(2 h), %IE	$g^{-1} \pm S.D$	(2 h), $\%ID g^{-1} \pm S.D$		$(2 h)$ , $\%ID g^{-1} \pm S.D$		(2 h), %ID g <sup>-1</sup> ± S.D		(2 h), %ID	$(2 h)$ , $\%ID g^{-1} \pm S.D$	
Tissue	Normal group	Blocking group	Normal group	Blocking group	Normal group	Blocking group	Normal group	Blocking group	Normal group	_	
1 100 000	(n = 4)	(n = 3)	(n = 3)	(n = 3)	(n = 4)	(n = 3)	(n = 3)	(n = 4)	(n = 4)		
Blood	$1.71 \pm 0.16$	$1.33\pm0.25$	$7.13 \pm 1.25$	$2.69 \pm 0.13$	$6.74 \pm 1.14$	$4.56\pm0.04$	$2.10 \pm 0.42$	$0.92 \pm 0.11$	$5.53 \pm 0.64$	-	
Tumour	$2.39 \pm 0.43$	$1.48\pm0.19$	$4.16\pm0.58$	$2.47\pm0.42$	$4.72\pm1.07$	$2.9\pm0.46$	$1.80 \pm 0.59$	$1.58 \pm 0.38$	$2.93 \pm 0.55$	-	
Heart	$0.91 \pm 0.07$	$0.67\pm0.13$	$2.81 \pm 0.81$	$1.94\pm0.08$	$2.01\pm1.07$	$1.98\pm0.23$	$1.20 \pm 0.24$	$0.70\pm0.05$	$1.94 \pm 0.11$	-	
Lungs	$1.28\pm0.12$	$0.91\pm0.18$	$3.69 \pm 0.87$	$2.24\pm0.18$	$4.18\pm0.67$	$2.98\pm0.20$	$1.59 \pm 0.54$	$0.71 \pm 0.14$	$3.90 \pm 1.10$	-	
Liver	$21.31\pm3.24$	$14.25\pm2.74$	$33.16\pm7.93$	$29.66 \pm 5.56$	$32.65 \pm 2.31$	$29.02 \pm 4.78$	$20.41 \pm 3.99$	$10.67 \pm 1.44$	$43.04 \pm 3.41$	-	
Spleen	$6.85\pm0.99$	$3.58\pm0.86$	$18.07 \pm 7.10$	$8.57 \pm 1.13$	$12.56 \pm 2.91$	$21.62 \pm 1.35$	$12.89 \pm 5.24$	$2.68 \pm 0.85$	$14.99 \pm 2.40$	-	
Stomach	$0.49\pm0.04$	$0.51\pm0.06$	$0.76 \pm 0.09$	$0.73\pm0.06$	$0.54\pm0.21$	$0.95\pm0.53$	$0.66 \pm 0.14$	$0.67 \pm 0.32$	$1.31 \pm 0.57$	-	
Pancreas	$0.85 \pm 0.08$	$0.57\pm0.09$	$1.80\pm0.34$	$1.64 \pm 0.04$	$1.74\pm0.08$	$1.34 \pm 0.18$	$1.39 \pm 0.31$	$0.58 \pm 0.10$	$1.17 \pm 0.66$	-	
Kidney	$139.28 \pm 65.59$	$116.50 \pm 19.07$	$4.44 \pm 0.66$	$148.80 \pm 25.75$	$36.13\pm4.28$	$80.85 \pm 10.27$	149.63 ± 25.00	$139.03 \pm 23.11$	$35.66 \pm 5.02$	-	
Sm. Int.	$1.03\pm0.36$	$1.65\pm1.03$	$2.74\pm1.09$	$1.47 \pm 0.11$	$2.12\pm0.67$	$2.75 \pm 1.68$	$1.31 \pm 0.39$	$1.70 \pm 1.22$	$3.05 \pm 1.85$	-	
Large	$0.60\pm0.04$	$0.59\pm0.15$	$1.08\pm0.09$	$1.19 \pm 0.08$	$1.47\pm0.43$	$1.16 \pm 0.08$	$0.70 \pm 0.22$	$0.54 \pm 0.18$	$1.55 \pm 0.59$	-	
Fat	$0.65 \pm 0.40$	$0.63\pm0.38$	$0.78 \pm 0.14$	$1.48\pm0.88$	$0.89 \pm 0.21$	$1.05 \pm 0.13$	$0.46 \pm 0.20$	$0.47\pm0.17$	$0.77 \pm 0.06$	-	
Muscle	$0.37 \pm 0.03$	$0.30 \pm 0.11$	$0.90\pm0.10$	$0.88 \pm 0.22$	$0.82\pm0.09$	$0.72 \pm 0.19$	$0.49 \pm 0.10$	$0.33 \pm 0.08$	$0.70 \pm 0.17$	-	
Bone	$1.26\pm0.17$	$0.83\pm0.12$	$2.64 \pm 0.81$	$2.32\pm0.73$	$3.15\pm0.64$	$2.77 \pm 1.43$	$2.35 \pm 0.39$	$1.11\pm0.56$	$2.32 \pm 0.94$	-	
Skin	$1.76\pm0.31$	$1.29 \pm 0.05$	$2.35\pm0.43$	$2.55 \pm 0.04$	$2.05\pm0.27$	$1.80 \pm 0.12$	$2.29 \pm 0.26$	$1.98 \pm 0.63$	$1.74 \pm 0.10$	-	

**Table S15.** Tumour-to-tissue contrast data following the biodistribution of  $[^{68}Ga]GaNODAGA$ -DARPins 3, 6, 9, 12, and 13 in the normal and blocking groups in athymic nude mice bearing s.c. LNCaP xenografts. Note: S.D. = 1 standard deviation, grey: n = 3, as one value had to be excluded.

	$\begin{bmatrix} ^{68}Ga]GaNODAGA-DARPin \ \textbf{3} \\ (2 \ h), tumour to tissue contrast \\ ratio \pm S.D \\ \end{bmatrix}$		[ $^{68}$ Ga]GaNODAGA-DARPin <b>6</b> (2 h), tumour to tissue contrast ratio $\pm$ S.D		$[^{68}Ga]GaNODAGA-DARPin$ 9 (2 h), tumour to tissue contrast ratio $\pm$ S.D		[68Ga]GaNODAGA-DARPin 12 (2 h), tumour to tissue contrast ratio ± S.D		$[^{68}$ Ga]GaNODAGA-DARPin 13 (2 h), tumour to tissue contrast ratio $\pm$ S.D	
Tissue	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 3)$	Normal group $(n = 4)$	Blocking group $(n = 3)$	Normal group $(n = 3)$	Blocking group $(n = 4)$	Normal group $(n = 4)$	-
Blood	$1.39 \pm 0.15$	$1.12 \pm 0.08$	$0.58 \pm 0.20$	$0.91 \pm 0.11$	$0.70 \pm 0.09$	$0.64 \pm 0.10$	$0.85 \pm 0.18$	$1.72 \pm 0.38$	$0.54 \pm 0.14$	-
Tumour	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	-
Heart	$2.63 \pm 0.43$	$2.23 \pm 0.18$	$1.48 \pm 0.47$	$1.27 \pm 0.20$	$3.17 \pm 2.28$	$1.47 \pm 0.22$	$1.51 \pm 0.43$	$2.24 \pm 0.44$	$1.50 \pm 0.21$	-
Lungs	$1.85 \pm 0.20$	$1.64 \pm 0.17$	$1.13 \pm 0.31$	$1.10 \pm 0.11$	$1.14 \pm 0.22$	$0.97 \pm 0.11$	$1.19 \pm 0.38$	$2.27 \pm 0.59$	$0.78 \pm 0.20$	-
Liver	$0.11 \pm 0.02$	$0.11 \pm 0.03$	$0.13 \pm 0.03$	$0.08 \pm 0.02$	$0.15 \pm 0.03$	$0.10 \pm 0.02$	$0.09 \pm 0.02$	$0.15 \pm 0.04$	$0.07 \pm 0.01$	-
Spleen	$0.35 \pm 0.06$	$0.43 \pm 0.15$	$0.23 \pm 0.10$	$0.29 \pm 0.01$	$0.38 \pm 0.08$	$0.13 \pm 0.01$	$0.15 \pm 0.04$	$0.63 \pm 0.20$	$0.20 \pm 0.05$	-
Stomach	$4.87 \pm 1.04$	$2.91 \pm 0.18$	$5.49 \pm 1.01$	$3.38 \pm 0.31$	$9.12 \pm 2.28$	$3.60 \pm 1.68$	$2.82 \pm 1.02$	$2.78 \pm 1.19$	$2.61 \pm 1.33$	-
Pancreas	$2.80 \pm 0.33$	$2.59 \pm 0.07$	$2.31 \pm 0.54$	$1.50 \pm 0.24$	$2.72 \pm 0.60$	$2.21 \pm 0.54$	$1.32 \pm 0.46$	$2.79 \pm 0.71$	$4.15 \pm 4.17$	-
Kidney	$0.02 \pm 0.01$	$0.01 \pm 0.00$	$0.94 \pm 0.19$	$0.02 \pm 0.00$	$0.13 \pm 0.04$	$0.04\pm0.00$	$0.01 \pm 0.00$	$0.01 \pm 0.00$	$0.08 \pm 0.02$	-
Sm. Int.	$2.43 \pm 0.53$	$1.21 \pm 0.79$	$1.52 \pm 0.64$	$1.67 \pm 0.18$	$2.38 \pm 0.84$	$1.26 \pm 0.56$	$1.41 \pm 0.36$	$1.31 \pm 0.79$	$1.28 \pm 0.79$	-
Large Int.	$3.98 \pm 0.67$	$2.58 \pm 0.49$	$3.84 \pm 0.63$	$2.06 \pm 0.24$	$3.30 \pm 0.65$	$2.49 \pm 0.25$	$2.80 \pm 1.40$	$3.22 \pm 1.54$	$2.14 \pm 0.89$	-
Fat	$4.46 \pm 2.14$	$2.81 \pm 1.20$	$5.33 \pm 1.22$	$2.13 \pm 1.24$	$5.39 \pm 0.91$	$2.82 \pm 0.77$	$5.06 \pm 4.05$	$3.70 \pm 1.62$	$3.68 \pm 1.03$	-
Muscle	$6.57 \pm 1.49$	5.15 ± 1.38	$4.64 \pm 0.84$	$2.97 \pm 0.99$	$5.84 \pm 1.76$	$4.29 \pm 1.53$	$3.87 \pm 1.78$	$4.98 \pm 1.12$	4.39 ± 1.44	-
Bone	$1.90 \pm 0.28$	$1.78 \pm 0.07$	$1.58 \pm 0.53$	$1.10 \pm 0.18$	$1.64 \pm 0.53$	$1.20 \pm 0.50$	$0.76 \pm 0.17$	$1.77\pm0.90$	$1.63 \pm 1.22$	-
Skin	$1.39 \pm 0.39$	$1.14 \pm 0.11$	$1.77 \pm 0.41$	$0.97 \pm 0.15$	$2.31 \pm 0.46$	$1.61 \pm 0.19$	$0.82 \pm 0.37$	$0.86 \pm 0.29$	$1.70 \pm 0.38$	-

**Table S16.** Statistical analysis of tumour uptake. Comparison of [ $^{68}$ Ga]GaNODAGA-DARPin **3**, **6**, **9**, **12**, and **13** *vs.* blocking studies and versus data obtained for the non-binding control DARPin **13**. *P*-values of unpaired, two-tailed Student's *t*-test, light blue P < 0.05, medium blue P < 0.01, dark blue P < 0.001.

	[ <sup>68</sup> Ga]G	[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>3</b> <i>P</i> -values		[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>6</b> P-values			[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>9</b> P-values			[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>12</b> P-values		
Tissue	Normal vs. Blocking Group	Normal group vs. DARPin 13	Blocking group vs. DARPin 13	Normal vs. Blocking Group	Normal group vs. DARPin 13	Blocking group vs. DARPin 13	Normal vs. Blocking Group	Normal group vs. DARPin 13	Blocking group vs. DARPin 13	Normal vs. Blocking Group	Normal group vs. DARPin 13	Blocking group vs. DARPin 13
Blood	0.09702	0.00078	0.00025	0.07526	0.34196	0.00224	0.03115	0.12729	0.05563	0.03455	0.00038	0.00055
Tumour	0.01709	0.17359	0.00834	0.01770	0.04219	0.26351	0.03494	0.03512	0.94992	0.59791	0.05810	0.00836
Heart	0.06954	0.00002	0.00019	0.20224	0.20172	0.97336	0.96178	0.90396	0.78926	0.06482	0.02125	0.00003
Lungs	0.04507	0.01690	0.01076	0.09533	0.78346	0.05382	0.03033	0.68382	0.19072	0.10138	0.01734	0.00957
Liver	0.02755	0.00009	0.00007	0.56930	0.15188	0.03244	0.32036	0.00339	0.01682	0.04143	0.00141	0.00006
Spleen	0.00610	0.00333	0.00100	0.14316	0.53572	0.00717	0.00399	0.24686	0.00627	0.07492	0.57054	0.00088
Stomach	0.73876	0.06555	0.06761	0.64946	0.15064	0.13569	0.31674	0.07049	0.43532	0.95300	0.10464	0.11210
Pancreas	0.01054	0.40687	0.16807	0.50175	0.16451	0.24656	0.04791	0.18370	0.65532	0.03533	0.57747	0.17065
Kidney	0.54947	0.05044	0.01440	0.01040	0.00093	0.01478	0.00991	0.89153	0.00809	0.59568	0.01359	0.00220
Sm. Int.	0.40708	0.11419	0.25938	0.17990	0.79220	0.18442	0.58626	0.39633	0.83095	0.57899	0.15430	0.27558
Large Int.	0.91093	0.04792	0.04207	0.20419	0.21244	0.31691	0.24881	0.83937	0.28324	0.36956	0.05622	0.03644
Fat	0.94203	0.59263	0.57902	0.30036	0.95868	0.29866	0.26168	0.35688	0.04516	0.91593	0.09864	0.03169
Muscle	0.42882	0.02525	0.01294	0.89979	0.11505	0.31192	0.44409	0.25741	0.89627	0.07747	0.08716	0.01294
Bone	0.01276	0.10899	0.04993	0.64339	0.64747	0.99076	0.70600	0.22263	0.66190	0.03928	0.94824	0.08940
Skin	0.05506	0.92484	0.00084	0.51372	0.12774	0.00012	0.15926	0.09973	0.55495	0.42602	0.05489	0.50551

**Table S17.** Statistical analysis of tissue uptake showing the comparison of [ $^{68}$ Ga]GaNODAGA-DARPin **3**, **6**, **9** for the normal groups vs. the other PSA-binding [ $^{68}$ Ga]GaNODAGA-DARPin. P-values of unpaired, two-tailed Student's t-test, light blue P <0.05, medium blue P <0.01, dark blue P <0.001.

	[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>3</b> <i>P</i> -values			[ <sup>68</sup> Ga]	GaNODAGA-DA P-values	RPin 6	[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>9</b> <i>P</i> -values		
Tissue	vs. DARPin 6	vs. DARPin 9	vs. DARPin 12	vs. DARPin 9	vs. DARPin 12	-	vs. DARPin 12	-	-
Blood	0.01591	0.00268	0.24203	0.79963	0.05598	-	0.00168	-	-
Tumour	0.00739	0.01587	0.22942	0.41247	0.00778	-	0.00659	-	-
Heart	0.01952	0.13031	0.16199	0.31077	0.06445	-	0.22822	-	-
Lungs	0.00666	0.00261	0.42988	0.46117	0.03160	-	0.00254	-	-
Liver	0.02570	0.00179	0.76375	0.92265	0.08995	-	0.01767	-	-
Spleen	0.04748	0.02386	0.18095	0.30989	0.37125	-	0.92771	-	-
Stomach	0.00227	0.66281	0.17180	0.13585	0.34894	-	0.43481	-	-
Pancreas	0.02825	0.000005	0.08375	0.77092	0.19686	-	0.18667	-	-
Kidney	0.02720	0.05113	0.78697	0.00050	0.00970	-	0.01424	-	-
Sm. Int.	0.06255	0.03838	0.38266	0.44223	0.14004	-	0.10325	-	-
Large Int.	0.00038	0.02622	0.52107	0.16799	0.07363	-	0.02947	-	-
Fat	0.62365	0.34850	0.43407	0.44159	0.08939	-	0.04145	-	-
Muscle	0.00036	0.00105	0.15258	0.39092	0.00753	-	0.00769	-	-
Bone	0.01973	0.03087	0.02674	0.43982	0.61878	-	0.15407	-	-
Skin	0.19665	0.20550	0.05877	0.35818	0.83524	-	0.29381	-	-

**Table S18.** Statistical analysis of tissue uptake showing the comparison of [ $^{68}$ Ga]GaNODAGA-DARPin **3**, **6**, **9** from the blocking groups vs. the other [ $^{68}$ Ga]GaNODAGA-DARPin (also blocking groups). P-values of unpaired, two-tailed Student's t-test, light blue P < 0.05, medium blue P < 0.01, dark blue P < 0.001.

	[ <sup>68</sup> Ga](	GaNODAGA-DA P- values	RPin 3	[ <sup>68</sup> Ga]	GaNODAGA-DA  P-values	RPin 6	[ <sup>68</sup> Ga]GaNODAGA-DARPin <b>9</b> <i>P</i> -values		
Tissue	vs. DARPin 6 - Block	vs. DARPin 9 - Block	vs. DARPin 12 - Block	vs. DARPin 9 - Block	vs. DARPin12 - Block	-	vs. DARPin 12 - Block	-	-
Blood	0.00360	0.00161	0.08937	0.00068	0.00003	-	0.000001	-	-
Tumour	0.03824	0.02115	0.66540	0.29570	0.04241	-	0.01674	-	-
Heart	0.00196	0.00245	0.72667	0.77004	0.0000003	-	0.00844	-	-
Lungs	0.00083	0.00021	0.18900	0.00958	0.00047	-	0.00025	-	-
Liver	0.02433	0.01653	0.13697	0.88617	0.02286	-	0.01679	-	-
Spleen	0.00464	0.00013	0.23376	0.00025	0.00242	-	0.00015	-	-
Stomach	0.01065	0.28610	0.39659	0.54458	0.73304	-	0.47184	-	-
Pancreas	0.00054	0.00809	0.94634	0.09375	0.00004	-	0.00775	-	-
Kidney	0.16186	0.06331	0.21852	0.03119	0.63059	-	0.00912	-	-
Sm. Int.	0.78852	0.39668	0.95085	0.31560	0.72504	-	0.41825	-	-
Large Int.	0.00784	0.00943	0.71762	0.66831	0.00201	-	0.00249	-	-
Fat	0.23106	0.18448	0.55098	0.48872	0.18223	-	0.00354	-	-
Muscle	0.02737	0.03935	0.79346	0.39475	0.03775	-	0.05193	-	-
Bone	0.06860	0.14282	0.48526	0.66259	0.08927	-	0.17179	-	-
Skin	0.00001	0.00924	0.11534	0.00459	0.16907	-	0.60288	-	-

## References

- (1) Zanzonico, P. Routine Quality Control of Clinical Nuclear Medicine Instrumentation: A Brief Review. *J. Nucl. Med.* **2008**, 49 (7), 1114–1131. https://doi.org/10.2967/jnumed.107.050203.
- (2) Institute for Laboratory Animal Research. *Guide for the Care and Use of Laboratory Animals: 8th Ed.*; 2011. https://doi.org/10.2307/1525495.
- (3) Fridman, R.; Benton, G.; Aranoutova, I.; Kleinman, H. K.; Bonfil, R. D. Increased Initiation and Growth of Tumor Cell Lines, Cancer Stem Cells and Biopsy Material in Mice Using Basement Membrane Matrix Protein (Cultrex or Matrigel) Co-Injection. *Nat. Protoc.* 2012, 7 (6), 1138–1144. https://doi.org/10.1038/nprot.2012.053.
- (4) Dreier, B.; Plückthun, A. Rapid Selection of High-Affinity Binders Using Ribosome Display. In *Methods in Molecular Biology*; Douthwaite, J. A., Jackson, R. H., Eds.; Springer New York: New York, NY, 2012; Vol. 805, pp 261–286. https://doi.org/10.1007/978-1-61779-379-0 15.
- (5) Binz, H. K.; Stumpp, M. T.; Forrer, P.; Amstutz, P.; Plückthun, A. Designing Repeat Proteins: Well-Expressed, Soluble and Stable Proteins from Combinatorial Libraries of Consensus Ankyrin Repeat Proteins. *J. Mol. Biol.* **2003**, 332 (2), 489–503. https://doi.org/10.1016/S0022-2836(03)00896-9.
- (6) Kramer, M. A.; Wetzel, S. K.; Plückthun, A.; Mittl, P. R. E.; Grütter, M. G. Structural Determinants for Improved Stability of Designed Ankyrin Repeat Proteins with a Redesigned C-Capping Module. *J. Mol. Biol.* **2010**, *404* (3), 381–391. https://doi.org/10.1016/j.jmb.2010.09.023.
- (7) Plückthun, A. Designed Ankyrin Repeat Proteins (DARPins): Binding Proteins for Research, Diagnostics, and Therapy. *Ann. Rev. Pharmacol. Toxicol.*, 2015, pp 489–511. https://doi.org/10.1146/annurev-pharmtox-010611-134654.
- (8) Brauchle, M.; Hansen, S.; Caussinus, E.; Lenard, A.; Ochoa-Espinosa, A.; Scholz, O.; Sprecher, S. G.; Plückthun, A.; Affolter, M. Protein Interference Applications in Cellular and Developmental Biology Using DARPins That Recognize GFP and MCherry. *Biol. Open* 2014, 3 (12), 1252–1261. https://doi.org/10.1242/bio.201410041.
- (9) Schilling, J.; Schöppe, J.; Plückthun, A. From DARPins to LoopDARPins: Novel LoopDARPin Design Allows the Selection of Low Picomolar Binders in a Single Round of Ribosome Display. *J. Mol. Biol.* **2014**, *426* (3), 691–721. https://doi.org/10.1016/j.jmb.2013.10.026.
- Zahnd, C.; Wyler, E.; Schwenk, J. M.; Steiner, D.; Lawrence, M. C.; McKern, N. M.; Pecorari, F.; Ward, C. W.; Joos,
   T. O.; Plückthun, A. A Designed Ankyrin Repeat Protein Evolved to Picomolar Affinity to Her2. *J. Mol. Biol.* 2007,
   369 (4), 1015–1028. https://doi.org/10.1016/j.jmb.2007.03.028.
- (11) Andres, F.; Schwill, M.; Boersma, Y. L.; Plückthun, A. High-Throughput Generation of Bispecific Binding Proteins by Sortase a-Mediated Coupling for Direct Functional Screening in Cell Culture. *Mol. Cancer Ther.* **2020**, *19* (4), 1080–1088. https://doi.org/10.1158/1535-7163.MCT-19-0633.
- (12) Forrer, P.; Jaussi, R. High-Level Expression of Soluble Heterologous Proteins in the Cytoplasm of Escherichia Coli by Fusion to the Bacteriophage Lambda Head Protein D. *Gene* **1998**, *224* (1–2), 45–52. https://doi.org/10.1016/S0378-1119(98)00538-1.
- (13) Lindmo, T.; Boven, E.; Cuttitta, F.; Fedorko, J.; Bunn, P. A. Determination of the Immunoreactive Function of Radiolabeled Monoclonal Antibodies by Linear Extrapolation to Binding at Infinite Antigen Excess. *J. Immunol. Methods* **1984**, 72 (1), 77–89. https://doi.org/10.1016/0022-1759(84)90435-6.