

Supporting information for

UHPLC-MS/MS combined with microdialysis for simultaneous determination of nicotine and neurotransmitter metabolites in the rat hippocampal brain region: application to pharmacokinetic and pharmacodynamic study

Mingyu Zhu^{a,b,c}, Lili Cui^{b,c}, Guanglin Liu^{b,c}, Pengpeng Yu^{b,c}, Qingyuan Hu^{b,c}, Huan Chen^{*b,c}, and Hongwei Hou^{*b,c}

^aInstitute of Biomedical Engineering, College of Life Sciences, Qingdao University, Qingdao, 266071, China

^bChina National Tobacco Quality Supervision and Test Center, key Laboratory of Tobacco Biological Effects, Zhengzhou, 450001, China

^c Beijing Life Science Academy, Beijing, 100101, China

*^{b,c} Huan Chen, E-mail: hunny_ch@163.com; Tel: +86-0371-67672597

*^{b,c} Hongwei Hou, E-mail: qsfctc@163.com; Tel: +86-0371-67672727

Table.S1 Comparison of the LOD and LOQ of the proposed method with other methods.

Analyte	LOD (ng/mL)	LOQ (ng/mL)	LOD ^[1] (ng/mL)	LOQ ^[1] (ng/mL)	LOD ^[2] (ng/mL)	LOQ ^[2] (ng/mL)	LOD ^[3] (ng/mL)	LOQ ^[3] (ng/mL)	LOD ^[4] (ng/mL)	LOQ ^[4] (ng/mL)
DA	0.05	0.2	-	-	0.05	0.25	0.20	0.50	0.05	0.25
DOPAC	0.02	0.05	-	-	--	--	2.50	5.00	-	-
L-DOPA	0.01	0.02	-	-	0.26	0.50	-	-	0.50	0.75
3-MT	0.01	0.04	-	-	0.25	0.75	-	-	0.10	0.25
HVA	0.01	0.05	-	-	17.1	50.0	15.0	450	-	-
5-HT	0.004	0.01	-	-	0.07	0.25	0.01	0.03	0.01	0.05
5-HIAA	0.004	0.01	-	-	0.16	0.50	1.00	3.00	0.10	0.25
NE	0.01	0.03	-	-	0.96	2.00	1.00	3.00	0.75	2.00
E	0.02	0.10	-	-	-	-	0.02	0.05	-	-
GABA	0.01	0.03	-	-	0.008	0.03	30.0	100	0.75	1.00
GA	0.10	0.20	-	-	0.02	1.25	-	-	-	-
ACh	0.004	0.01	-	-	0.03	0.10	-	-	0.03	0.05
Nic	0.01	0.03	0.01	0.04	-	-	-	-	-	-
Cot	0.01	0.02	0.02	0.05	-	-	-	-	-	-
NorNic	0.02	0.10	0.03	0.09	-	-	-	-	-	-
NorCot	0.01	0.03	0.02	0.07	-	-	-	-	-	-
NNO	0.01	0.02	0.02	0.05	-	-	-	-	-	-
CNO	0.02	0.10	0.03	0.10	-	-	-	-	-	-
3OH-Cot	0.03	0.10	0.04	0.14	-	-	-	-	-	-
NicGlu	0.01	0.04	0.08	0.27	-	-	-	-	-	-




Fig.S1 Chemical formula of the reaction of DA with benzoyl chloride




Fig.S2 (a)Chromatograms of 3Bz-DA and 2Bz-DA under the condition of 2:1:1;

(b)Chromatograms of 3Bz-DA and 2Bz-DA under the condition of 1:2:2.

References:

- [1] XU Y, ZHANG Q, LI P, et al. Nicotine Pharmacokinetics in Rat Brain and Blood by Simultaneous Microdialysis, Stable-Isotope Labeling, and UHPLC-HRMS: Determination of Nicotine Metabolites [J]. Analytical Chemistry, 2019, 91(4): 2916-22.
- [2] BECKER S, SCHULZ A, KREYER S, et al. Sensitive and simultaneous quantification of 16 neurotransmitters and metabolites in murine microdialysate by fast liquid chromatography-tandem mass spectrometry [J]. Talanta, 2023, 253: 123965.
- [3] WANG Y C, ZHENG R J, WU P P, et al. Determination of Multiple Neurotransmitters through LC-MS/MS to Confirm the Therapeutic Effects of Althaea rosea Flower on TTX-Intoxicated Rats [J]. Molecules, 2023, 28(10).
- [4] HELMSCHRODT C, BECKER S, PERL S, et al. Development of a fast liquid chromatography-tandem mass spectrometry method for simultaneous quantification of neurotransmitters in murine microdialysate [J]. Analytical and Bioanalytical Chemistry, 2020, 412(28): 7777-87.