

Supporting Information

Convenient in situ self-assembled formation of dual-functional Ag/MXene nanozymes for efficient chemiluminescence sensing

Xiangjuan Zheng^{b,d}, Zhiying Shi^b, Chun Han^c, Hongyi Mu^b, Shiyun Cheng^b, Xiluan Yan^{*,a,b}

^a College of Pharmacy, Nanchang University, Nanchang 330031, China. ^b School of Chemistry and

^b Chemical Engineering, Nanchang University, Nanchang 330031, China.

^c The Collaboration Unit for Field Epidemiology of State Key Laboratory for Infectious Disease

Prevention and Control, Nanchang Centre for Disease Control and Prevention, Nanchang, P. R.

China 330038

^d Chongqing Research Institute of Nanchang University, Chongqing 402660, China

corresponding author: Xiluan Yan

Yanxiluan@ncu.edu.cn

EXPERIMENTAL SECTION

Reagents


Ti₃AlC₂ was purchased from Lai Zhou Kaiene Ceramic Material Co., LTD. HF, H₂O₂, NaOH were obtained from Xi long Chemical Co., LTD. Citric acid, glucose and urea were purchased from Shanghai Sinopharm Group Chemical Reagent Co., LTD. Luminol was purchased from Beijing Balingway Technology Co., LTD. Reduced Glutathione, Uricase (20 U/mg), Calcium acetate monohydrate (C₄H₈CaO₅), uric acid, cysteine, lysine, tryptophan, threonine, glutamic acid, arginine, alanine, valine, glycine, isoleucine, aspartic acid, leucine, methionine, hydroxylamine hydrochloride, thiourea, ascorbic acid were purchased from Shanghai Maclin Biochemical Technology Co., LTD. NaCl, KCl, MgCl₂, AgNO₃ were purchased from Shanghai Aladdin Biochemical Technology Co., LTD. All reagents are analytically pure and used directly. The ultrapure water used in all experiments was purified by Millipore Synergy System (18.24 MΩ cm⁻¹, Milli-Q, Millipore).

Apparatus and Characterization

All chemiluminescence (CL) signals were collected using an ultra-weak luminescence analyzer (BPCL-1-TGC, Institute of Biophysics, Chinese Academy of Sciences, China). UV-vis absorption analysis was performed on the UV-vis spectrometer (TU-1901, Puxi General Instrument Co., Ltd., China Beijing). A scanning electron microscope (SEM,

TESCAN MIRA LMS, Czech Republic) was used to characterize the formation of Ag/MXene. A transmission electron microscope (TEM, JEM-2010HR, JEOL, Japan), with a 200 kV acceleration voltage was used to measure the size. X-ray powder diffraction spectra (XRD) of the Ag/MXene powders were characterized by Cu K α radiation ($\lambda = 1.5418 \text{ \AA}$) using an x-ray diffractometer (D8-FOCUS, Bruker, Karlsruhe, Germany) with diffraction angles ranging from 5° to 90°. The elemental valence states of Ag/MXene were determined by X-ray photoelectron spectroscopy (XPS) using Thermo Scientific K-Alpha (EscaLab250, Thermo Fisher Scientific, America). The 300E electron spin resonance (EPR) spectrometer (Bruker EMXplus, Germany) was employed to record the EPR spectra at room temperature.

Supplementary Figures



Scheme S1. Preparation process of Ag/MXene composites.




Figure S1. The EDS images of Ag/MXene composites.




Figure S2. XRD image for Ag/MXene composites (A). Full-scan XPS spectrum of Ag/MXene (B). High-resolution XPS spectra of C 1s (C), O 1s (D), Ti 2p (E) and Ag 3d (F) for Ag/MXene.




Figure S3. Optimization of Ag/MXene synthesis conditions, including MXene concentration (A), AgNO_3 concentration (B), the churning time (C) and the mixing volume ratio between MXene solution and AgNO_3 solution (D). Reaction condition: luminol, 1 μM , 50 μL ; H_2O_2 , 10 mM, 50 μL ; Ag/MXene, 10 times diluted, 50 μL .




Figure S4. Ag/MXene-luminol-H₂O₂ measurement condition optimization, including substrates luminol and H₂O₂ volume ratio (A), pH of luminol (B) and luminol-H₂O₂ reaction time (C).




Figure S5. Sensing analysis of H₂O₂ concentration by Ag/MXene composites in chemiluminescence luminol-H₂O₂ system. The error bars represented the standard deviation of three replicates.




Figure S6. The time stability (A) and temperature stability (B) tests of Ag/MXene composites.




Figure S7. Raman spectra of Ag/MXene composites.




Figure S8. (A) Zeta potential and (B) particle size analysis between different concentrations of AgNO₃ and MXene.




Figure S9. UV-vis absorption spectra of luminol, luminol-H₂O₂, Ag/MXene and Ag/MXene-luminol-H₂O₂ (A). Scavenging experiment of reactive oxygen species with ascorbic acid (B).

Table S1. Application with the proposed method to detect human Serum spiked with varied GSH amounts

samples	Spiked (μM)	Found (μM)	Recovery (%)	RSD (% , n=3)
1	0.5	0.56	10.68	1.94
	5	5.59	111.80	3.08
	10	10.68	106.80	5.55
2	0.5	0.48	95.72	1.56
	5	4.94	98.70	5.68
	10	10.42	104.20	5.28

Table S2. Application with the proposed method to detect human Serum spiked with varied UA amounts

samples	Spiked (μM)	Found (μM)	Recovery (%)	RSD (% , n=3)
1	10	10.49	104.90	1.88
	15	14.63	97.55	5.87
	20	19.52	97.59	5.12
2	10	10.6	106.00	4.07
	15	16.25	108.32	4.56
	20	21.67	108.34	1.23