

Supporting Information

Multipurpose Mitochondrial NIR Probe for Imaging Ferroptosis and Mitophagy

Deeksha Rajput^a, Nachiket Pradhan^a, Shabnam Mansuri^a, Virupakshi Soppina^{c*} and Sriram Kanvah^{a*}

^aDepartment of Chemistry, Indian Institute of Technology Gandhinagar, Palaj, Gandhinagar, 382355

^bDepartment of Biological Engineering, Indian Institute of Technology Gandhinagar, Palaj, Gandhinagar 382355

sriram@iitgn.ac.in

Contents	
S1.	Table S1: List of a few previously reported probes for similar function
S2.	Absorption and emission spectra
S3.	Lifetime data: Figure S2 and Table S2
S4.	Viscosity sensing in solution
S5.	Influence of added analytes on the fluorescence properties of PP and NP
S6.	Live cell imaging
S7.	Cytotoxicity assay
S8.	Live Cell Imaging with probe NP with Nystatin treatment
S9.	Live Cell Imaging with probe NP during Ferroptosis
S10.	Synthetic scheme and synthetic procedures
S11.	Characterization data
S12.	Spectral copies

Previously reported mitochondrial probes for ferroptosis/mitophagy:

Table S1: Details of Previously reported probes

Probe	Emission Maxima	MMP Independent	Ferroptosis/ Mitophagy	Viscosity Sensing (Solution/Endogenous in cells)	Ref.
1.	668	Yes	No/Yes	Yes/ Yes	1
2.	675	Yes	Yes/No	Yes/Yes	2
3.	670	No	No/Yes	Yes/Yes	3
4.	704	No	Yes/No	Yes/Yes	4
5.	795	No	Yes/No	Yes/Yes	5
6.	602	No	Yes/No	Yes/Yes	6
7.	553	No	Yes/No	No/No	7
8.	700	Yes	Yes/Yes	Yes/Yes	<i>This Work</i>




Figure S1: Molecular structures of previously reported probes which detect ferroptosis/mitophagy process but none utilizing this architecture are known.

S2: Absorption and emission spectra:





Figure S2: A) Absorption and B) Emission spectra of PP in all solvents. Concentration: 10 μ M.

S3: Lifetime decay




Figure S3: A) & B) Lifetime decay for NP & PP in different solvents.

Table S2: Lifetime decay for the compounds in different solvents.

NP						
Components	τ_1 (ns)	A1	τ_2 (ns)	A2	Average Lifetime (τ) (ns)	χ^2
Dioxane	0.09	47.4%	0.30	52.5%	0.15	1.03
CH ₃ CN	0.10	100%	-	-	0.10	0.96
DMSO	0.26	100%	-	-	0.26	1.18
Water	0.07	68.0%	0.41	31.9%	0.07	1.04
PP						
Components	τ_1 (ns)	A1	τ_2 (ns)	A2	Average Lifetime (τ) (ns)	χ^2
Dioxane	0.07	63.75%	0.25	36.25%	0.10	1.32
CH ₃ CN	0.03	100%	-	-	0.03	1.22
DMSO	0.11	100%	-	-	0.11	1.01
Water	0.03	100%	-	-	0.03	1.16

*Most of the decays are very fast and overlap with the instrument response function (IRF). The FWHM of IRF is 70 ps, and the resolution of the experiment is 13 ps. For such data, it is difficult to resolve with our current instrument setup.

S4: Viscosity sensing in solution:




Figure S4: Emission spectra of A) PP in different viscous medium; B) Viscosity vs Log I with linear relationship for PP (Inset – Colour of the fluorophore in 0% glycerol and 100 % glycerol)

S5: Ion selectivity of the probe:




Figure S5: Fluorescence response of A) NP and B) PP to different competitive analytes (A–K: F⁻, Cl⁻, Br⁻, I⁻, NO₂⁻, NO₃⁻, S²⁻, CH₃COO⁻, Ca²⁺, Mg²⁺, Fe³⁺. Concentration of probes and analytes: 10 μM.

S6: Live Cell imaging:




Figure S6: Colocalization of PP in COS-7 cells. CLSM image of A) PP, B) MitoView 405, C) Merged image, D) Bright field image and E) Intensity profile.

S7: Cytotoxicity assay:




Figure S7: Cytotoxicity assay for PP & NP

S8: Live Cell Imaging with probe NP with Nystatin treatment:




Figure S8: CLSM Images of NP with before (first column) and after Nystatin treatment (second column). (Images with multiple cells)

S9: Live Cell Imaging with probe NP during Ferroptosis:




Figure S9: CLSM Images of NP with no treatment (first column), erastin treatment (second column), erastin + ferrostatin-1 treatment (third column), and ferrostatin-1 treatment (fourth column) after 0 min (first row) and 30 min (second row) incubation with NP. (Images with multiple cells)

S10: Synthetic procedures:

General procedures:

Styrene 3 (Heck Coupling):

In a sealed reactor, 1 molar equivalent of N, N-dimethylamino bromonaphthalene derivative **2** was dissolved in dry DMF, followed by the addition of 3 molar equivalents of triethylamine, a pinch of tris (o-tolyl) phosphine ligand and Pd(OAc)₂. The contents were stirred for 10 min at room temperature. After the addition of 3 molar equivalents of vinyl pyridine, the reaction was tightly sealed and allowed to heat at 100° C. After completion of the reaction, as monitored by TLC, the reaction mixture is triturated with n-pentane and filtered off to get yellowish orange coloured powder. (64% yield)

General procedure for N-Alkylation: The styryl pyridyl derivatives (**3** or **5**) were dissolved in dry acetonitrile and allowed to stir at room temperature for 15 mins, followed by the addition of 0.51 molar equivalents of 1,12-dibromo dodecane. The reaction mixture was refluxed for 24

h. The solvent was evaporated using a rotatory evaporator and the obtained precipitate was washed with diethyl ether and pentane and filtered to get pure dark red colored solid powder. (80% yield for PP and 70% yield for NP)

S11: Characterization data:

PP: 1,1'-dodecane-1,12-diylbis(4-((E)-4-(dimethylamino)styryl)pyridin-1-ium bromide

^1H NMR (500 MHz, *DMSO-d6*) δ 8.79 (d, $J = 6.7$ Hz, 4H), 8.08 (d, $J = 6.7$ Hz, 4H), 7.94 (d, $J = 16.0$ Hz, 2H), 7.60 (d, $J = 8.8$ Hz, 4H), 7.19 (d, $J = 16.1$ Hz, 2H), 6.79 (d, $J = 8.9$ Hz, 4H), 4.41 (t, $J = 7.2$ Hz, 4H), 3.03 (s, 12H), 1.87 (m, 4H), 1.25 (d, $J = 17.4$ Hz, 16H). ^{13}C NMR (126 MHz, *DMSO*) δ 153.12, 151.31, 142.87, 141.57, 129.58, 121.86, 121.77, 116.49, 111.34, 58.49, 39.51, 39.42, 39.34, 39.25, 39.17, 39.08, 39.01, 38.91, 38.75, 38.58, 38.41, 29.86, 28.27, 28.21, 27.81, 24.84. HRMS: Observed Mass: 308.2207; Calculated Mass: 308.2247

NP: 1,1'-dodecane-1,12-diyl bis(4-((E)-2-(6-(dimethylamino) naphthalen-2-yl) vinyl) pyridin-1-ium bromide

^1H NMR (500 MHz, *DMSO-d6*) δ 8.91 (d, $J = 6.7$ Hz, 4H), 8.21 (d, $J = 6.8$ Hz, 4H), 8.13 (d, $J = 16.2$ Hz, 2H), 8.01 (s, 2H), 7.80 (dd, $J = 8.8, 5.9$ Hz, 4H), 7.73 (d, $J = 8.8$ Hz, 2H), 7.50 (d, $J = 16.2$ Hz, 2H), 7.29 – 7.21 (m, 2H), 6.97 (d, $J = 2.4$ Hz, 2H), 4.47 (t, $J = 7.2$ Hz, 4H), 3.07 (s, 12H), 1.89 (m, 4H), 1.26 (d, $J = 19.8$ Hz, 16H). ^{13}C NMR (126 MHz, *DMSO*) δ 153.72, 150.06, 144.40, 142.28, 136.43, 130.77, 130.15, 128.86, 127.27, 125.93, 124.40, 123.69, 121.19, 116.99, 105.88, 59.93, 40.59, 40.50, 40.42, 40.33, 40.25, 40.16, 40.09, 40.00, 39.83, 39.66, 39.50, 35.71, 32.69, 30.98, 29.37, 29.30, 28.91, 28.56, 27.97, 27.11, 25.93, 20.31. HRMS: Observed mass: 358.2392; Calculated mass: 358.2404

S12: Spectral copies:




Figure S12a: ¹H NMR of PP




Figure S12b: ¹³C NMR of PP




Figure S12c: ^1H NMR of Styrene 3




Figure S12d: ^1H NMR of NP





Figure S12e: ¹³C NMR of NP

SKG NP NME2== IIT, GN
IITGN_Pos_15_MARCH_23_24_355 (8.372) AM2 (Ar,30000.0,566.28,0.00,LS 7); ABS

17-Mar-2023 14:10:26
1: TOF MS ES+
9.66e5



SKG DR NAPTH ++
IITGN_Pos_29_May_23_76_393 (9.252) AM2 (Ar,30000.0,556.28,0.00,LS 10); ABS

IIT, GN

03-Jun-2023 12:50:48
1: TOF MS ES+
1.81e6




Figure S12g: HRMS of NP

References:

1. J. Hong, X. Guan, Y. Chen, X. Tan, S. Zhang and G. Feng, *Anal. Chem.*, 2023, **95**, 5687-5694.
2. X. Wu, R. Zhang, Y. Li, Y. Gai, T. Feng, J. Kou, F. Kong, L. Li and B. Tang, *Anal. Chem.*, 2023, **95**, 7611-7619.
3. P. Lei, M. Li, C. Dong and S. Shuang, *ACS Biomaterials Science & Engineering*, 2023, **9**, 3581-3589.
4. L. Fan, Q. Yang, Q. Zan, K. Zhao, W. Lu, X. Wang, Y. Wang, S. Shuang and C. Dong, *Anal. Chem.*, 2023, **95**, 5780-5787.
5. J. Yin, Q. Xu, X. Mo, L. Dai, M. Ren, S. Wang and F. Kong, *Dyes and Pigments*, 2022, **200**, 110184.
6. F. Zheng, J. Ding, S. Huang, A. Bi, S. Liu, K. Zhang, F. Chen and W. Zeng, *Dyes and Pigments*, 2023, **217**, 111424.
7. X. Xie, Y. Liu, G. Liu, Y. Zhao, J. Bian, Y. Li, J. Zhang, X. Wang and B. Tang, *Anal. Chem.*, 2022, **94**, 10213-10220.