

Supplementary information for

Comparison of neutron activation analysis, X-ray fluorescence spectrometry and inductively-coupled plasma mass spectrometry for the determination of rare earth element concentrations in Jordanian monazite ore

Ahmad Alsabbagh^{1*}, Kamal Harahshah², Kafa Al-Khasawneh³, Elizabeth J. New⁴

1. Department of Nuclear Engineering, Jordan University of Science and Technology, Irbid 22110, Jordan
2. Jordan Atomic Energy Commission, P.O. Box 70, Amman, 11934, Jordan
3. Jordan Research and Training Reactor, Irbid 22110, Jordan
4. School of Chemistry, The University of Sydney, NSW 2006, Australia

* ahalsabbagh@just.edu.jo

Table of Contents

Supplementary Figures	2
Figure S1: Publication data (<i>Web of Science</i> , accessed 1 st October 2024) for different analytical techniques applied to monazite: (a) “neutron activation analysis” and “monazite”, (b) “X-ray fluorescence” and “monazite” and (c) “ICP-MS” and “monazite”	2
Figure S2: Percentage recoveries for certified reference materials (a) REE-1 and (b) REE-2 for INAA, XRF, and ICP-MS.....	3
Figure S3: Comparison of value deflections (Def) for La, Ce, Pr, Nd, and Sm for INAA, XRF and ICP-MS measurements of certified reference materials REE-1 and REE-2. Error bars represent standard error of the mean.	4
Figure S4: Total rare earth element content of different monazite samples, determined by ICP-MS.	4
Figure S5: Bland-Altman plots for analysis of La, Ce, Nd, Pr, and Sm in monazite samples using ICP-MS, INAA, and XRF. Black line indicates bias, gray lines indicate limits of agreement. (a) and (d) are low concentration monazite samples (Bulk, Tail 1, Tail 2, Tail 3); (b) and (e) are medium concentration samples (Tail 4, Tail 5, Tail 7, Conc. 1); (c) and (f) are high concentration samples (Tail 6, Tail 8, Conc. 2, Conc. 3, Conc. 4).....	5
Supplementary Tables	6
Table S1: Statistical analysis of differences between Def values for the three techniques for analysis of REE-1 and REE-2.....	6
Table S2: Elemental concentrations of monazite samples determined by ICP-MS, INAA and XRF	7

Supplementary Figures




Figure S1: Publication data (*Web of Science*, accessed 1st October 2024) for different analytical techniques applied to monazite: (a) “neutron activation analysis” and “monazite”, (b) “X-ray fluorescence” and “monazite” and (c) “ICP-MS” and “monazite”




Figure S2: Percentage recoveries for certified reference materials (a) REE-1 and (b) REE-2 for INAA, XRF, and ICP-MS.




Figure S3: Comparison of value deflections (Def) for La, Ce, Pr, Nd, and Sm for INAA, XRF and ICP-MS measurements of certified reference materials REE-1 and REE-2. Error bars represent standard error of the mean.




Figure S4: Total rare earth element content of different monazite samples, determined by ICP-MS.




Figure S5: Bland-Altman plots for analysis of La, Ce, Nd, Pr, and Sm in monazite samples using ICP-MS, INAA, and XRF. Black line indicates bias, gray lines indicate limits of agreement. (a) and (d) are low concentration monazite samples (Bulk, Tail 1, Tail 2, Tail 3); (b) and (e) are medium concentration samples (Tail 4, Tail 5, Tail 7, Conc. 1); (c) and (f) are high concentration samples (Tail 6, Tail 8, Conc. 2, Conc. 3, Conc. 4).

Supplementary Tables

Table S1: Statistical analysis of differences between Def values for the three techniques for analysis of REE-1 and REE-2. Student's t-test (two-tailed, heteroscedastic).

Sample	INAA vs XRF	<i>p</i> -values	
		INAA vs ICP-MS	XRF vs ICP-MS
REE-1	0.047*	0.43	0.017*
REE-2	0.78	0.86	0.86

* significance determined as $p < 0.05$

Table S2: Elemental concentrations of monazite samples determined by ICP-MS, INAA and XRF

Sample	Ce			La			Nd		
	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]
Bulk	3631.12	3595.67	3767.5	1709.09	1765.62	1881.82	1440.26	1277.59	1455.76
Tail 1	1391.04	1419.86	1400.2	650.04	716.55	679.57	569.75	653.9	602.71
Tail 2	1702.77	2006.58	1895.96	949	1014.2	872.27	746.55	783.89	722.74
Conc.1	24200.2	22809.1	21028.1	11712.2	11166.4	11695.0	9619.67	10432.5	9264.4
Tail 3	3829.27	3970.14	4032.07	1789.68	1986.91	1853.68	1624.47	1894.41	1477.19
Tail 4	9843.81	9880.14	9746.83	4767.13	5105.66	4842.26	4029.88	4305.33	4039.78
Conc.2	59045.4	53281.0	48590.8	28484.8	27280.0	29899.3	23617.9	21627.6	22281.3
Tail 5	11876.7	11283.2	11890.2	5973.27	6004.23	6066.68	4951.74	4631.82	5033.44
Tail 6	37324.6	34988.5	30345.1	18260.2	17059.0	17971.5	14837.4	20301.4	13786.8
Conc.3	69652.5	60157.1	59185.9	34991.3	30965.9	37083.0	27235.9	24603.1	27240.2
Tail 7	22637.3	23468.9	22221.5	10910.6	7397.6	12180.2	9714.67	10912.2	9612.48
Tail 8	46045.2	44778.6	39452.1	21832.0	22092.3	24378.4	17854.6	23150.3	18340.1
Conc.4	74036.3	69700.0	64710.1	35255.1	35245.0	40373.4	28773.89	29427.16	29558.47

Sample	Pr			Sm			Y		
	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]
Bulk	398.22	465.05	361.7	248.73	243.97	232.84	365.56	-	507.13
Tail 1	149.41	143.48	101.8	91.57	89.21	BDL	161.45	-	214.19
Tail 2	229.95	254.44	178.78	114.81	132.97	BDL	225.38	-	326.8
Conc.1	2629.62	2315.87	2636.15	1619.77	1479.7	1771.3	2129.85	-	2321.44
Tail 3	432.96	362.23	320.31	267.46	275.73	257.85	640.66	-	863.85
Tail 4	1100.01	1173.4	1099.98	695.23	668.15	726.97	1320.63	-	1707.22
Conc.2	6564.72	5547.66	6499.75	4086.83	3464.59	4401.52	4480.79	-	3996.38
Tail 5	1375.18	1306.19	1508.03	841.72	807.75	998.62	1442.09	-	1723.76
Tail 6	4105.39	3404.45	4081.28	2585.87	2085.45	2796.65	3360.7	-	3167.97
Conc.3	7718.68	7128.71	8069.03	4641.49	4040.4	5401	4893.88	-	4490.12
Tail 7	2734.76	1953.63	2715.61	1698.38	843.36	1864.44	2175	-	2336.4
Tail 8	5009.76	4987.82	5447.77	3081.55	2971.17	3652.12	3844.82	-	3708.95
Conc.4	8083.75	7613.73	8670.75	4914.35	4658.13	6102.97	4985.58	-	4703.52

Sample	Dy			Er			Eu		
	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]
Bulk	94.85	88.5	-	94.85	88.5	-	13.19	12.32	-
Tail 1	37.35	39.55	-	37.35	39.55	-	-	6.75	-
Tail 2	50.15	49.78	-	50.15	49.78	-	-	7.77	-
Conc.1	585.5	684.65	-	585.5	684.65	-	75.43	70.69	-
Tail 3	130.88	144.07	-	130.88	144.07	-	16.6	17.93	-
Tail 4	289.22	286.9	-	289.22	286.9	-	36.46	36.51	-
Conc.2	1359.97	1303.85	1025.53	1359.97	1303.85	1025.53	175.58	159.37	-
Tail 5	333.94	342.72	-	333.94	342.72	-	43.61	43.82	-
Tail 6	902.67	817.72	709.24	902.67	817.72	709.24	113.81	109.38	-
Conc.3	1586.91	1405.38	1241.61	1586.91	1405.38	1241.61	198.76	173.98	-
Tail 7	557	554.47	-	557	554.47	BDL	74	73.46	-
Tail 8	1045.95	1145.14	879.15	1045.95	1145.14	879.15	136.4	137.13	-
Conc.4	1650.66	1455.77	1315.67	1650.66	1455.77	1315.67	218.76	200.42	-

Sample	Gd			Yb			Ho		
	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]
Bulk	190.32	-	-	37.69	36.99	-	14.45	14.58	-
Tail 1	71.09	-	-	18.65	29.07	-	-	6.47	-
Tail 2	93.66	-	-	31.16	36.7	-	-	9.41	-
Conc.1	1243.72	-	1466.25	196.98	243.98	-	79	59.36	-
Tail 3	219.56	-	-	97.49	107.46	-	24.18	23.59	-
Tail 4	538.26	-	600.38	168.35	204.71	-	48.82	36.42	-
Conc.2	3018.26	-	3693.39	234.36	175.52	-	132.52	102.69	-
Tail 5	642.19	-	754.82	160.7	198.71	-	54.96	35.95	-
Tail 6	1919.14	-	2313.9	310.76	345.82	318.78	127.66	92.23	-
Conc.3	3565.3	-	4537.57	323.06	405.99	-	192.17	174.01	-
Tail 7	1390.55	-	1496.62	190.66	156.96	-	84.68	33.92	-
Tail 8	2206.25	-	3173.69	332.01	403.1	374.11	151.49	107.22	-
Conc.4	3667.96	-	4908.9	313.1	412.18	320.54	217.23	151.71	-

Sample	Lu			Tb			Sc		
	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]	ICP-MS [µg/g]	INAA [µg/g]	XRF [µg/g]
Bulk	-	6.18	-	22.35	18.64	-	-	16.54	-
Tail 1	-	4.08	-	8.69	8.68	-	-	17.91	-
Tail 2	-	6.07	-	11.04	12.29	-	-	12.94	-
Conc.1	36.9	32.78	-	141.4	118.76	-	-	54.59	-
Tail 3	17.07	18.52	-	27.36	25.4	-	-	38.69	-
Tail 4	32.14	35.54	-	65.36	62.38	-	-	53.04	-
Conc.2	41.66	38.7	-	256.35	282.23	-	-	63.7	-
Tail 5	31.06	31.99	-	77.55	70.56	-	-	56.03	-
Tail 6	58.72	74.53	-	219.86	184.59	-	-	78.17	-
Conc.3	54.37	49.76	-	385.41	316.23	-	-	56.96	-
Tail 7	34.62	36.39	-	142.73	132.42	-	-	56.22	-
Tail 8	70.34	74.54	-	259.52	229.87	-	-	80.46	-
Conc.4	54.05	47.97	-	432.95	356.13	-	-	56.76	-