

Supporting information

Table S1. Oligonucleotide sequences used in this work.

Oligonucleotide	Sequence (5'→3')
SL-probe	CGC GTT ACC CGA TGT TTG GGT TAG TCC ATG ACG TGC ATC ATT CAT TTG TTT CAC CTC AGC TGA TGC ACG TCA TGG GAC GTG
PER template probe	GGG TTA CCC GGT TTT CCG GGT AAC CCT AAC CCA AAC ATC GGG TAA CGC G
Target miRNA	CUC AGU GAU UUU AGA GAG AGG AU
“4*”	GGG TTA
miRNA-149	UAA CAC UGU CUG GUA AAG AUG G




Figure S1. Selectivity of the method to homological miRNAs.