

Dynamic biomimicry in skin-on-a-chip: Multi-scale construction to translational dermatology, drug screening and cosmetic evaluation

Linwei Sang,[†] Ajing Liu,[†] Haochen Wang, Wanyi Li, Qin Tu and Jinyi Wang**

College of Chemistry and Pharmacy, Northwest A&F University, Yangling, Shaanxi 712100, P. R. China

E-mail: tuqin@nwsuaf.edu.cn, jywang@nwsuaf.edu.cn

Table S1. Summary of representative skin-on-a-chip studies across different dynamic levels

Model type /dynamic level	Cells	Matrix / scaffold	Materials	Flow rates	Culture duration	Throughput	Barrier Metrics	Primary application	Limitations / future need	Ref.
Dynamic barrier (Epidermis)	NHKs	PET membrane	PMMA	Medium flow 1 $\mu\text{L}/\text{min}$; air flow 1 $\mu\text{L}/\text{min}$; shear stress $3 \times 10^{-4} \text{ dyn}/\text{cm}^2$	2 d submerged + 14 d ALI; monitored to 21 d ALI	24-unit chip	TEER; Cascade Blue/Texas Red permeability; epidermal differentiation markers	Irritation testing; Barrier function	Epidermis-only model with limited multicellular complexity; flow profile and shear stress require further investigation, and more sophisticated biomimetic models with additional irritation biomarkers are needed.	1
Dynamic barrier (Epidermis)	HaCaT	Col I + sponge scaffold	PDMS + PMMA	41.67 $\mu\text{L}/\text{h}$ submerged; 30 $\mu\text{L}/\text{h}$ lower-channel perfusion during ALI	7 d submerged + 35 d ALI	Single unit	H&E-based stratification and stratum corneum; no TEER/permeability	Barrier function; Tissue engineering	Mainly demonstrated culture feasibility and epidermal morphology; no clear advantage was observed for low-density 2D submerged cultures, and further validation of functional barrier readouts and broader cell/tissue applications is needed.	2
Dynamic barrier (Epidermis)	HaCaT	PET membrane	PDMS	NR	15 d for permeation assay; 5 d for viability assay	24-site format	FITC-dextran permeability across HaCaT monolayer	Drug permeation / PK; Irritation testing	Limited to a HaCaT monolayer-based barrier without full epidermal differentiation; long-term culture robustness and cell-layer integrity need improvement for more reliable kinetic permeation assays.	3
Dynamic barrier (Epidermis); with pigmentation	NHKs + NHMs	PET membrane	PMMA	Medium flow 1 $\mu\text{L}/\text{min}$; air flow 1 $\mu\text{L}/\text{min}$	2 d submerged + 14 d ALI	3-unit chip	Cascade Blue/Texas Red permeability; prednisone acetate permeation; ZO-1/differentiation markers	Barrier function; Cosmetic safety	Epidermis-based model remains limited in multicellular tissue complexity; broader validation with more chemicals, inter-laboratory testing, and expanded high-throughput formats is needed for pharmaceutical and cosmetic applications.	4
Dynamic mechanics (Full-thickness)	HEKs + HDFs	PS scaffold + fibroblast-derived matrix	polycarbonate + PDMS	Dermal perfusion 1.5 $\mu\text{L}/\text{min}$; epidermal perfusion 2 $\mu\text{L}/\text{min}$; air flow 1.5 $\mu\text{L}/\text{min}$	6 d dermis + 2 d submerged + 11 d ALI	Single unit	TEER; FITC-dextran permeability; epidermal stratification and barrier function	Barrier function; Drug testing	Dermal maturation and FDM deposition require careful optimization to avoid keratinocyte infiltration or channel clogging; future work should integrate additional cell types and expand disease- or patient-specific validation.	5
Dynamic mechanics (Full-thickness)	HEKs + HDFs	PS scaffold + fibroblast-derived matrix	PS + PDMS	All at 1.5 $\mu\text{L}/\text{min}$	12 d dermis + 2 d submerged + up to 14 d ALI	Single unit	TEER; lucifer yellow penetration	Barrier function; Irritation testing	Mainly validated as a proof-of-concept barrier-sensing platform; TEER sensitivity depends on electrode geometry and chamber design, and broader validation with diverse barrier tissues and testing scenarios is needed.	6
Dynamic mechanics (Full-thickness)	HaCaT + HDFs	Fibrin gel + PC membrane	vinyl + PMMA + PDMS	lower-channel flow 40 $\mu\text{L}/\text{h}$; keratinocyte seeding 50 $\mu\text{L}/\text{min}$	24 h fibroblast spreading + 6 h keratinocyte attachment	Single unit	Multilayer tissue formation and cell viability by fluorescence/confocal imaging; no TEER/permeability assay	Tissue engineering; Multilayer tissue generation	Proof-of-concept platform mainly demonstrates multilayer tissue generation rather than mature skin barrier function; further epidermal maturation, functional readouts, and mechanical stimulation strategies are needed.	7
Dynamic mechanics (Full-thickness)	HEKs + HDFs	Col I + PET membrane	PDMS	Gravity-driven rocker flow; average shear stress $\sim 0.06 \text{ Pa}$	3 d dermis + 3 d submerged + 11 d ALI	24-unit plate; scalable to 128 units	TEER; dextran permeability; LOR and collagen IV staining	Barrier function; Toxicity screening	Highly scalable platform, but TEER sensitivity remains non-uniform and requires further optimization of electrode integration, frequency range, and electrode geometry; tissue-specific maturation and application validation should be expanded.	8

Model type /dynamic level	Cells	Matrix / scaffold	Materials	Flow rates	Culture duration	Throughput	Barrier Metrics	Primary application	Limitations / future need	Ref.
Dynamic mechanics (Full-thickness)	HaCaT + fibroblasts	Col I + PET membrane	PDMS	Medium flow 1 $\mu\text{L}/\text{min}$; air flow 1 $\mu\text{L}/\text{min}$; shear stress $3 \times 10^{-4} \text{ dyn}/\text{cm}$	2 d submerged + day 3–16 ALI	Single unit	TEER; fluorescent molecule permeability; stratum corneum thickness; collagen IV/FLG/LOR	Inflammation / infection; Drug efficacy	The platform remains based on a simplified HaCaT/fibroblast collagen model and a P. acnes/SLS-induced injury setting; future work should incorporate more native skin cell types, optimize dermal stiffness, and expand multi-chamber or high-throughput drug-screening validation.	9
Dynamic mechanics (Full-thickness); with stretch	HEKs + HDFs	Col I	PDMS	10 $\mu\text{L}/\text{min}$; tensile stimulation: 0.01 Hz, 5.3 mm/min	28 d	Single unit	H&E morphology and barrier-related proteins; no TEER/permeability assay	Mechanobiology; Wound-healing-like modeling	Mainly demonstrates tensile-stimulation-induced skin responses rather than a validated disease or drug-testing model; loading parameters and long-term tissue damage need further optimization for damaged-barrier or wound-healing applications.	10
Dynamic mechanics (Full-thickness); with stretch	HEKs + HDFs	Col I + sulfo-SANPAH	PDMS	Gravity-driven operation; mechanical stimulation: 5% compressive strain @ 0.01 Hz, 12 h on/12 h off	5 d dermis + 3 d epidermis + 28 d ALI/mechanical aging	Single unit	H&E epidermal thickness and COL4A/fibronectin/FLG/K10/IVL markers; no TEER/permeability assay	Skin aging model; Cosmetic evaluation	The ageing model remains simpler than native skin because it is limited by current fibroblast–keratinocyte co-culture conditions; future studies should add immune cells, capillaries, peripheral nerves, hair follicles, and sweat glands for more physiologically complete ageing models.	11
Dynamic mechanics (Full-thickness); with stretch	HEKs + HDFs	Col I + porous membrane	PDMS	NR; Stretching: 10% strain @ 0.01 or 0.05 Hz; 12 h/day	4 d fibroblasts + 4 d keratinocytes + 7 d ALI/stretch	Single unit	H&E morphology, stratum corneum thinning, and collagen IV/fibronectin/K10 markers; no TEER/permeability assay	Anti-wrinkle evaluation; Mechanobiology	The model establishes wrinkle-like morphology under rapid cyclic stretching but remains an early ageing-mimic platform; stretch frequency and strain require careful optimization because higher-frequency stimulation can damage epidermal and dermal structures.	12
Dynamic mechanics (Full-thickness); with stretch	HEKs + HDFs + HUVECs	Col I with hollow vascular channels	PDMS + Ecoflex + 3D-printed resin	Perfusion 2–3 mL/h; Stretching: 10% strain @ 0.05 Hz	2 d dermis + 2 d endothelialization + 3 d ALI + 4 d stretch	Single unit	Epidermal capacitance measurement; epidermal thickness/stratification; collagen IV and Ki67 staining	Mechanobiology; Transdermal testing	Mainly demonstrates that combined perfusion and dynamic stretching can improve skin-equivalent morphology; further optimization of perfusion/stretching parameters and validation of barrier function, disease modelling, and drug-testing applications are needed.	13
Dynamic circulation (immune)	Dermal fibroblasts + HUVECs + M1/M2 macrophages	Matrigel + collagen/fibronectin coating	PDMS	NR	24-48 h	Single unit	NR	Wound healing; Immune response modeling	The model captures the paracrine component of early wound inflammation but does not reconstruct a full-thickness epidermal barrier or all wound-healing stages; future work should incorporate additional cell types and extend the platform to later proliferation and remodelling phases	14
Dynamic circulation (Vascularized & immune)	HDFs + HaCaT + HUVECs + HL-60	Col I + porous membrane	PDMS	10 $\mu\text{L}/\text{min}$	7 d dermis + 3 d HaCaT stabilization + 7–10 d differentiation	Single unit	Epidermal differentiation: K5/FLG; endothelial barrier: 70 kDa FITC-dextran permeability and CD31	Immune response modeling; Irritation testing	The vascularized skin chip recapitulates leukocyte recruitment after UV stimulation, but vascular permeability was higher than in Transwell and endothelial differentiation was partly affected by chip culture conditions; further optimization of medium compatibility, endothelial stability, and disease-specific immune validation is needed.	15
Dynamic circulation (Vascularized & immune)	HEKs + HDFs + dermal microvascular endothelial cells + primary neutrophils	Col I	PMMA/glass	Gravity-driven microvascular perfusion; flow rate NR	2 d endothelialization + 2 d keratinocyte confluence + 7–10 d ALI;	Single unit	Epidermal markers: K14/K10/FLG/LOR/IVL; endothelial barrier: 40 kDa FITC-dextran permeability and CD31/VE-cadherin	Inflammation / infection; Immune response modeling	The model is limited by mixed-origin commercial primary cells and lacks skin-resident immune components; next-generation platforms should use autologous cellular components and incorporate tissue-resident immune cells to model adaptive immunity and personalized HSV responses.	16

Model type /dynamic level	Cells	Matrix / scaffold	Materials	Flow rates	Culture duration	Throughput	Barrier Metrics	Primary application	Limitations / future need	Ref.
Dynamic circulation (Vascularized & immune)	HEKs + HDFs + HUVECs + HL-60	Porous scaffold + Col I	NR	Gravity-driven perfusion; flow rate NR	Long-term ALI culture	Single unit	No TEER/permeability assay; stratum corneum/epidermal thickness and FLG/IVL/LOR markers; CD31/CD144 and microsphere perfusion	Inflammation / immune response modeling; Drug testing	Although the platform improves contraction control and vascularized AD-like modelling, pericyte incorporation, stepwise higher-shear maturation, broader cytokeratin/stromal marker panels, and validation with primary immune subsets and disease-specific therapeutic benchmarks are still needed.	17
Dynamic circulation (Vascularized)	HaCaT + HS27 + HUVECs	Fibronectin-coated PET membranes (0.4 μ m)	PDMS	Gravity-driven perfusion; flow rate NR	3 d to confluence; up to 3 weeks culture stability	Single unit	Endothelial barrier: ZO-1 and 4 kDa FITC-dextran permeability	Inflammation / infection; Drug testing	The model simulates TNF- α -induced inflammation and edema but remains based on simplified monolayer cell layers using HaCaT, fibroblasts, and HUVECs; future work should improve epidermal maturation, immune-cell integration, and broader drug or irritant validation.	18
Dynamic circulation (Vascularized)	HDFs + HEKs / HaCaT; HUVECs	Col I + PC membrane	PDMS	Gravity-driven perfusion; 10 μ L/min	5-7 d fibroblasts + 2 d keratinocytes + 5-10 d ALI	Single unit	H&E stratification and K5/IVL/FLG markers; no TEER/permeability assay	Barrier function; Tissue engineering	Mainly establishes on-chip skin construction and perfusion-dependent cell viability; collagen contraction, gel thickness-dependent nutrient limitation, and lack of immune-function validation remain to be further improved for physiological skin modelling.	19
Dynamic circulation (Vascularized)	HaCaT + Hs27 fibroblasts + HUVECs	Porous membranes	PDMS	Gravity-driven perfusion; flow rate NR	3 d to confluence + 42 h post-exposure	Single unit	MTT viability and ZO-1 tight-junction ratio; no TEER/permeability assay	Irritation testing; Toxicity screening	The dual-parameter model improves irritation prediction by combining cell viability and endothelial tight-junction disruption, but validation remains limited to 20 reference substances; more chemicals, human patch-test datasets, and standardized inter-laboratory testing are needed.	20
Dynamic circulation (Vascularized)	HDFs + HEKs + HUVECs	Fibrin gel	PDMS	NR	6 d total (2 d pre-culture + 4 d treatment)	Single unit	Claudin-1/K14 and live/dead staining for keratinocyte integrity; no TEER/permeability assay	Irritation testing; Cosmetic safety	Proof-of-concept irritation model based on angiogenic responses rather than a mature full-thickness epidermal barrier; broader validation with diverse irritants, allergens, corrosives, and standardized quantitative biomarkers is needed for cosmetics and drug-safety testing.	21
Dynamic circulation (Vascularized)	N/TERT keratinocytes + HCA2 fibroblasts + HUVECs + pericytes	Fibrin + Col I/Matrigel	PDMS	NR	4 d vascularization + 14 d skin culture	Single unit	No TEER/permeability assay; epidermal stratification and laminin 332/K14/IVL/TGM1 markers; 70 kDa FITC-dextran vascular perfusability	Drug testing; Vascularized skin equivalent	The model improves vascularized skin-equivalent construction, but epidermal maturation under ALI remains limited and the tissue is less stratified than native skin; future work should optimize nutrient transport, topical/systemic delivery assays, and scalable multi-chip formats for higher-throughput therapeutic testing.	22
Systemic dynamics (Multi-organ)	HEKs + hNSC-derived neurons / hiPSC-HEPs	Col I	PDMS	NR	3 d stabilization + ALI to day 7	6-unit chip	Cornified epidermis with IVL/CK10/CK5 markers; no TEER/permeability assay	Toxicity screening; multi-organ toxicity	The platform enables skin-nerve and skin-liver coupling, but remains a proof-of-concept system with a simplified epidermal barrier and limited test compounds; broader chemical validation, quantitative barrier/permeation assessment, and standardized screening workflows are needed.	23

Model type /dynamic level	Cells	Matrix / scaffold	Materials	Flow rates	Culture duration	Throughput	Barrier Metrics	Primary application	Limitations / future need	Ref.
Systemic dynamics (Multi-organ)	iPSC-cardiomyocytes + primary hepatocytes	Strat-M membrane	PMMA + PDMS	Rocker-driven; flow rate NR	Acute 3 d; chronic topical 10 d	Single unit	Strat-M membrane drug permeation by LC-MS/MS; no living skin barrier metrics	Drug permeation/PK; multi-organ toxicity	The system links dermal absorption to heart–liver toxicity, but uses a synthetic Strat-M skin surrogate rather than a living skin equivalent; future work should incorporate biological skin barriers, skin heterogeneity, and broader pharmacokinetic validation of topical drug delivery.	24
Systemic dynamics (Multi-organ)	Caco-2 + HEKs + gut microbes	PET membrane + PC membrane	PC	Gravity-driven perfusion; 15° tilt, 500 s interval	14 d gut differentiation + skin ALI differentiation	Single unit	Gut barrier: TEER, 70 kDa FITC-dextran permeability, and ZO-1; reconstructed epidermis: H&E and K5/IVL/FLG/Ki67 markers.	Gut–skin axis modeling; Inflammation / infection	The chip captures gut microbiota–skin crosstalk, but relies on Caco-2 cells and lacks immune components; future iterations should incorporate primary or patient-derived intestinal/skin cells, immune cells, and broader microbiome or disease-specific validation.	25
Systemic dynamics (Multi-organ)	<i>Ex vivo</i> human epidermis + KGN or TM3 aggregates	<i>Ex vivo</i> epidermis	PDMS	NR	48 h aggregate formation + 24 h stimulation	Single unit	H&E epidermal structure and stratum corneum thickness; no TEER/permeability assay	Cosmetic safety; Hormone crosstalk model	Sex-hormone regulation is modelled mainly through local hormone gradients without vascular delivery, limiting systemic hormone distribution and vascular–epidermal interactions; future work should incorporate vascular structures and broader donor or hormone-condition validation.	26

PMMA: polymethyl methacrylate; PDMS: polydimethylsiloxane; PS: polystyrene; PC: polycarbonate; PET: polyethylene terephthalate; HEKs: human epidermal keratinocytes; HDFs: human dermal fibroblasts; HUVECs: human umbilical vein endothelial cells; Col I: type I collagen; ALI: air-liquid interface; TEER: transepithelial electrical resistance; FLG: filaggrin; IVL: involucrin; LOR: lorixin; NR: not reported.

References

1. J. Zhang, Z. Chen, Y. Zhang, X. Wang, J. Ouyang, J. Zhu, Y. Yan, X. Sun, F. Wang, X. Li, H. Ye, S. Sun, Q. Yu, J. Sun, J. Ge, Q. Li, Q. Han, Y. Pu and Z. Gu, *Lab Chip*, 2021, **21**, 3804-3818.
2. Q. Lang, Y. Ren, Y. Wu, Y. Guo, X. Zhao, Y. Tao, J. Liu, H. Zhao, L. Lei and H. Jiang, *Rsc Adv*, 2016, **6**, 27183-27190.
3. N. Sasaki, K. Tsuchiya and H. Kobayashi, *Sensor Mater*, 2019, **31**, 107-115.
4. Q. Li, C. Wang, X. Li, J. Zhang, Z. Zhang, K. Yang, J. Ouyang, S. Zha, L. Sha, J. Ge, Z. Chen and Z. Gu, *J Tissue Eng*, 2023, **14**, 1-13.
5. P. C. A. Zoio, S. Lopes-Ventura and A. Oliva, *Micro*, 2022, **2**, 191-211.
6. P. Zoio, S. Lopes-Ventura and A. Oliva, *Micromachines-Basel*, 2021, **12**, 20.
7. L. Valencia, V. Canalejas-Tejero, M. Clemente, I. Fernaud, M. Holgado, J. L. Jorcano and D. Velasco, *Sci Rep-Uk*, 2021, **11**, 13160.
8. H. T. Nguyen, S. L. Rissanen, M. Peltokangas, T. Laakkonen, J. Kettunen, L. Barthod, R. Sivakumar, A. Palojärvi, P. Junttila, J. Talvitie, M. Bassis, S. L. Nickels, S. Kalvala, P. Ilina, P. I. Tammela, S. Lehtonen, J. C. Schwamborn, S. Mosser and P. Singh, *Tissue Barriers*, 2024, **12**, 2315702.
9. Q. Quan, D. Weng, X. Li, Q. An, Y. Yang, B. Yu, Y. Ma and J. Wang, *Front Bioeng Biotech*, 2022, **10**, 13.
10. K. Kim, S. Jeong and G. Y. Sung, *Biochip J*, 2022, **16**, 501-514.
11. S. Jeong, J. Kim, H. M. Jeon, K. Kim and G. Y. Sung, *International Journal of Molecular Sciences*, 2021, **22**, 12788.
12. H. Y. Lim, J. Kim, H. J. Song, K. Kim, K. C. Choi, S. Park and G. Y. Sung, *J Ind Eng Chem*, 2018, **68**, 238-245.
13. N. Mori, Y. Morimoto and S. Takeuchi, *Biofabrication*, 2019, **11**, 11001.
14. S. Biglari, T. Y. L. Le, R. P. Tan, S. G. Wise, A. Zambon, G. Codolo, M. De Bernard, M. Warkiani, A. Schindeler, S. Naficy, P. Valtchev, A. Khademhosseini and F. Dehghani, *Adv Healthc Mater*, 2019, **8**, e1801307.
15. B. S. Kwak, S. P. Jin, S. J. Kim, E. J. Kim, J. H. Chung and J. H. Sung, *Biotechnol Bioeng*, 2020, **117**, 1853-1863.
16. S. Sun, L. Jin, Y. Zheng and J. Zhu, *Nat Commun*, 2022, **13**, 15.
17. K. Kim, H. Jang, E. Kim, H. Kim and G. Y. Sung, *Lab Chip*, 2025, **25**, 5936-5949.
18. M. Wufuer, G. Lee, W. Hur, B. Jeon, B. J. Kim, T. H. Choi and S. Lee, *Sci Rep-Uk*, 2016, **6**, 37471.
19. S. Lee, S. P. Jin, Y. K. Kim, G. Y. Sung, J. H. Chung and J. H. Sung, *Biomed Microdevices*, 2017, **19**, 22.
20. B. Jeon, G. Lee, M. Wufuer, Y. Huang, Y. Choi, S. Kim and T. H. Choi, *Toxicol in Vitro*, 2020, **68**, 104955.
21. N. Jusoh, J. Ko and N. L. Jeon, *Apl Bioeng*, 2019, **3**, 36101.
22. C. F. E. Jones, S. Di Cio, J. T. Connelly and J. E. Gautrot, *Front Bioeng Biotech*, 2022, **10**, 915702.
23. J. S. Lee, J. Kim, B. Cui, S. K. Kim, S. A. Cho, S. An and S. W. O. Cho, *Lab Chip*, 2022, **22**, 343-353.
24. C. P. Pires De Mello, C. Carmona-Moran, C. W. McAleer, J. Perez, E. A. Coln, C. J. Long, C. Oleaga, A. Riu, R. Note, S. Teissier, J. Langer and J. J. Hickman, *Lab Chip*, 2020, **20**, 749-759.
25. B. Ko, J. Son, J. In Won, B. M. Kang, C. W. Choi, R. Kim and J. H. Sung, *Lab Chip*, 2025, **25**, 2609-2619.
26. X. Wu, Y. Zhang, Y. Xu, X. Cao and Y. Zhao, *Bioact Mater*, 2025, **53**, 240-252.