

Deep Eutectic Solvent-Driven Synthesis of La_2NiO_4 @rGO Ruddlesden–Popper Hybrids for Ethyl Parathion Sensing Supported with Theoretical Studies

Farhana Yasmin Rahman^{a,b,c,d}, Rajendran Surya^{d,e}, Subramanian Sakthinathan^{d,e}, Sangita Majumder^a, Syed Arshad Hussain^{a*}, Te-Wei Chiu^{d,e*} Mani Govindasamy^{b,c*}

^aThin Film and Nanoscience Laboratory, Department of Physics, Tripura University (Suryamaninagar-799022), Tripura, India.

^bInternational Ph. D Program in Innovative Technology of Biomedical Engineering & Medical Devices, Ming Chi University of Technology, New Taipei City 243303, Taiwan.

^cResearch Center for Intelligence Medical Devices, Ming Chi University of Technology, New Taipei City 243303, Taiwan.

^dDepartment of Materials and Mineral Resources Engineering, National Taipei University of Technology, No. 1, Sec. 3, Zhongxiao E. Rd., Taipei 10608 Taiwan.

^eInstitute of Materials Science and Engineering, National Taipei University of Technology, No. 1, Sec. 3, Zhongxiao E. Rd., Taipei 10608 Taiwan.

Corresponding Author

Syed Arshad Hussain: sa_h153@hotmail.com, sahussain@tripurauniv.ac.in

ManiGovindasamy: govindasamy420700@gmail.com,
manymany2025@mail.mcut.edu.tw

Te-Wei Chiu: tewei@ntut.edu.tw

S.I. 1 Chemical reagents:

Lanthanum (III) nitrate hexahydrate ($\text{La}(\text{NO}_3)_3 \cdot 6\text{H}_2\text{O} \geq 99.0\%$), nickel (II) nitrate hexahydrate ($\text{Ni}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O} \geq 99.0\%$), (ACS grade), ethanol ($\text{C}_2\text{H}_5\text{OH} \geq 99.8\%$), sodium hydroxide ($\text{NaOH} \geq 99.0\%$), sodium phosphate monobasic (NaH_2PO_4), sodium phosphate dibasic (Na_2HPO_4), alumina slurry, and ethyl parathion (EP) were procured from Sigma-Aldrich, Taiwan. Additionally, Choline chloride ($[(\text{CH}_3)_3\text{NCH}_2\text{CH}_2\text{OH}]^+\text{Cl}^-$), fructose ($\text{C}_6\text{H}_{12}\text{O}_6$) and citric acid ($\text{C}_6\text{H}_8\text{O}_7$), and potassium hydroxide (KOH) were procured from Sigma-Aldrich, Taiwan. All chemicals were of analytical grade with the highest purity available and did not require additional purification. A phosphate buffer solution (PBS) at a concentration of 0.1 M was prepared using Na_2HPO_4 and NaH_2PO_4 , and pH adjustment was achieved using 0.1 M HCl/NaOH solutions. Distilled water served as the solvent for formulation. Finally, carrot juice, apple juice, cabbage, and orange juice have been purchased in the Taiwan market for the real sample analysis.

S.I. 2 Chemical Methods and Instruments:

High-Resolution Transmission Electron Microscopy (HR-TEM; JEM 2100F) and Field Emission Scanning Electron Microscopy (FESEM, JEOL JSM-6500F) coupled with Energy-Dispersive X-ray Spectroscopy (EDX) were employed to identify the surface morphology of nanomaterials. Furthermore, X-ray Photoelectron Spectroscopy (XPS) with the Thermo Scientific MultiLab 2000 instrument and X-ray Diffraction (XRD, PANalytical X'PERT PRO, Netherlands) using $\text{Cu-K}\alpha$ as the source radiation ($\lambda = 1.5417 \text{ \AA}$), Raman analysis was performed within the frequency range of 100 cm^{-1} to 3000 cm^{-1} utilizing a Dong Woo 500i system from Korea and Fourier transform infrared spectroscopy (FTIR) has been studied by the FTIR-6600 spectrometer in the wavenumber range between 4000 and 400 cm^{-1} were utilized to scrutinize the elemental composition and surface characteristics. In the realm of electrochemical studies, a three-electrode system was utilized, comprising a glassy carbon electrode (GCE, 0.071 cm^2) as the working electrode, a platinum electrode as the counter electrode, and an Ag/AgCl electrode as the reference electrode (ALS Co., Ltd, Japan). Cyclic Voltammetry (CV) and Differential Pulse Voltammetry (DPV) were conducted in a nitrogen gas atmosphere using phosphate buffer solution, with a potential range of 0.6 to -1.2 V (vs Ag/AgCl) and a scan rate of 100 mV/s . Additionally,

Electrochemical Impedance Spectroscopy (EIS) was carried out over a frequency range of (10^2 - 10^5) Hz in a solution containing potassium chloride and $[\text{Fe}(\text{CN})_6]^{3-/4-}$. These electrochemical methods provided valuable insights into the material's electrochemical behavior and properties. Real sample validation carried out by Ultimate 3000 RS UHPLC Systems with Ultimate 3000 Analytical column.




Figure S1. Elemental mapping of La₂NiO₄@rGO nanocomposite. (a) La₂NiO₄@rGO nanocomposite, (b) whole elemental distribution, (c) La, (d) C, (e) Ni, (f) O.




Figure S2. Anti-interference study with five-fold excess concentration




Figure S3. Calibration plot for the real samples, carrot juice (a), cabbage juice (b), apple juice (c), orange juice (d).